User Guide

Revision 1 (August 2024)

Document No.: ME-UG-7022

CONTENTS

02	English	354	日本語
24	العربية	376	한국어
46	Български	392	Latviešu
68	中文(简体)	414	Lietuvis
90	Hrvatski	436	Norsk
112	čeština	458	Polskie
134	Dansk	480	Português
156	Dutch	502	Română
178	Eesti keel	524	русский
200	Suomi	546	Slovenčina
222	Français	568	Slovenščina
244	Deutsch	590	Español
266	Ελληνικά	612	Svenska
288	magyar	634	Türk
310	íslenska	656	українська
332	Italiano		

English

English

About This Guide	4
1 General	4
1.1 Intended Use	4
1.2 Workflow	4
1.3 Scanner & Software	4
1.4 Qualification of Operating User	5
1.5 Symbols	5
1.6 Components Overview	6
1.6.1 Additional Components (Sold Separately)	7
1.6.2 How to Use Components	8
2 Medit Scan for Labs Overview	9
2.1 Introduction	9
2.2 Installation	9
2.2.1 System Requirements	9
2.2.2 Software Installation Guide	10
2.2.3 Hardware Configuration	10
2.3 Medit Scan for Labs User Guide	11
2.4 Scanner Calibration	11
3 Safety Guide	14
3.1 General Precautions	14
3.2 Environmental Conditions	16
3.3 Electrical Safety	16
3.4 Electromagnetic Compatibility Information	17
3.4.1 Electromagnetic Emissions	17
3.4.2 Electromagnetic Immunity	18
4 Specifications	22

About This Guide

Convention in This Guide

This guide uses various symbols to highlight important information to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

WARNING

The WARNING symbol indicates information that, if ignored, could result in a medium risk of personal injury.

CAUTION

The CAUTION symbol indicates safety information that, if ignored, could result in a slight risk of personal injury, property damage, or damage to the system.

1. General

1.1 Intended Use

The Medit T-Series system is a table top dental 3D scanner intended to record topographical characteristics of teeth models and restorations. The Medit T-Series system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations.

1.2 Workflow

The workflow is designed to provide high-quality scan data in the dental clinic or laboratory for any shape and size.

- Model or Impression Scan
Medit Scan for Labs will scan the model according to the information entered in the order form in Medit Link.
This enables you to create a prosthesis directly by scanning the impressions compared to the conventional methods of prosthesis manufacturing.
- CAD Processing
Design the prosthesis using a CAD program.
- CAM Processing
Convert the designed prosthesis into NC data using a CAM program.
- Manufacturing
Fabricate the prosthesis using a machine according to the NC data.
- Finishing
Perform the finishing on the prosthesis.

1.3 Scanner & Software

The scanner comes equipped with the accompanying software.

- Scanner: Medit Table Top Scanner (Medit T-Series)
The scanner is designed to acquire scan data from a variety of dental models and impressions in a convenient way.
A full arch scan takes only 8 seconds (T500 takes 12 seconds).
- Software: Medit Scan for Labs
The accompanying software is designed to be user-friendly, making it easy to acquire scanned data.

1.4 Qualification of Operating User

The system can only be used by trained dental professionals or technicians.

You are solely responsible for the accuracy and completion of all the data acquired using your 3D scanner system.

The user should verify the accuracy of each scan result and use it to evaluate the applicability of each treatment.

The scanner system must be used in accordance with the accompanying user manual.

Improper use or handling of the scanner system will void your warranty.

If you need additional information or assistance in using the equipment, please contact your local service provider.


You cannot modify or change the device of the software system on your own.


1.5 Symbols

No.	Symbol	Description
1		Serial number
2		Date of manufacture
3		Manufacturer
4		Caution
5		Warning
6		Caution. Hand hazard and optical hazard
7		Read the user guide
8		The official mark of the Europe Certificate
9		Authorized representative in the European community
10		WEEE mark
11		MET mark
12		AC
13		DC
14		Protective Earth
15		Temperature limitation: -5 – 45°C (23 – 113°F)
16		Humidity limitation

17		Atmospheric pressure limitation
18		Fragile
19		Keep dry
20		This way up
21		Stacking more than three layers is prohibited
22		Consult instructions for use

1.6 Components Overview

No.	Item	Qty	Appearance
1	3D Table Top Scanner	1ea	
2	USB Cable	1ea	
3	Power Cable & External Adapter	1ea	
4	LED Protector	1ea	
5	Calibration Panel	1ea	
6	Articulator Plate	1ea	


7	USB Flash Drive (Installer of Medit Scan for Labs included)	1ea	
8	Blu Tack	1ea	
9	Flexible Multi-Die	1ea	
10	Single Die	16ea	
11	Spray Supporter	2ea	
12	Impression Jig (Not available for T510, T310, and T4)	1ea	

1.6.1 Additional Components (Sold Separately)


No.	Item	Qty	Appearance
1	KAS Jig	1ea	
2	AM Jig	1ea	
3	Screw Jig	1ea	

1.6.2 How to Use Components


- Fill each of the spray supporters with Blu Tack.


- Prepare Blu Tack to apply on the surface of the Flexible Multi-die.

- Apply the scan spray onto an object on the spray supporter and fix it to the Flexible Multi-Die for scanning.

- Use the LED protector to cover the light of the scanner projector.

2. Medit Scan for Labs Overview

2.1 Introduction

Medit Scan for Labs is a software program that allows users to perform model and impression scans using Medit's scanner series. Users can edit data, complement it with data from the intraoral scanner and prepare for CAD/CAM processes. Explicit explanations and guide messages for each step can be found on the left side of the window. Medit Scan for Labs is to be run only on computers that meet the specifications outlined in the System Requirements. Otherwise, the device may not function properly. In case Windows is not updated before the installation, USB 3.0 will not work properly.

CAUTION

- This device is designed for USB 3.0 port only. Please make sure to connect it to a USB 3.0 port on your computer.
- This device is compatible only with Windows 10 and later. It does not work with Mac operating systems.
- Before installing the scanning S/W, please make sure that the Windows version in use, the mainboard, the VGA card, and the USB drivers are up to date.

2.2 Installation

2.2.1 System Requirements

Minimum System Requirements

	Laptop	Desktop
CPU	Intel Core i7-8750H or higher	Intel Core i7-8700K or higher
RAM	16 GB or higher	
Graphic	NVIDIA GeForce GTX 1060 6 GB or higher	
OS	Windows 10 64-bit / Windows 11 64-bit	

Recommended System Requirements

	Laptop	Desktop
CPU	Intel Core i7-8750H or higher	Intel Core i7-8700K or higher
RAM	32 GB or higher	
Graphic	NVIDIA GeForce GTX 1060 6 GB or higher	
OS	Windows 10 64-bit / Windows 11 64-bit	

2.2.2 Software Installation Guide

- ① Connect the enclosed USB flash drive to a PC.
- ② Run the installation file.
- ③ Select the setup language and click "Next."
- ④ Choose the installation path.
- ⑤ Read the "License Agreement" carefully, check "I agree to the License terms and conditions," and then click "Install."
- ⑥ The installation process may take several minutes. Please do not shut down the PC until the installation is complete.
- ⑦ After the installation is complete, restart the PC to ensure optimal program operation.

The installation will not be processed if the Medit T-Series is connected to a PC. Please turn off the scanner before starting the installation.

2.2.3 Hardware Configuration

Once the software installation is complete, reboot your PC before installing the hardware.


CAUTION


The package includes a power cable and a USB cable. All cables used with the scanner must be properly connected to the PC.

* Use only a USB 3.0 port when connecting the scanner to your PC.


How to Connect T710/T510/T310

- ① Plug in the scanner power cable and connect the USB 3.0 cable to a USB 3.0 port.

- ② Press the power button to turn on the scanner.

2.3 Medit Scan for Labs User Guide

Please refer to the User Guide of Medit Scan for Labs: Medit Scan for Labs > Menu > User Guide.

2.4 Scanner Calibration

Calibration is recommended for proper scanning and performance of the device. Please calibrate the scanner when:

- The quality of scan data has decreased as compared to the previous scans.
- The external conditions, such as device temperature, changed during use.
- It is already past the configured calibration period.

It is recommended to calibrate the device periodically. Go to Menu > Settings > Table Top Scanner, and configure the calibration period in Calibration Period (Days) option. The default calibration period is 30 days.


CAUTION


The calibration panel is a delicate component. Please do not touch it.
If the calibration fails, inspect the panel and contact the service provider if it is contaminated.

Calibration of T710/T510/T310



- ① Turn on the scanner and connect the scanner to the software.
- ② Click on the scanner icon at the bottom left to run the Calibration Wizard.
- ③ Prepare and place the calibration panel.
- ④ Select one of the two calibration options and click Next.
 - Auto calibration: Automatic calibration is conducted with the QR code on the back side of the calibration panel.
 - Manual calibration: The corresponding PNL file is required to conduct manual calibration.

- ⑤ Please enter the serial number of the calibration panel according to the option you chose above.
 - Auto calibration
 - » The scanner will scan the QR code on the back side of the calibration panel, and the calibration process starts automatically.
 - Manual calibration
 - » Check the serial number on the calibration panel and select the corresponding PNL file from the file list.
 - » If you cannot find the serial number on the list, please check if you have a PNL file on the PC or installation USB.
 - If you have a PNL file, click to search for it.
 - If you do not have a PNL file, click and enter the serial number.

- ⑥ The calibration process might take a few minutes. Please do not touch the scanner.
- ⑦ Wait until the calibration is completed successfully.

Intraoral Scanner Calibration

- ① Turn on the intraoral scanner and launch the Medit Scan for Clinics.
- ② Run the Calibration Wizard at the bottom of the main toolbar panel in Medit Scan for Clinics.
- ③ Prepare the calibration tool and the handpiece.
- ④ Turn the dial of the calibration tool to the starting position.
- ⑤ Put the handpiece into the calibration tool.
- ⑥ Click "Next" to start the calibration process.
- ⑦ If the calibration tool is mounted properly in the correct position, the system automatically acquires data.
- ⑧ When data acquisition is complete at the starting position, turn the dial to the next position.
- ⑨ Repeat the steps to the last position.
- ⑩ When data acquisition is complete at the last position, the system automatically calculates and shows the calibration results.

3. Safety Guide

3.1 General Precautions

- This scanner should only be used by trained professionals.
- Before operating the scanner, users should receive training on how to use the system or thoroughly review and understand this User Guide.
- Do not spray any liquid inside the scanner.
- Do not touch the lenses.
- Do not operate the scanner with the cover on.
- Do not store items on the top of the scanner.

During Installation

- Install the device in a dust-free environment with proper ventilation and minimum changes in air pressure, temperature, and humidity.
- Please be aware of the safety conditions, such as tilting the device, vibration, or shock.
- Do not install the device in a place with chemical storage or near gas-generating points.
- Do not install the device in a place with poor ventilation.
- Please consider the power requirements and consumption.
- Ensure that all packaged components are free from physical damage. Safety cannot be guaranteed if there is any physical damage to any component of the units.
- Install and use only approved programs to ensure proper functionality of the 3D scanner system.
- Hand Hazard

- » Refrain from touching while in operation.
- » To prevent access to the moving part, maintain a safe distance of 1.5 m from the scanner to keep you and nearby equipment, including the PC.

Before Using Equipment

- Please make sure that all the components and cables are properly connected.
- Ensure that the device appears in the Device Manager.
- Exercise caution and refrain from applying force when manually adjusting the scanner's axis.

While Using Equipment

- Make sure that the scanner is not subjected to shock during operation.
- Do not turn off the scanner while in use.
- Ensure that none of the openings on the device are blocked during operation.
- If you detect any smoke or unusual odor from the scanner, immediately turn it off, disconnect the cables, and contact the manufacturer for assistance.
- Unplug the power cable if the scanner will be unattended or unused for a long period of time.
- If there is a problem, do not attempt to disassemble and assemble the scanner on your own. Once a problem is identified, please contact your local authorized resellers for repair assistance.

- Please refrain from altering the product in any way.
- If the scanner fails to function properly (e.g., experiencing accuracy issues), stop using the product and contact the manufacturer or authorized resellers.
- Manipulation of the scanner should only occur while the moving part is not in operation.
- Do not manipulate the scanner during scanning or while moving parts are in operation, and maintain a distance of 1.5 m from the scanner.
- Position the device in a manner that facilitates easy management and detachment of accessories.
- Use the scanner strictly in accordance with the manufacturer's specifications to ensure the integrity of its protection features.
- Do not place objects weighing over 1.5 kg on the moving part.
- In case there is an error during the scanning process, the scanner will automatically stop, and the LED indicator on the scanner will switch to red.
 - » If the LED indicator on the scanner turns red, click the "Initialize Axis" button in the software to fix the problem.
- LED Caution (Risk Group 2)

- » Avoid staring at the lamp for long periods during operation, as it can be harmful to the eyes.

Maintenance

- If the device was not in use for a long time, please make sure it is properly installed and calibrated before reuse. (The recommended calibration period is one month.)
- Please check if the scanner is acquiring scan data properly.
- Clean the scanner regularly with a dry cloth to prevent moisture from entering it.
- Make sure to unplug the power cable before examining the scanner.
- Do not apply the scanning spray to the scanner or any objects on the moving part.
- Replacement of any parts should only be performed by service personnel with appropriate tools.

Disposal

- All components are designed to conform to the following directives:
 - RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)
 - WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

Storage

- Gently wipe the scanner surface with a dry cloth, ensuring no foreign objects or liquids come into contact with it.
- In case of such contact, immediately wipe the surface without applying force.
- Store the scanner in a safe place to prevent any damage.

3.2 Environmental Conditions

WARNING

Observe the following environmental conditions.

Operating Conditions

- Temperature: 18 – 28°C (64.4 – 82.4°F)
- Humidity: 20 – 75% (Non-condensing)
- Atmospheric pressure: 800 – 1,100 hPa
- Altitude: Up to 2,000 m
- Pollution degree 2
- Indoor use

Storage & Transportation Conditions

- Temperature: -5 – 45°C (23 – 113°F)
- Humidity: 20 – 80% (Non-condensing)
- Atmospheric pressure: 800 – 1,100 hPa (No condensation)

3.3 Electrical Safety

WARNING

- The table top 3D scanner system is a Class I device. It includes the scanner and its components that are described in chapter 1.6 Components Overview altogether.
- To prevent electric shock, the table top 3D scanner system must only be connected to a power source with a protective earth connection. If you cannot insert the plug provided with your table top 3D scanner into a main outlet, contact a qualified electrician to replace the plug or outlet. Do not try to circumvent these safety guidelines.
- Do not use a grounding-type plug connected to the table top 3D scanner system for any other purpose than its intended use.
- Do not install the table top scanner system in a place where it is difficult to unplug the power cord.
- There is a risk of electric shock if you attempt to access the inside of a table top 3D scanner system. Only qualified service personnel should access the system.
- Do not connect your table top 3D scanner system to a regular power strip or extension cord, as these connections are not as safe as grounded outlets.
- Do not place liquids such as beverages near the table top 3D scanner system, and avoid spilling liquid of any kind on the system.
- Condensation due to changes in temperature or humidity can cause moisture buildup inside the table top 3D scanner system, which may damage the system. Before connecting the table top 3D scanner system to a power supply, be sure to keep the table top 3D scanner system at room temperature for at least two hours to prevent condensation. If you see any condensation on the product surface, the table top 3D scanner should be left at room temperature for more than 8 hours.
- You should only disconnect the table top 3D scanner system from the power supply via its power cord.
- When disconnecting the power cord, hold the surface to remove it.
- Before disconnecting, make sure to turn off the power on the device using the power switch on the scanner.
- Avoid pulling on the communication cables, power cables, etc., used with the table top 3D scanner system.
- Use only the adapter provided for the table top 3D scanner. Using other adapters can damage the table top 3D scanner system.
- Do not touch the connectors of the device.

3.4 Electromagnetic Compatibility Information

3.4.1 Electromagnetic Emissions

The Medit T-Series system is intended for use in the electromagnetic environment as specified below. The customer or the user of the Medit T-Series system should ensure that it is used in such an environment.

Guidance and Manufacturer’s Declaration - Electromagnetic Emission		
Emission Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The Medit T-Series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations / Flicker Emissions	Complies	The Medit T-Series is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

WARNING

This Medit T-Series system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the Medit T-Series or shielding the location.

3.4.2 Electromagnetic Immunity

Guidance 1

The Medit T-Series system is intended for use in the electromagnetic environment as specified below.
The customer or the user of the Medit T-Series system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with a synthetic material, relative humidity of at least 30% is recommended.
Electrical Fast Transient / Burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	The mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode	±0.5 kV, ±1 kV differential mode	The mains power quality should be that of a typical commercial or hospital environment.
	±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV, ±2 kV common mode	
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	0% U _T (100% dip in U _T) for 0.5/1 cycles 70% U _T (30% dip in U _T) for 25/30 cycles 0% U _T (100% dip in U _T) for 250/300 cycles	0% U _T (100% dip in U _T) for 0.5/1 cycles 70% U _T (30% dip in U _T) for 25/30 cycles 0% U _T (100% dip in U _T) for 250/300 cycles	The mains power quality should be that of a typical commercial or hospital environment. If the user of the Medit T-Series system requires continued operation during power mains interruptions, it is recommended that the Medit T-Series system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Fields (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a location in a typical commercial or hospital environment.

Proximity Magnetic Fields in the Frequency Range 9 kHz to 13.56 MHz Immunity IEC 61000-4-39	8 A/m 30 kHz CW modulation	8 A/m 30 kHz CW modulation	Resistance to magnetic fields was tested and applied only to surfaces of enclosures or accessories accessible during intended use.
	65 A/m 134.2 kHz PM 2.1 kHz	65 A/m 134.2 kHz PM 2.1 kHz	
	7.5 A/m 13.56 MHz PM 50 kHz	7.5 A/m 13.56 MHz PM 50 kHz	

NOTE: U_T is the main voltage (AC) prior to the application of the test level.

Guidance 2

Recommended Separation Distances Between Portable and Mobile Communication Equipment and the Medit T-Series		
Rated Maximum Output Power of the Transmitter [W]	Separation Distance According to the Frequency of Transmitter [M]	
	IEC 60601-1-2:2014	
	150 kHz to 80 MHz d = 1.2 √ P	80 MHz to 2.7 GHz d = 2.0 √ P
0.01	0.12	0.20
0.1	0.38	0.63
1	1.2	2.0
10	3.8	6.3
100	12	20

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

- Guidance 3**
The Medit T-Series system is intended for use in the electromagnetic environment specified below.
The customer or the user of the Medit T-Series system should ensure that it is used in such an environment.

Guidance and Manufacturer’s Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz Outside ISM Bands amateur	3 Vrms	Portable and mobile RF communications equipment should not be used closer to any part of the Ultrasound System, including cables, than the recommended separation distance. This is calculated using the equation applicable to the frequency of the transmitter. Recommended Separation Distance (d): $d = 1.2 \sqrt{P}$ IEC 60601-1-2:2007 $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80 MHz to 2.5 GHz IEC 60601-1-2:2014 $d = 2.0 \sqrt{P}$ 80 MHz to 2.7 GHz
	6 Vrms 150 kHz to 80 MHz In ISM Bands amateur	6 Vrms	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.
NOTE 3: The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

- Guidance 4**
The Medit T-Series system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the Medit T-Series system. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer’s Declaration - Electromagnetic Immunity					
Immunity Test	Band ¹⁾	Service ¹⁾	Modulation	IEC 60601 Test Level	Compliance Level
Proximity Fields from RF Wireless Communications IEC 61000-4-3	380 – 390 MHz	TETRA 400	Pulse Modulation 18 Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460; FRS 460	FM ±5 kHz Deviation 1 kHz sine	28 V/m	28 V/m
	704 – 787 MHz	LTE Band 13, 17	Pulse Modulation 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM 800:900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	Pulse Modulation 18 Hz	28 V/m	28 V/m
	1700 – 1990 MHz	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217 Hz	28 V/m	28 V/m
	2400 – 2570 MHz	Bluetooth; WLAN 802.11b/g/n; RFID 2450; LTE Band 7	Pulse Modulation 217 Hz	28 V/m	28 V/m
	5100 – 5800 MHz	WLAN 802.11a/n	Pulse Modulation 217 Hz	9 V/m	9 V/m

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

**WARNING**

- Use of the Medit T-Series adjacent to or on other equipment must be avoided as it may result in improper operation. If this use is necessary, it is advisable that this and the other equipment be observed to verify that they are operating normally.
- The use of accessories, transducers, and cables other than those specified or provided by the Medit of the Medit T-Series could result in high electromagnetic emissions or reduced electromagnetic immunity from this equipment and result in improper operation.

¹ For some services, only the uplink frequencies are included.

4. Specifications

Brand Name (Model Name)	T710 (MD-ID0400), T510/T310 (MD-ID0420)
Resolution of Camera	Mono 5.0 (MP)
Point Spacing	0.040 mm
Scan Area	100 x 73 x 60 mm (W x H x D)
Scan Principle	Phase-shifting optical triangulation
Size	271 x 340 x 505 mm (W x H x L)
Weight	15 kg
Light Source	LED, 150 ANSI-lumens
Light Color	Blue LED
Texture Scan	Provide texture scanning
Connection	USB 3.0 B Type
Protection Against Electric Shock	Class I
Power	AC 100 – 240 V, 50 / 60 Hz
Main Supply Voltage Fluctuation	±10%
Power Consumption	Standby: 20 W (Idle: 35W, Scan: 48 W)

Power Requirements	Power supply: 100 – 240 VAC, 50 / 60 Hz
EMC	CE Class A, Passed conducted and radiated test
Protection	OVP (Over Voltage Protection)
	OCP (Over Current Protection)
Mode of Operation	Continuous
DC Adapter (for MD-ID0400, MD-ID0420)	
Model Name	ATM120T-P240
Input Voltage	Universal 100 – 240 VAC, 50 / 60 Hz input, without any slide switch
Output	24 V $\overline{\text{---}}$, 5 A
Case Dimension	168.1 x 65.9 x 39 mm
Protection	OVP (Over Voltage Protection)
	SCP (Short Circuit Protection)
	OCP (Over Current Protection)

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Meditrial Srl
Medit Corp.