

MEDIT Link



Revision 6 (January 2023)



Medit Link

Revision 6 (January 2023)

Table of Contents

Introduction and Overview	1
Cloud Service	2
Workflow Management Using Ordering System	3
Utilizing the Dashboard	3
Installation	4
System Requirements	4
Installation on Windows	5
Installation on macOS	7
Registration	11
Landing Page	11
Registration	11
Login and Welcome Page	13
Medit Link Web	15
System Requirements for Medit Link Web	15
Dashboard	16
Case Box/Work Box	19
Case Talk	21
Creating a Case Talk	22
Web Viewer	26
Data Tree	29
Partners	30
Partner Search	30
Featured Partners	30
My Partner List	31
Partner Invitation	32
Membership	33
Membership	34
Billing Info	39
Redeem	43

Payment	45
Payment Details (Clinics, Labs accounts)	45
Payment account information (Lab account only)	45
My Information	46
My Account	46
My Organization (Admin account only)	46
My Price List (Lab: Admin account only)	47
Notice and Updates	48
Help Center	48
Download the App	48
Medit Link App	49
Title Bar	49
User Information	50
Notifications	51
Case Synchronization Manager	52
Network Status	53
Main Functions	54
Box Overview	54
Case Box/Work Box	55
Order Box/In Box	58
Case Talk	59
App Box	59
Trash Box	60
Settings	60
Patient Management	69
Notices and Updates	71
Go to Medit Link Web	71
Help Center	71
File Viewer	72
Data Tree	72
Model Control	73
Capture and Memo	74
Orthodontic Capture	75
Fill Holes	77
Trimming	77
Optimize Occlusion Alignment	81
Mandibular Movement	83
Manage the Multi Occlusion Data in File Viewer	84

View Options	85
Status Filters	85
Table Items	87
View Type	88
Grouping	89
Sorting	90
Filtering and Search	92
Filtering by Date	92
Advanced Search Capabilities	93
Detailed Explanation of the Process	94
Managing Cases: Creating, Editing, Cloning, Deleting	94
Creating Cases	94
Editing cases	97
Cloning cases	99
Deleting cases (Lab and Clinic)	100
Inputting form information	100
Inputting Form Information in Case Box	100
Inputting Form Information in Work Box	101
Inputting Bridge Information	102
Setting Scanning Options	103
Copying Form Information	104
Deleting treatment information	104
Adding Product Information to “Favorite”	105
Adding to Favorites	105
Managing Favorite Templates	105
Copying form information between Scan and CAD tabs	106
From Scan tab to CAD tab	106
From CAD tab to Scan tab	107
Scanning process	108
CAD Process	108
Get Result Files	109
Export scan result files	109
Export selected files	111
Export the Margin line	112
Export the case list as CSV	113
Download the files from the Medit Link Web	114
Order Process and Order Management	117
Ordering	117
Checking the ordered case	119
Managing Delivery Date	120
Job Completed – Work Box, In Box	122
Job Complete – Case Box, Order Box	124

Introduction and Overview

Medit Link is a collaborative tool designed to facilitate communication between clinics and labs to support patient-centered workflows. It includes tools according to your business type (Clinic or Lab) and specific features can be used according to the account type. Medit Link helps users manage a series of digital dentistry workflows including scanning with Medit's intraoral and tabletop scanners, prosthesis order and manufacturing process for prosthesis through a user-friendly system.

Medit Link consists of a web-service and a computer application.

Medit Link Overview

		Web	App
		<ul style="list-style-type: none"> • Register • Search • Search and manage partners • Share cases via Web Viewer with a sharable link • Access cases made from the App. 	<ul style="list-style-type: none"> • Enter form information • Scan • Place/Receive orders and make payments • Link with CAD • Manage patients
Clinic	Admin	Search and manage partners and your organization, manage payment information, and share cases. <ul style="list-style-type: none"> • Dashboard • Case Box • Case Talk • Partners • Membership • Payment • My Info 	Enter form information, scan, link with CAD, place orders, make payments, and manage patients. <ul style="list-style-type: none"> • Dashboard • Case Box • Order Box • Case Talk • App Box • Trash Box • Patient Management
	Member	<ul style="list-style-type: none"> • Case Box • Case Talk • My Info 	<ul style="list-style-type: none"> • Case Box • Case Talk • App Box • Trash Box

Lab	Admin	<p>Search and manage partners and your organization, manage payment information, and share cases.</p> <ul style="list-style-type: none"> • Dashboard • Work Box • Case Talk • Partners • Membership • Payment • My Info 	<p>Enter form information, scan, link with CAD, receive/ manage orders.</p> <ul style="list-style-type: none"> • Dashboard • In Box • Work Box • Case Talk • App Box • Trash Box
	Member	<ul style="list-style-type: none"> • Work Box • Case Talk • My Info 	<ul style="list-style-type: none"> • Work Box • Case Talk • App Box • Trash Box



The following service is not available on Medit Link China.

- Payment

Cloud Service

Medit Link supports the automatic upload/download of databases and files through AWS (Amazon Web Services). The data is automatically synchronized to the computers of multiple users in the organization so that they can access the same data from any computer at any time. Cloud-based backup and recovery keeps all your data safe. It keeps your data safe even in case of data loss due to sudden computer failure.

The Cloud is also useful for dentists and labs to quickly and easily share data, creating a working environment that minimizes the hassle of sharing data using email. All data with patient information is securely encrypted, stored, and transmitted, and complies with medical law, HIPAA regulations, and GDPR.

Workflow Management Using Ordering System

Medit Link provides order and payment system.

You can quickly and easily share your data by signing a contract between the clinic and the lab.

You can view/manage the status of your orders.

Credit card-backed system works with payment modules of online payment platforms, such as Stripe (www.stripe.com) and I'mPort (<http://www.iamport.kr>), making payment services easy and convenient.




The following service is not available on Medit Link China.

- Payment

Utilizing the Dashboard

You can view the statistics for the entire work done on Medit Link in visualized graphs and charts – both in the Medit Link App and Web.



Installation

System Requirements

System Requirements for Windows

	Laptop	Desktop
CPU	Intel Core i5 2.6 GHz or higher	
RAM	16 GB or higher	
Graphics	NVIDIA GeForce GTX 760 (2 GB) or higher / or equivalent AMD video card	
OS	Windows 8 64 Bit (unavailable in 32 Bit) or higher	

System Requirements for macOS

	Mac
Processor	M1 with 8-core, 7-core GPU, 16-core Neural Engine
RAM	16 GB or higher
OS	Monterey 12


Installation on Windows

Medit Link is installed as a package with Medit Scan for Clinics and Medit Scan for Labs.

- ① Run Medit_Link_ChinaSetup.exe.
- ② Select the language for installation.




- ③ Select the folder for installation of Medit Link. Read “License terms and conditions” carefully and check the box “I agree to the License terms & Conditions.” Click “Install” to proceed to the installation.




- ④ Make sure that the scanner cable is disconnected from the computer before proceeding to the next step.



The installation may not proceed if the scanner is connected to your PC.



- ⑤ From this step, the installation will continue automatically and may take a few minutes. Do not turn off the computer until the installation is complete.



- ⑥ After the installation is complete, restart the computer to use Medit Scan Application with Medit Link. You can skip this step If you do not have a scanner.




The scanners may not work without restarting the computer.


Installation on macOS

Medit Link is installed as a package with Medit Scan for Clinics.


- ① Run Medit_Link_x.x.x.pkg and click “Continue.”




- ② Read the End User License Agreement and click “Continue.”



- ③ Select the disk to install and click “Continue.”




- ④ Click the “Change Install Location...” to change the installation destination.




- ⑤ Click “Install” and type your Mac password to allow installation if requested.



- ⑥ Wait until the installation completes.




- ⑦ Click “Close” to exit the installation window.



Registration

Landing Page

Register for Medit Link at www.meditlink.com.




Register for Medit Link China at www.meditlink.cn.

Registration


- Medit Link requires an initial registration as Clinic (administrator & member) or Lab (administrator & member).
- Administrator is usually a representative of your organization and has unlimited access to all the services, whereas members have a limited access to Medit Link services.

- ① Click the “Don’t have an account? Sign up.” on the login page.



The image shows a login form for MEDIT Link. At the top, the text "MEDIT Link" is displayed in blue. Below this, there are two input fields: "ID" with a placeholder "Email Address" and "Password" with a placeholder "Password" and an eye icon to toggle visibility. Under the password field, there is a checkbox labeled "Keep me logged in." and a link "Need help?" with a dropdown arrow. A large blue button labeled "Log In" is positioned below these fields. Underneath the button is a link that says "Don't have an account? Sign up." At the bottom of the form is a language selection dropdown menu showing "English (United States)" with a small American flag icon and a downward arrow.

- ② Click the “Register” button of the desired account type among the four types below.



The image displays two columns of registration options. The left column is titled "Clinic" and features an illustration of a doctor with a tablet. The right column is titled "Lab" and features an illustration of a lab technician with a tablet. Each column contains four registration cards. The top two cards are for "Administrator" roles, and the bottom two are for "Member" roles. Each card lists available features and includes a blue "Register" button.

Category	Role	Available Functions	
Clinic	Administrator Dentist (Manager)	Order/Payment, Member Management	
	Member Dental Hygienist and Dental Staff	Oral Scanning	
	Lab	Administrator Representative/Leader (Manager)	Price Management for Dental Prostheses, Order Management, Member Management
		Member Laboratory Employees	Receiving Data for Making Dental Prostheses, Sharing work files and production processes.

- ③ After specifying the name of your clinic or lab, conduct the following action according to your account type.
 - If you are registering as a clinic or lab, click “Check Availability” to see if the clinic name or lab name is available.
 - If you are registering as a member, click Search.” Your organization should be activated prior to your registration.
- ④ Once you complete all the fields and click “Register,” a verification email will be sent to you.
- ⑤ To complete the registration process, verify your email address with the link in the verification.

Login and Welcome Page

- ① Enter your email address and password and click “Login” on the login page.

MEDIT Link

ID

Email Address

Password

Password 

Keep me logged in. [Need help?](#) ▾

Log In

[Don't have an account? Sign up.](#)

 English (United States) ▾

- ② The following dialog appears if you have not filled organization information.



To use Medit Link services, information about your organization should be entered using the administrator account.

Welcome to Medit Link

Grow your organization through close collaboration with clinics or labs.

Address Information Organization Information Additional Information Scanner Information

[Here are four simple steps to get you started with your personal order system.](#)

You will be able to access all of Medit Link's services after you have completed filling out your information.

Step 1
Fill in your organization information.

Step 2
Select a partner clinic or lab you want to work with.

Step 3
Create a partnership and contract with your selected partner.

Step 4
Download Medit Link on your computer and start sending and receiving cases.

[Set Up Later](#) [Set Up Now](#)

- ③ Click the “Set Up Now” and enter the Address Information, Organization Information, Additional Information and Scanner Information step by step.




You can also register and manage intraoral and tabletop scanners and milling machines.



If you agreed to receive marketing information during registration, you will receive useful information for the product and service from the manufacturer.

- ④ The information entered on the **Welcome Page** (My Organization etc.) can be edited anytime on “My Organization” page.

Medit Link Web



System Requirements for Medit Link Web


Recommended Browsers

PC	Chrome, Firefox, Safari, Edge (Chromium), 360 Safe
Android	Chrome, Firefox
iOS	Safari, Chrome, Firefox

Minimum Specifications for Mobile Devices

Device	Android Same or above the specs of Galaxy S7	iOS Same or above the specs of iPhone 7
OS	Android 7.1 Nougat	iOS 11.0
RAM	3 GB or higher	2 GB or higher
AP	Above Snapdragon 710	A10 Fusion

Dashboard



Graphs and charts in the dashboard show statistics for the entire work done in Medit Link.

Scanner	Provide information about the active scanners, including scanner usage time, calibration interval, and calibration elapsed time.
Cloud Storage	Display the total capacity of the cloud and the amount in use.
Statistics	Provide statistics on cases based on job status.
Order Status	Provide statistics on cases based on order status.
Graph	Graph job information by period.
Tooth Type	Provide statistics on the number of scans and orders by treatment information.


Medit Insights

Medit provides the reporting service and helps you analyze your usage statistics. You can receive a report email on the second day of every month. No additional sign-up process is required. Just enable “Subscribe Medit Insights.”



The following service is not available on Medit Link China.

- Medit insights



If you don't use the email address you used for the Medit Link account, you can set another email address to receive Medit Insights.

Settings

Receive the report via

Medit Link account


Another email address

Enter another email address to receive the report email.

i The report will be sent on the 2nd day of every month.

Cancel **OK**

You will receive the statistics email like below on the second day of every month.




- **Banner**
Medit provides useful tips and information here. If you don't want to see the banner, uncheck the "Marketing and Promotion Agreement" check box in My Info > My Account.
- **Statistics Information**
The statistics information is displayed, such as scanning time, created cases, most ordered lab, etc.
- **Medit Link Download**
The button provides the download link of the Medit Link.
- **Dashboard**
All statistics information is based on Dashboard. You can see more details in Dashboard.

How to unsubscribe the email

- Disable "Subscribe Medit Insights" in Dashboard.
- Click "Unsubscribe" in the received email.


Case Box/Work Box

It shows all cases which the user has created on the Medit Link application. See files via a web browser without installing a separate program.



Filtering by status

You can filter cases based on their status.



- The user can click the status they want to view. The browser will display only cases in the chosen status.
- The user can select multiple statuses at once.

Case Box – Clinic	All/Form Only/Need Processing/Scan Completed/CAD/Ordered/Case Completed
Work Box – Lab	All/Form Only/Need Processing/Scan Completed/CAD/Done


Filtering using period settings

You can filter and display cases that were last modified within the period you set.

- ① Click the duration filter.



- ② Select the range.




- All: Shows all cases.
- Custom Range: Specify start and end dates to display cases within the period.
- Today: Show cases last modified today.
- Yesterday: Show cases last modified with yesterday's date.
- Past 7 days: Show cases last modified during the last 7 days, excluding today.
- Past 30 days: Show cases last modified during the last 30 days, excluding today.
- This Month: Show this month's cases.
- Previous Month: Show previous month's cases.

The case list will be updated in real-time with the cases created or last modified within the selected duration.

Filtering with search capabilities

You can leverage search capabilities to search for the required information easily and quickly.



A search bar interface with a dropdown menu on the left showing 'Case Name' with a downward arrow, a vertical separator line, the text 'Case Name / Contains', and a magnifying glass search icon on the right.

Each box provides the following search options:

Case Box – Clinic	Case Name, Patient Name, Patient ID, Lab Name
Work Box – Lab	Case Name, Patient Name, Clinic Name, Order ID

Case Talk

Cases created in the Medit Link can be shared with restrictions, without restrictions, with specific people, and registered members of Medit Link. The user can see all shared lists through “Case Talk” in the Medit Link application and Medit Link web service.

[App] Create a case in Medit Link App

[App] Create a form and scan

[App/Web] Share the case

[App/Web] Access the shared case in Case Talk menu




Shared cases can be viewed via a web browser without installing a separate program, enabling easy communication over the 3D data on various devices.


Creating a Case Talk

You can create a shared link from a case created in Case Box/Work Box.

- ① Select the case you want to share.



- ② Click the “Share” button.
- ③ Select the files you want to share.



- ④ Set the sharing permission type and sharing options.
- ⑤ Click “Get Link” to create a shared link.

Set the sharing permission type and sharing options

Share through Case Talk

Share With:

Share via:

Expiration Date:

PIN Code: PIN Code

Share the patient's name

Grant Download Permission **Premium Feature**

Allow comments

Comment:

0/1024

By pressing the "Get Link" button, you agree that you understand personal information can be collected and provided to a third party (or shared in any public domain) ONLY with the consent of the person whose images on which you have worked. You also agree the consent of that person was truthfully obtained, and Medit is not responsible for any problems which may arise there after.

[Close](#) [Get Link](#)

Share With

- Share with anyone
 - Creates a shared link. You can share the files with anyone.
- Share with specific people
 - Creates a shared link so you can send the link to specific people.

Share via

If you select “Share with specific people,” you can share the link to specific people with various way.

- Email Address
- SMS
- Partners
- NFT



The following service is not available on Medit Link China.

- NFT
- SMS

Expiration Date

- Set the expired date of the link. After the link expired, the people you shared the link with will no longer be able to access it. You can extend the expiration period up to 90 days from today.

 If you select NFT, this option is unavailable to control. The expiration date will be set as unlimited.

PIN Code

- Set a 6-digit access code to ensure data security.

 If you select NFT, you should set the PIN code to protect your personal information.

Share the patient's name

- If you check this option, all people who are connected to link will see the patient's name.

Grant Download Permission

- If you check this option, all users who have access to the shared link can download the files.

 If you select NFT, this option is unavailable to control. The option is always on.

Allow comments

- If you check this option, you will be able to communicate through the shared link.
- If you uncheck the option, the recipient can only view your comments, but will not be able to add the ones of their own.


 If you select NFT, this option is unavailable to control. The option is always on.

Comment

You can leave a message as a comment.

Case Talk Details

You can see files in the shared case. In the case talk details page, the below functions are provided.



	Download	Download all the files in this case to the local PC.
	Case Information	Show the detailed information of the case.
	Share	Modify the sharing setting.
	Copy link	Copy the shared link to the clipboard. The copied link shows the shared files on your web browser.
	Live	Share your screen with the connected users.
	Attendees	Show all users who are accessing the shared link.
	Delete	Delete the shared link. All shared items will be deleted and the people who try to access to the link will no longer be able to access it.




The following services are not available on Medit Link China.

- Live
- Attendees

Web Viewer

You can view shared data in Case Talk in a web browser.

It is possible to communicate with the people who have access to the link using annotations and comments.




	Download	Download all the files in this case to the PC.
	Case Information	Show detailed information of the case.
	Share	Modify share setting.
	Live	Share your screen with the connected users. The screen sharing is synchronized with the connected users.
	Attendees	Show all users who are connected to the shared link.
	Tooth Labeling	Set the way of tooth labeling

 The following services are not available on Medit Link China.

- Live
- Attendees

Annotation


You can add annotations to a specified location on the 3D model. The annotations you add can be viewed by anyone with the link.




The following service is not available on Medit Link China.

- Annotation

① Click on “New Annotation” on an any specific spot.



② Type what you want to deliver, then click “OK” to create an annotation.



③ You can use annotations as a simple communication channel.

Comments

Communicate with users who have access to the link.



The following service is not available on Medit Link China.

- Comments

Real-time screen sharing


Share your screen in real time.



The following service is not available on Medit Link China.

- Live



- ① Click the “LIVE”  to share your Web Viewer screen with the connected users.



- ② Click “OK” to start sharing your screen.
- ③ Other people who have accessed the link accept screen sharing.
- ④ The screen of the user who requested screen sharing is synchronized with the screen of the person who received the request.

Data Tree

You can group and display data according to attributes, show/hide data, adjust opacity, and so on.

Show/Hide	Click on the  icon to adjust the display status of your data.
Opacity	<p>When you mouse over an element, a slider bar is displayed to adjust the opacity.</p> <p>Use the slider bar to adjust the opacity of your data.</p> 
Extended menu	<p>The extension menu is provided on the right side of an element.</p>  <ul style="list-style-type: none"> • Show This Only: Hide all other data except for current data. • Hide: Hide the current data.

3D Model Control

	Pan	Move the model.
	Rotate	Rotate the model.
	Zoom in/out	Zoom in/out the model.
	Zoom Fit	Fit the model to the screen.
	Model Display Mode	Change the color display mode for the model.

Partners


You can search for potential partners to request a contract and see your current partner list.

 Services using Google Maps are not available on Medit Link China.

Partner Search


You can search for Medit Link users based on their business name and/or address.

In the search results, you can send a partnership request to an organization you would like to work with.



Featured Partners

Medit Link shows the potential partners close to you based on your address information.



My Partner List

My Partner List displays a list of the companies you currently partner with or have a contract with.



On Partnership vs. Contract

- Once you partner with a business, you can request to make a contract with the organization.
- If your partner accepts your contract request, you can proceed with payment for the placed orders.

The screenshot shows the 'My Partner List' page in the MEDIT Link system. The page has a top navigation bar with 'Search for Partners', 'Featured Partners', and 'My Partner List'. Below this, there are filter tabs: 'All 3', 'Partnership 2', 'Contract 1', and 'Pending 0'. The main content area is divided into two sections, each showing a partner card for 'Wonhoonchoi-Lab' and 'Wonhoonchoi_Lab3'. Each card includes a profile picture, name, location (Seoul, Korea, Republic of), and a 'Contract Details' button. To the right of the partner list is a map showing the locations of these partners in Seoul, Korea. The map includes various landmarks and streets, with markers indicating the locations of the partners.




The following service is not available on Medit Link China.

- Contract

Partner Invitation

You can proceed with a partnership request with a subscription request to a user who is not a member of Medit Link.

- ① Click “Invite Partner” to display the following screen:





Invite Partner

Invite your desired partner to join Medit Link. When the lab accepts your invitation, a partnership is automatically formed.

Cancel

Send Invitation

- ② Enter the email address of the vendor you want to request a partnership from and send the invitation mail to proceed with the invitation.
- ③ If the invited company confirms the membership and proceeds with the membership, the partnership will be established, and the Medit Link service based on this will be available.

Membership

Medit Link’s membership policy aims to establish a more extensive digital environment, gradually expanding to provide various services based on the cloud platform. We aim to provide our customers a more consistent and personalized experience with our new membership policy.

The available membership plans are

- Standard
 - Free plan for every member
 - 1GB cloud storage
 - Efficient for new users
- Premium
 - Paid plan
 - 10 TB of cloud storage
 - Efficient for users with many cases and various data
 - One-month free trial included


	Standard	Premium
Cloud Storage*	1GB	10TB
Price (Monthly)	Free	\$0.99/Month

*Raw data will be excluded from the cloud storage.




The membership plans will be available on Medit Link China. You can get the Premium plan for free until the official release of the feature.

Membership




My Plan

You can check the membership plan you are using now. You can also subscribe to or unsubscribe from the membership when you want.



Cloud Storage

You can manage your cloud storage.




<p>Set Period</p>	<p>Set the period. You can select the files to remove in the period you set.</p> <ul style="list-style-type: none"> • All: See all files. • 1 Week Ago: See the files that have passed a week after creating. • 1 Month Ago: See the files that have passed a month after creating. • 6 Months Ago: See the files that have passed six months after creating. • 1 Year Ago: See the files that have passed a year after creating.
<p>Raw Data</p>	<p>The raw data does not take up the cloud storage, so it will not be in the pi-chart. You can manage the raw data as well for a specific purpose, such as data security.</p>
<p>Result & Attachment</p>	<p>The result and attachment files take up the cloud storage. If the cloud storage is full, you can secure the storage space by removing unnecessary files.</p>
<p>Delete</p>	<p>Remove selected files from the cloud storage to free up the space.</p>

Subscribe Membership

If you are using the Standard plan or using the Premium plan with a redeem code, you can subscribe to the Premium membership.


 The payment system is based on Stripe. (www.stripe.com)

① Click the “Subscribe” button.



The screenshot shows the Medit Link membership page. The left sidebar contains navigation options: Dashboard, Case Box, Case Task, Partners, Membership, Payment, and My info. The main content area is titled "Membership" and includes "Billing Info" and "Redeem" links. Under "My Plan", the "Standard" plan is displayed as "Free". It lists benefits: "Free plan for every member", "1 GB of cloud storage", and "Efficient for new users". A "Subscribe" button is visible at the bottom right of this section. To the right, the "Cloud Storage" section shows a donut chart indicating "54.68 MB / 1 GB" used. Below the chart, there are checkboxes for "Total Selected" (0.00 KB), "Raw Data" (117.51 GB), and "Result Files & Attachments" (54.68 MB). The footer contains links for Home, Medit Notice, Terms & Conditions, Privacy Policy, and Copyright © Medit Corp.

② Check the description for the plan and click “Subscribe.”




The screenshot shows a modal window titled "Subscribe Membership". It features an illustration of a person using a smartphone connected to a cloud. The "Premium" plan is highlighted, with a price of "0.99\$ USD / Monthly" and a note: "1 Month Free - for the first subscription!". A list of benefits includes: "A subscription plan", "10TB cloud storage", "Efficient for users with many cases and various data", and "One-month free trial after subscription". A note at the bottom states "* VAT not included". At the bottom right, there are "Cancel" and "Subscribe" buttons, with the "Subscribe" button highlighted by a red box.

- ③ Check the price and conditions and input necessary information for the payment. Then, click “Subscribe.”

The screenshot shows a payment interface for a 'Premium membership for Clinic'. On the left, the price is listed as \$0.00, with a note 'Then \$0.99 per day'. A table shows the breakdown: Premium membership for Clinic (\$0.99), Subtotal (\$0.99), and a 1-month free trial discount (-\$0.99), resulting in a total due today of \$0.00. On the right, the 'Pay with card' section includes an email field, card information (number, expiry, CVC), name on card, and country (South Korea). A 'Subscribe' button is highlighted with a red border. Below the button, a disclaimer states: 'By confirming your subscription, you allow MEDIT SG PTE. LTD. to charge your card for this payment and future payments in accordance with their terms.'


- ④ The payment will take several seconds. When it's done, you will see the below message.




Unsubscribe Membership

You can unsubscribe from the membership whenever you want. After unsubscribing, you can still use the Premium membership for the remaining days. The plan will be changed to the Standard plan from the next payment day, and the cloud storage will be limited to 1 GB. If the cloud storage capacity is exceeded, the files will not be backed up, and the order and share functions will not be available.


- ① Click “Unsubscribe.”



- ② Click “Cancel Subscription.”




- ③ Click “Cancel Subscription.”




- ④ The subscription has been canceled.

The plan will be changed to the Standard plan from the next payment day, and the cloud storage will be limited to 1 GB. If the cloud storage capacity is exceeded, the files will not be backed up, and the order and share functions will not be available.



Billing Info

You can check the payment history for membership subscriptions, including an invoice for each payment. You can also add or change the payment method whenever you want. The new payment method will be applied from the next payment.



The screenshot shows the "MEDIT Link" interface. The top navigation bar includes "Membership", "Billing Info", and "Redeem". The "Billing Info" section is active, displaying "My Plan and Billing Information".

My Plan and Billing Information

Plan	Premium
Amount	0.99\$ USD / Monthly
Next Payment Date	2021-09-09
Payment Method	Visa *****4242 (11/21) Stripe Edit Payment Method

If you haven't subscribed yet?
Experience a differentiated premium.


Billing History

Payment Date	Membership	From	To	Payment Amount	Status
2021-09-09	Premium	2021-09-09	2021-09-10	0.99 \$ (USD)	Payment Completed
2021-09-09	Premium	2021-09-09	2021-09-09	0 \$ (USD)	Payment Completed


The footer contains links for Home, Medit Notice, Terms & Conditions, Privacy Policy, and Copyright © Medit Corp. All rights reserved.

Add Payment Method

- ① Click “Edit Payment Method.”



- ② Click “Add payment method.”



- ③ Input the required information and click “Add.”

The screenshot shows a form titled "Add payment method". Under the heading "Card information", there is a text input field for "Card number" and a smaller input field for "MM / YY". Below this is a checked checkbox labeled "Use as default payment method". There are two buttons: a blue "Add" button and a white "Go back" button. At the bottom, a small text line reads: "By adding a payment method, you agree to MEDIT SG PTE. LTD.'s Terms of Service and Privacy Policy."

- ④ The newly added payment method will be set as default and used from the next payment date.

The screenshot shows a "Billing" page. Under "CURRENT PLANS", there are two identical entries for "Premium membership for Clinic" at "\$0.99 per day" with a renewal date of "September 10, 2021". Under "PAYMENT METHODS", there are two listed methods: a VISA card ending in "5556" marked as "Default" with an expiration date of "Expires 11/2021", and another VISA card ending in "4242" with an expiration date of "Expires 11/2021". A red box highlights the first payment method. At the bottom, there is a "+ Add payment method" link.

Change Payment Method

- ① Click “xxx.”

Billing

CURRENT PLANS

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

PAYMENT METHODS

 **** 5556 Default	Expires 11/2021	<input type="checkbox"/>
 **** 4242	Expires 11/2021	<input type="checkbox"/>

+ Add payment method

- ② Click “Make default” to change the default payment method.

Billing

CURRENT PLANS

Premium membership for Clinic

\$0.99 per day
Your plan renews on September 10, 2021.

Premium membership for Clinic

\$0.99 per day
Your plan renews on September 10, 2021.

PAYMENT METHODS

 **** 5556 Default	Expires 11/2021	<input type="checkbox"/>
 **** 4242	Expires 11/2021	<input type="checkbox"/>

+ Add payment method

Make default

Delete

BILLING INFORMATION

Redeem

Register the redeem code provided as a part of the promotion. You can use premium membership with redeem code.



The following service is not available on Medit Link China.

- Redeem

Register and Activate Redeem Code

- ① Input redeem code and click “Register.”

- ② Click “Register” or “Apply.”
- ③ Click “Register” or “Apply” on the pop-up dialog.

Register: Register the redeem code. You can activate and use the redeem code you registered.

Apply: Register and activate the redeem code. The promotion starts immediately.

Activate Redeem Code

You can have multiple redeem codes as “Registered” status. In this case, you can activate the code you want.

① Click “Activate.”

If there isn't any activated code, the one you choose will be started immediately. If there is an activated code, the other code you choose will be in “Waiting” status, and it will be activated automatically right after the current one expires.



Only one code can be in “Waiting” status.

Promotion	Code Number	Dealer Name	Registration Date	From	To	Status
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Yoshida	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate

Payment



The following service is not available on Medit Link China.

- Payment

Payment Details (Clinics, Labs accounts)

You can check all payment history and transactions through the Payment Detail page.

Payment account information (Lab account only)

This feature is only available to the administrator of Lab account.

- You must enter your payment account information so that your account can use the online payment system when ordering.



If you select a country that can support Stripe, you will be notified of the Stripe connection, and you can create a Stripe account by pressing the appropriate button.

My Information

My Account

You can view/modify the information you entered during the registration and check your membership.

The screenshot shows the 'My Account' page with a navigation bar containing 'My Account', 'My Organization', and 'My Price List'. The main content is divided into two columns. The left column, titled 'Account Information', features a profile picture placeholder with a plus sign and a note: 'JPG, GIF, or PNG only. Max size of 10MB'. Below this are fields for 'Email Address' (coldgreentea@naver.com) with a 'Change Email Address' button, 'Password' with a 'Change Password' button, 'Name' (WonhoonChoi_Lab_Stage), 'Language' (English(United States)), and 'Country' (United States). A note states: 'The language is applied to the email sent by Medit Link and the notification.' At the bottom of this column is an 'Organization Information' section with an 'Edit Information >' button. The right column, titled 'Group Information', shows a list of groups with a filter 'All 1'. The first group is 'WonhoonChoi_Lab_Stage' with a profile icon, email 'coldgreentea@naver.com', and 'Registration 2020-07-03'.

My Organization (Admin account only)

View/modify information of your organization.

The screenshot shows the 'My Organization' page with a navigation bar containing 'My Account', 'My Organization', and 'My Price List'. Below the navigation bar are four tabs: 'Address Information' (checked), 'Organization Information', 'Additional Information', and 'Scanner Information' (checked). The main content area is titled 'Address Information Mandatory *'. It contains several fields: 'Address1 *' with the value '23 Incheon-ro 22-gil, Seongbuk-gu, Seoul, Korea' and a 'View Location On Map' button; 'Address2' with the value 'Apartment, Suite, Unit, Building, Floor, etc.'; 'Zip/Postal Code *' (02835) and 'City *' (Seoul); 'Country *' (Korea, Republic of) and 'State/Province *' (Seoul). A note at the bottom states: 'Please select your country information carefully because it is important for completion of your transactions.' At the very bottom, there are labels for 'Office Phone *' and 'Mobile Phone *'.

My Price List (Lab: Admin account only)

Enter the products offered by your lab and the price information for each product.

The screenshot displays the 'My Price List' interface. At the top, there are navigation tabs: 'My Account', 'My Organization', and 'My Price List'. Below the navigation, there are filter buttons for various dental product categories: 'All (0/99)', 'Inlay (0/8)', 'Onlay (0/8)', 'Coping (0/10)', 'Denture Coping (0/5)', 'Crown (0/15)', 'Pontic (0/23)', 'Post & Core (0/3)', 'Veneer (0/8)', 'Custom Abutment (0/2)', 'Surgical Guide (0/1)', and 'Offset substructure (0/4)'. Below the filters, there is a table with columns: Name, Type, Method, and Material. The table lists several inlay products:

<input type="checkbox"/>	Name	Type	Method	Material
<input type="checkbox"/>	Inlay Zirconia	Inlay	-	Zirconia
<input type="checkbox"/>	Inlay Gold	Inlay	-	Gold
<input type="checkbox"/>	Inlay Press Ceramic	Inlay	-	Press Ceramic
<input type="checkbox"/>	Inlay Resin	Inlay	-	Resin
<input type="checkbox"/>	Inlay Hybrid Ceramic	Inlay	-	Hybrid Ceramic



You can make a transaction based on the product list and price information if you have a contract relationship with the client.



Medit Link does not support price-based transactions for partnerships.

Notice and Updates

Displays notifications about using and updating the Medit Link service.

Help Center

Once you press “Help Center” you get transferred to Medit’s Help Centre which contains information about the company’s products and their functions.

Download the App


Download the installer for the latest version of the program.

- You can install the program on another computer if you have it already installed after the purchase of the product.



Program updates are automatic when the Internet is connected.

Medit Link App




1. Title Bar
2. Main Functions
3. Detailed Functions
4. Case List

Title Bar

The title bar displays icons for network status, cloud synchronization, and notifications.

Windows Title Bar




macOS Title Bar



User Information

User information displays basic information such as profile image, organization name, email address, and membership name. It also provides logout option and home page shortcuts.



My Information

1 

2 Organization Name
Wonhoonchoi_Clinic_Dev

Email Address
wonhoon.choi@medit.com


Name
Wonhoonchoi_Clinic

Close Log Out Edit >

1	Displays the profile image.
2	Displays the organization name, email address, and membership name.
3	Logs out of the currently connected account.
4	Logs in with the currently accessed account information and access the homepage. Check, modify, and manage registered information

Notifications

Provides notifications related to contracts/orders and work processes.



1	Shows only the unread messages.
2	Quickly navigate to the specific case or webpage by selecting the corresponding notification.

Download Case

Total Count / Item Selected	Check the number of all/selected cases that need to be downloaded.
Download All	Download all cases to the PC.
Cancel All	Cancel all requests.
	Download selected cases to PC.
	Cancel downloading of the selected cases.
	Retry downloading of the selected cases.

Network Status

Medit Link enables you to use all the features of the service seamlessly when you are connected to a network.

The Title Bar displays the network status for the smooth use of the service.

If the network connection is not good, the following restrictions occur:

- Some functions, such as ordering a case, creating a case, sharing via Case Talk, and cloud synchronization, are restricted.
- The overall performance of the service can be reduced due to the impact of the network environment.

In this case, you can enter the Offline Mode to perform basic tasks, and then use the features that require internet once you are connected.




The network speed symbol indicates the response speed between the current area and AWS (Amazon Web Service). It does not indicate your network speed. Sometimes it may show the status as poor or bad, even though your network speed is fast enough.

If the other process uses the same network, the speed may show poor or bad, and the overall performance of Medit Link could be bad. In this case, please do not share the network with other processes. For example, the “Upload Raw Data” feature takes up a huge part in the network process. If the symbol shows poor or bad status, go to Settings > Upload Raw Data and turn “Upload Automatically” off.

Main Functions


Box Overview



The table below shows which Boxes are available for Clinic and Lab accounts.

		Case Box	Order Box	In Box	Work Box	Case Talk	Trash Box	App Box
Clinic Account	Admin	O	O	X	X	O	O	O
	Member	O	X	X	X	O	O (Restoring cases)	O
Lab Account	Admin	X	X	O	O	O	O	O
	Member	X	X	X	O	O	O (Restoring cases)	O

Case Box/Work Box



1	Display the case name.
2	Display the key features available: Scan, CAD, Order, and Apps
3	Display available tabs: Scan, CAD, and File Viewer.
4	Functions for managing cases and files: Delete Case, Share Case, Attach File, and Export File
5	Display the form for the selected tab, which can be used in a scan or CAD program.
6	Select the type/method/material/shade information for the form based on the CAD settings.
7	Display the type/method/material/shade information applied to the selected teeth.



Tabs in Case Box/Work Box


Scan	CAD	File Viewer
Enter information for basic scan tasks (form and scanning options). The information is passed to the connected scan program and used as the order information.	Set CAD options. The information you set up is passed to the CAD program.	View result files generated by scanning programs and CAD. Provides the ability to effectively display and manage 3D and 2D data, both acquired and attached.

The following functions are provided for Case Box/Work Box.

	Status	Display the case status.
	Go to Detail	Provide a shortcut for the same case in another box (Case Box ↔ Order Box; In Box ↔ Work Box).
	Case Information	Utilize this tool to save case information as a PDF or print-out.
	Memo	Show the notes for the case. You can add memos when ordering to communicate with your partner with detailed information about the case.
	Case History	Provide detailed information about the case history. <div data-bbox="730 1482 1252 1921" data-label="Complex-Block"> <p>History</p> <p>Last Modified Date 3/29/2022 1:57 PM</p> <p>Scanning Date 4/12/2021 2:06 PM</p> <p>Creation Date 3/22/2022 1:17 PM</p> <hr/> <p>Case Size 435MB</p> <p style="text-align: right;">Close</p> </div>

	Delete	<p>Delete a case.</p> <p>You cannot delete ordered cases or cases created offline.</p>
	Clone	<p>Clone a case on the Case Box.</p> <p>This feature is helpful for replicating treatment information for the same patient or making an additional order with the same information or files in the case.</p>
	Attach	<p>Attach files in various formats.</p>
	Export	<p>Export scan or CAD result files to the specified file format.</p> <p>You can use the exported files for future processes such as CAD, CAM, etc.</p>
	Share	<p>Create a link for sharing.</p> <p>View and check the data through a browser without installing additional programs.</p> <p>Share the link freely or to specific people exclusively. Set a PIN code to access the link through the Medit Link website.</p>
	Support Request	<p>Send a request for technical support to the Medit Technical Support Center.</p>

Order Box/In Box



1	Display the case name and patient name.
2	Provide functions for order management (accepting/cancelling; delivering/receiving).
3	Tabs for each of the tasks.
4	Show case status and shortcuts. <ul style="list-style-type: none"> Go to Detail: A shortcut to check the details of ordered cases. Click the button to go to the Case Box or Work Box case. Status: Displays the job status of the case.
5	Display order number, patient name, partner name, scanning Date, order date, delivery date, memo.
6	Display the product information and delivery information.




Tabs in Order Box/In Box

Scan	File Viewer
Check the detailed information related to the order, such as order number, patient name, account name, scanning date, order date, delivery date, and memo.	View files generated by scan/CAD operations. The file viewer displays and manages the acquired or attached 3D models and 2D images.

Trash Box


The Trash Box is a space for keeping deleted cases before they are entirely removed or recovered.

 The data deleted from the Trash Box will be permanently deleted from your PC and the server and cannot be recovered.




Settings

General: It provides features for overall system setup.



Save Folder	Specify a folder to store the results generated from Medit Link and connected programs.
Languages Settings	Select the program language. Improve your current language and update in real time when a new language is added.
Teeth Labeling	Select a labeling system.
Notifications	Set whether to receive notifications for program updates.
Upload Raw Data	Determines whether to automatically upload raw scan data to the cloud. It helps to back up the original data and perform seamless scanning when using Medit Link with the same account on another PC. The original data cloud backup only applies to data acquired by Medit Scan for Clinics.
Quick Scan	Choose if you would like to have the scan stages below to automatically appear in Medit Scan for Clinics without having to fill in the form information. <ul style="list-style-type: none"> • Pre-op for Maxilla • Pre-op for Mandible • Face Scan
Viewing Options	Advanced Rendering Displays the 3D model more vivid and realistic rendering. Adjust Color Texture Adjusts the brightness of the 3D model. When viewing data using other software, the resulting colors may be slightly different from Medit Link.
Medit User Experience Improvement Program	Join in the Medit User Experience Improvement Program. We use the collected information to improve our product and service. We do not collect any information that identifies individuals.

Storage Management: You can check and manage the local and cloud storage space.



Set Period

Set the period. You can select the files to remove in the period you set.

- All: See all files.
- 1 Week(s) Ago: See the files that have passed a week after creating.
- 1 Month(s) Ago: See the files that have passed a month after creating.
- 6 Month(s) Ago: See the files that have passed six months after creating.
- 1 Year(s) Ago: See the files that have passed a year after creating.

Local Storage

Manage the PC storage space.

Raw Data	Remove the selected raw data from the PC. If the raw data is already uploaded to the cloud, you can download and use it whenever you want.
Result & Attachment	Remove the result and attachment files such as mesh, image, video from the PC. If these files are already uploaded to the cloud, you can download and use them whenever you want.
Others	See the files used to operate the system and other programs. It shows the storage of your PC. You cannot remove these files directly in Medit Link.
Delete	Remove selected files from the PC.

Cloud Storage

Manage the cloud storage space.


Raw Data	The raw data does not take up the cloud storage, so it will not be in the pi-chart. You can manage the raw data as well for a specific purpose, such as data security.
Result & Attachment	The result and attachment files take up the cloud storage. If the cloud storage is full, you can secure the storage space by removing unnecessary files.
Buy Storage	If you use Standard membership, you can subscribe to the Premium plan to expand to more cloud storage.
Delete	Remove selected files from the cloud storage to free up the space.



The following service is not available on Medit Link China.

- Buy Storage

Set Connected Program: Set/manage programs that work with Medit Link.



You can connect the following programs:

	Clinic Account	Lab Account
Scanning S/W	Medit Scan for Clinics	Medit Scan for Clinics Medit Scan for Labs
CAD	exocad	exocad



The following services are not available for Medit Link for macOS.

- Medit Scan for Labs
- exocad integration

Export Option: Set how to export files from Medit Link.

Settings

General
Storage Management
Set Connected Program
Export Options
Case Converting Tool
Proxy Option

Export Options

CAD Type
Maestro

Export Files To
C:/Users/Medit-Wonhoonchoi/Desktop/test/

Mesh Format
Convert mesh files to the below formats
 meditMesh OBJ PLY STL
 Do not convert the attached files

Image Format
Convert images to the below format
 JPG JPEG PNG BMP

Compress & Combine
 Compress Files
 Combine Individual Mesh

Options when same folder or file already exists.
 Always ask when exporting
 Export a new file and overwrite the existing file. The file currently stored will be lost.
 Create a folder and export the files

Rule for Folder Name
 Creation Date Scanning Date Order Date
 Patient Name Case Name Partner Name
 Patient ID Order ID
 Folder Name
 Case Name


Rule for File Name
 Creation Date Scanning Date Order Date
 Patient Name Case Name Partner Name
 Patient ID Order ID
 File Name
 Case Name

Apply

Close


CAD Type	Select the CAD program you want to use for your design. The program will adjust the files based on this selection, allowing you to proceed with the design without any additional axis alignment.
Export Files To	Specify the path where you want to store the results.
Mesh Format	Specify the format of the final mesh file.
Image Format	Specifies the format of the exported image files.
Compress & Combine	Compress Files: Compresses all results into a single file. Combine Individual Mesh: Combines the result of maxilla or mandible with individual tooth. Useful when importing several files into a program is not possible.
Options when same folder or file already exists	Select an action when same folder or file name exists. <ul style="list-style-type: none"> • Always ask when exporting. • Export a new file and overwrite the existing files. • The file currently stored will be lost. • Create a folder and export the files.
Rule for Folder Name	Define the rules for the folder name. Creates a new folder with a defined name and saving the result files in it. Below terms is available for the rule. Clinic account Date Created, Scanning Date, Order Date, Patient Name, Case Name, Partner Name, Patient ID, Order ID Lab account Date Created, Scanning Date, Order Date, Patient Name, Case Name, Partner Name, Order ID
Rule for File Name	Define the rules for the output file name. Below terms is available for the rule. Clinic account Date Created, Scanning Date, Order Date, Patient Name, Case Name, Partner Name, Patient ID, Order ID Lab account Date Created, Scanning Date, Order Date, Patient Name, Case Name, Partner Name, Order ID

Case Converting Tool: Used to move cases between different accounts.




How to use

- ① Open Medit Link and press the “Settings” button located in the left bottom corner of the window.




- ② Choose the “Case Converting Tool” option from the left side menu as shown on the picture below:




- ③ Click on the  icon to select the file to import from your PC.

- ④ Press “Import” button.




- ⑤ The case will appear in your Case Box (Clinic account) or Work Box (Lab account).

Proxy Options: In some limited network environments, you can change network settings.



About: Displays the details of Medit Link.



User Guide: Gets the user guide for Medit Link. Requires an Internet connection.

Patient Management

Patient management provides functions for adding, editing, deleting, and merging patients.

Patient Management

Allows to select patients to delete or merge their information from different cases.
Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date	
<input checked="" type="checkbox"/>	5-9_JScan_CAMOrder_Teeth2				3	2020-09-18	2020-09-18	
<input checked="" type="checkbox"/>	5-9_JScan_CAMOrder_Teeth2				1	2020-09-18	2020-09-18	
<input checked="" type="checkbox"/>	test				2	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	test				22	2019-07-26	2020-09-17	
<input checked="" type="checkbox"/>	test				1	2020-09-17	2020-09-17	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16	
<input type="checkbox"/>	gttt				1	2020-09-15	2020-09-15	
<input type="checkbox"/>	Medit Link 2.3.0				2	2020-09-14	2020-09-14	
<input type="checkbox"/>	test				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	복사태스트				2	2020-09-13	2020-09-13	
<input type="checkbox"/>	test3				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	test1				1	2020-09-13	2020-09-13	
<input type="checkbox"/>	testq1231				26	2019-07-25	2020-09-11	

Close

Adding Patients

① Click “New Patient.”

New Patient

Patient Name*
Add new patient

General
 Female Male

Birth Date
yyyy mm dd

Patient ID

Phone Number

Email Address

Memo

Creating patient information involves consent from the patient.
Patient information can be shared with 3-parties on order.

[Terms & Conditions](#)

Cancel Save

② Fill in the required information in the New Patient dialog, then click “Save.”

Deleting Patients



You can delete a patient if there are no cases for the patient.

- 1 Select the patient to delete from the list.

Patient Management

Allows to select patients to delete or merge their information from different cases.
Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date		
<input type="checkbox"/>	test				0	2019-08-28	2019-09-02		
<input type="checkbox"/>	test				0	2019-08-28	2019-08-28		
<input type="checkbox"/>	test				0	2019-08-28	2019-08-28		
<input type="checkbox"/>	A.J Abutment2				0	2019-08-28	2019-08-28		
<input type="checkbox"/>	Scanbody				0	2019-08-28	2019-08-26		
<input type="checkbox"/>	Margin Line Creation				0	2019-08-28	2019-08-22		
<input type="checkbox"/>	Abutment Library				0	2019-08-28	2019-08-22		
<input type="checkbox"/>	GuideMessage				0	2019-08-28	2019-08-16		
<input checked="" type="checkbox"/>	margin				0	2019-08-28	2019-08-07		
<input checked="" type="checkbox"/>	test2				0	2019-08-28	2019-08-07		
<input checked="" type="checkbox"/>	salkjadsklf				0	2019-08-28	2019-08-07		
<input checked="" type="checkbox"/>	Test				0	2019-08-28	2019-08-07		
<input checked="" type="checkbox"/>	test				0	2019-08-28	2019-07-26		
<input checked="" type="checkbox"/>	test				0	2019-08-28	2019-07-25		

Close

- 2 Click "Delete."

Merging Patients

Manage all cases belonging to one patient all together by merging patient information.

- 1 Select the patient you want to merge from the list.

Patient Management

Allows to select patients to delete or merge their information from different cases.
Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date		
<input type="checkbox"/>	5-9 JScan_CAMOrder_teeth2				3	2020-09-18	2020-09-18		
<input type="checkbox"/>	5-9 JScan_CAMOrder_teeth2				1	2020-09-18	2020-09-18		
<input checked="" type="checkbox"/>	test				2	2020-09-17	2020-09-17		
<input checked="" type="checkbox"/>	test				22	2019-07-25	2020-09-17		
<input checked="" type="checkbox"/>	test				1	2020-09-17	2020-09-17		
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16		
<input checked="" type="checkbox"/>	ASD Demo				1	2020-09-16	2020-09-16		
<input checked="" type="checkbox"/>	gfff				1	2020-09-15	2020-09-15		
<input type="checkbox"/>	Medit Link 2.3.0			2020-9-24	2	2020-09-14	2020-09-14		
<input type="checkbox"/>	test				1	2020-09-13	2020-09-13		
<input type="checkbox"/>	测试测试				2	2020-09-13	2020-09-13		
<input type="checkbox"/>	test3				1	2020-09-13	2020-09-13		
<input type="checkbox"/>	test11				1	2020-09-13	2020-09-13		
<input type="checkbox"/>	testq1231				26	2019-07-25	2020-09-11		

Close

- ② Click “Merge Patients.”
- ③ Select the patient you would like to merge all cases to.
- ④ All cases are merged based on the selected information.

Merge all patient information to

Please select the patient you'd like to merge all data to.
After merging all data together, all information in other patient files will be deleted.

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date
<input checked="" type="radio"/>	5-9_jScan_CAMOrd...				1	2020-09-18	2020-09-18
<input type="radio"/>	test				2	2020-09-17	2020-09-17
<input type="radio"/>	test				22	2019-07-26	2020-09-17

[Cancel](#) [Apply](#)

- ⑤ Click “Apply” to proceed with the merge.

Notices and Updates

Displays program updates and critical alerts.

Go to Medit Link Web

Transfers you to Medit Link homepage under the account that you are currently logged in with.

Help Center

Goes to the Help Center page. You can get various information about using Medit Link.

File Viewer

File Viewer groups and shows various data, including the data acquired from programs linked to Medit Link and data attached by users.

It provides basic management functions such as showing/hiding data, adjusting opacity, and deleting data.

It supports previews of files in the following formats:

- 3D data: STL, OBJ, PLY
- 2D images: PNG, JPG, JPEG, BMP, GIF, TIF, TIFF
- Video files: MKV, MP4, AVI, WMV

Data Tree

You can group and display all data for the case in the data tree. You can also show/hide data, adjust transparency, and delete data from the data tree.

	Show/Hide	You can show or hide the data.
	Opacity	<p>The slide bar appears when you click one of the data on the tree. You can adjust the opacity of the shown data by moving the slider bar.</p> 
	More	<p>Clicking the More icon shows additional features available.</p>  <p>Show This Only: Hide all other data but display only the corresponding data.</p> <ul style="list-style-type: none"> • Hide: Hide the data. • Duplicate: Copy the data and add it to the attachment group. • Export: Export the selected files. • Rename: Change the name of the data. • Delete: Deletes the data.



You can rename data only when online.

Model Control

3D Model Control

	Pan	Move the model.
	Rotate	Rotate the model.
	Zoom In/Out	Zoom in/out the model.
	Zoom Fit	Fit the model to the screen.
	Model Display Mode	Switch the model display mode between the “Texture On” and “Texture Off.”
	Grid Settings (mm)	Show or hide the grid and control its position in relation to the model (overlay on/off).


2D Images Control

	Next/Previous Image	Move among images.
	Edit	Edit the image.
	Pinch/Paste Window	<p>You can detach the image window and adjust it to the desired position. You can resize the window.</p> <p>When the window is restored to its original location, it returns to its default size.</p>

Capture and Memo


Capture the current screen as a 2D image, draw on it, or add notes.

You can save it as a file with your notes and send it as a reference.



Orthodontic Capture


Capture the images in current viewing direction and top/bottom/left/right from the current viewing direction.



Sync View	<ul style="list-style-type: none"> • When the Sync View option is turned on, you can rotate or zoom in/out the center image. The other images will be adjusted simultaneously as fixed viewpoints as you move the center image. • When the Sync View option is turned off, you can pan, rotate, or zoom in/out all images individually.
Save Each Image Separately	Save each image as high-quality resolution.
Show Texture	Turn on/off the texture for the model

Reassign Data


Change the data for maxilla/mandible with the Reassign Data feature. You can use the changed data to create and utilize the result images of Ortho Capture.



- **Scan Data**
Scan data for the case are all shown.
- **Maxilla**
The data you selected for the maxilla are displayed.
- **Mandible**
The data you selected for the mandible are displayed.


Fill Holes

Fill holes in 3D mesh data.



Trimming

Trimming provides various tools to edit data.



The selected 3D mesh data can be edited using different tools located on the bottom of the window.

Toolbox

Selection

	<p>Smart Teeth Selection</p>	<p>Automatically select all teeth of the arch leaving out gingiva parts.</p> <div style="background-color: #e1f5fe; padding: 5px; border: 1px solid #ccc;">  This function is only available for the scan data that has been acquired by Medit Scan for Clinics with the “Use GPU” option on. </div>
	<p>Smart Single Tooth Selection</p>	<p>Automatically select the area of a single tooth leaving out gingiva parts. Click, press, and drag the mouse on the tooth.</p>
	<p>Polyline Selection</p>	<p>Select all entities within a polyline shape drawn on the screen.</p>
	<p>Brush Selection</p>	<p>Select all entities on a freehand-drawn path on the screen. Only the front face will be selected. The brush comes in three different sizes.</p>

After selecting an area


	<p>Shrink Selected Area</p>	<p>Reduce the selected area each time you press the button.</p>
	<p>Expand Selected Area</p>	<p>Expand the selected area each time you press the button.</p>
	<p>Invert Selected Area</p>	<p>Invert the selection.</p>

Deselection

	<p>Selection/ Deselection</p>	<p>When on, deselect the area using various tools.</p>
	<p>Clear All Selection</p>	<p>Clear all selected areas.</p>
	<p>Delete Selected Area</p>	<p>Delete the data from selected area.</p>


How to edit data

Using various tools to select the area you would like to edit.




“Polyline” selects all entities within a polyline shape drawn on the screen, while “Brush” selects only the front face data.

① Revert the selection by pressing the “Invert Selected Area” button.




– Reduce/expand the selected areas using the “Shrink Selected Area” / “Expand Selected Area” tools.




② Click the “Delete Selected Area” button.


Optimize Occlusion Alignment

Optimize the occlusion alignment. You can choose the target data and adjust the occlusion alignment.




Compensation Level



<p>Mid</p> <p>Compensation Level <input checked="" type="checkbox"/></p> <p>  </p>	
<p>High</p> <p>Compensation Level <input checked="" type="checkbox"/></p> <p>  </p>	
<p>Off</p> <p>Compensation Level <input type="checkbox"/></p> <p>  </p>	

Mandibular Movement

You can play the mandibular movements recorded in Medit Scan for Clinics to see the dynamic occlusal contacts between the upper and lower dentitions.



	Free Movement	Play the recording of free mandibular movement.
	Left Lateral Movement	Play the recording of the mandible moving to the left.
	Right Lateral Movement	Play the recording of the mandible moving to the right.
	Protrusive Movement	Play the recording of the mandible moving forward.

Manage the Multi Occlusion Data in File Viewer

You can get multiple results of occlusion data from Medit Scan Application.

Medit Link provides the function to manage the results through the “File Viewer.”


How to manage multiple results

When the user gets multiple results of occlusion from Medit Scan for Clinics, Medit Link provides the option to select a result.

The list of results is listed by name, and you can select one item. The selected model will be shown in the “Data Tree” and “Model Display Area.”



You can collapse the list and change the item using .



View Options

Change the View Options to customize the layouts of the case list in the Case Box/ Work Box/In Box/Order Box.

View Options

-  **Status Filters** 6 in use >
-  **Table Items** 6 in use >
-  **View Type** List View >
-  **Grouping** None >
-  **Sorting** Last Modified Date >

Status Filters

Show or hide each filter representing the status of the case. The hidden status filters are not displayed on the screen in the Case Box/Work Box/In Box/Order Box. You can easily find the case you need with the shown filters.

View Options

< **Status Filters** 


Shown [Hide All](#)

- Form Only 
- Need Processing 
- Scan Completed 

Hidden [Show All](#)

- CAD 
- Ordered 
- Case Completed 

[Cancel](#) [Apply](#)



- ① You can click the status you want to view to display only cases in that state.
- ② You can select multiple statuses at a time.

Available Filtering Options for Each Box:

Case Box	All/Form Only/Need Processing/Scan Completed/CAD/Ordered/Case Completed
Order Box	All/Pending/Canceled/Rejected/Accepted/Shipped/Case Completed
In Box	All/Pending/Canceled/Rejected/Accepted/Ready/Shipped/Case Completed
Work Box	All/Form Only/Need Processing/Scan Completed/CAD/Done
Trash Box	Not supported

Table Items

Select the items to display in the case list in the Case Box/Work Box/In Box/Order Box. You can change the order of table items or adjust the table column width in the list view to best optimize for the user environment.

View Options

< **Table Items**
↻

Shown Hide All

- ⋮ Status 👁
- ⋮ Case Name 👁
- ⋮ Patient Name 👁
- ⋮ Form Information 👁
- ⋮ Last Modified Date 👁
- ⋮ Scanning Date 👁

Hidden Show All



- ⋮ Birth Date 👁
- ⋮ Order Date 👁
- ⋮ Lab Name 👁

Cancel
Apply

Status	Case Name	Patient Name	Form Information	Last Modified Date	Scanning Date
Form Only	test's Case	test	-	2022-09-29 오전 8:45	
Ordered	test's Case	test	14-Implant Crown / 15-Implant Crown / 16-Implant Crown	2022-09-29 오전 8:13	2022-09-29 오전 7:56
Form Only	test's Case	test	-	2022-09-28 오후 5:56	
Form Only	test's Case	test	-	2022-09-28 오전 10:41	
Scan Completed	test's Case	test	24-Inlay / 25-Inlay / 26-Inlay	2022-09-28 오전 9:27	
Scan Completed	Multi Occlusion - Clone - Clone - Clone ...	test	24-Inlay / 25-Inlay / 26-Inlay	2022-09-27 오후 7:03	2021-05-14 오전 9:27
Scan Completed	Preop 사용한 케이스 - Clone	straumann	15-Inlay / 26-Pontic	2022-09-27 오후 5:09	2021-11-18 오후 6:53

View Type


Change the view type between the Gallery View and List View on the Case Box.



1	Case name
2	Patient name and birth date
3	Form information
4	Scanning Date
5	Case status

Grouping

Select a grouping option to classify and rearrange cases in the list by the selected option.



For each box, the following grouping options are available:

Box	Grouping
Case Box	None
	Patient Name
	Scanning Date
	Order Date
Order Box	None
	Scanning Date
	Order Date
	Requested Delivery Date
In Box	None
	Order Date
	Requested Delivery Date
Work Box	None
	Requested Delivery Date



The Grouping and Sorting options are not available for the Trash Box and Case Talk.

Sorting

Select a sorting option to change the order of cases in the list according to the selected option.

View Options

 Status Filters	3 in use >
 Table Items	6 in use >
 View Type	List View >
 Grouping	None >

 Sorting	Last Modified Date ▾
<input type="radio"/> Patient Name	
<input checked="" type="radio"/> Last Modified Date	
<input type="radio"/> Scanning Date	
<input type="radio"/> Order Date	
<input type="radio"/> Lab Name	

For each box, the following sorting options are available:

Box	Grouping	Sorting
Case Box	None	Patient Name / Last Modified Date / Scanning Date/ Order Date / Lab Name
	Patient Name	Last Modified Date / Scanning Date/ Order Date / Lab Name
	Scanning Date	Patient Name / Last Modified Date / Order Date / Lab Name
	Order Date	Patient Name / Last Modified Date / Scanning Date/ Lab Name
Order Box	None	Patient Name / Scanning Date / Order Date/ Requested Delivery Date / Lab Name
	Scanning Date	Patient Name / Order Date / Requested Delivery Date / Lab Name
	Order Date	Patient Name / Scanning Date / Requested Delivery Date / Lab Name
	Requested Delivery Date	Patient Name / Scanning Date / Order Date / Lab Name
In Box	None	Patient Name / Order Date / Requested Delivery Date / Clinic Name
	Order Date	Patient Name / Scanning Date / Requested Delivery Date / Clinic Name
	Requested Delivery Date	Patient Name/ Order Date / Clinic Name
Work Box	None	Patient Name / Last Modified Date / Requested Delivery Date / Clinic Name
	Requested Delivery Date	Patient Name / Last Modified Date / Clinic Name




The Grouping and Sorting options are not available for the Trash Box and Case Talk.

Filtering and Search


Filtering by Date

Filter and display cases last modified within the set duration.

- ① Click the duration filter.



- ② Select the duration from the list.



- All: Shows all cases.
- Custom Range: Specify the start and end dates to display cases last modified within the period.
- Today: Show the cases last modified today.
- Yesterday: Show the cases last modified yesterday.
- Past 7 days: Show the cases last modified during the previous 7 days excluding today.
- Past 30 days: Show the cases last modified during the previous 30 days excluding today.
- This Month: Show the cases last modified during this month.
- Previous Month: Show the cases modified during last month.

The case list will be updated in real-time with the cases created or last modified within the selected duration.

Advanced Search Capabilities

Leverage search capabilities to quickly find the information you want.

A search bar with a dropdown menu showing 'All' and a search logic dropdown showing 'Search / Contains'. A magnifying glass icon is on the right.

The following Searching Criteria are available for each box:

Case Box	All, Case Name, Patient Name, Patient ID, Form Information, Lab Name
Order Box	All, Case Name, Patient Name, Patient ID, Lab Name, Order ID
In Box	All, Case Name, Patient Name, Clinic Name, Order ID
Work Box	All, Case Name, Patient Name, Form Information, Clinic Name, Order ID
Case Talk	Case Talk Name (Contains), Patient Name (Equals)
Trash Box – Clinic	All, Case Name, Patient Name, Lab Name, Deleted By
Trash Box – Lab	All, Case Name, Patient Name, Clinic Name, Deleted By

You can select the Searching Criteria and select between the Contains and Equals logic.

A search dropdown menu. At the top, it shows 'All' and 'All / Contains'. Below is a section titled 'Searching Criteria' with radio buttons for: All (selected), Case Name, Patient Name, Patient ID, Form Information, and Lab Name. Below that is another section with radio buttons for: Contains (selected) and Equals.

Detailed Explanation of the Process

Managing Cases: Creating, Editing, Cloning, Deleting

Creating Cases

Creating cases in Case Box (Clinic)

- ① Click “New Case” button in the upper right corner.

New Case Information

Patient Name*


Case Name*

i Creating patient information involves consent from the patient.
Patient information can be shared with 3-parties on order.

[Terms & Conditions](#)

[Cancel](#)

- ② Enter basic information such as patient name.
 - Entering a name that is not in the patient list creates a new patient.
 - You can enter a name in the patient list to select from the list.
 - In case you are registering a case for an already existing patient, click .
 - You can select a new patient after registering or editing information for an existing patient.



- ③ Case name will be filled in automatically, however you can still change it.
- ④ Press “Register & Scan” to complete the case registration and proceed to scanning immediately.
 - Otherwise, press “Register” which is useful when registering multiple cases at the same time.
 - You will be redirected to the main page of the Case Box.

Creating cases in Work Box (Lab)

- ① Click “New Case” button in the upper right corner.

New Case Information

Clinic Name*

Patient Name*

Case Name*

Requested Delivery Date

2020-10-03

Cancel Register

- ② Enter or select the clinic name.
- ③ Enter the patient’s name and case name.
 - Case name will be filled in automatically, however you can still change it.
- ④ Enter the expected delivery date considering the completion date of the work.
- ⑤ After registering the case, go to the job information entry screen.

Editing cases

Editing cases in Case Box (Clinic)

- 1 Click on the  icon next to the case name.

Edit Case Information

Patient Name*



 'test' is selected.

Case Name*

- 2 Click on the  icon to change the patient's name. Note that it cannot be edited directly.; instead, you can find a patient name you would like to replace the current one with and can edit the patient information from the list.

Patient List

Select a patient in the list to set the patient of the case. You can add or manage the patient information in Medit Link.

All  [+ New Patient](#)

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Registered Date	Recent Visit Date	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input type="radio"/>	test				0	2020-08-28	2020-08-28	
<input checked="" type="radio"/>	test				13	2020-08-23	2020-08-28	
<input type="radio"/>	test				85	2020-08-23	2020-05-08	
<input type="radio"/>	Wonhoon Choi				1	2020-04-01	2020-04-01	
<input type="radio"/>	ttttt				0	2020-08-18		
<input checked="" type="radio"/>	13				0	2020-08-20		
<input type="radio"/>	14				0	2020-08-20		
<input type="radio"/>	15				0	2020-08-20		

[Close](#) [Select](#)

- ③ Check the patient's name you would like to replace the current one with.
- ④ To create a new patient, press the "New Patient" button.
- ⑤ Press  to edit the registered information about the patient.
- ⑥ Press "Apply" to finish the process.

Editing cases in Work Box (Lab)

- ① Click on the  icon next to the case name.

Edit Case Information

Clinic Name*

Patient Name*

Case Name*

Requested Delivery Date

[Cancel](#) [Apply](#)

- ② Press "Apply" to after the editing is done.




You cannot edit the information of the cases received through the order in In Box.

Cloning cases

You can use the Clone feature to clone a case to rescan with the same treatment information. You can duplicate the entire data to order a case that has already been ordered or choose other cloning options.

Cloning cases in Case Box (Clinic)

- ① Click on the  icon.
- ② You will be asked what kind of information you would like to clone to the new case.



- ③ Choose whether to replicate only treatment information or data all together.
- ④ When cloning is complete, a case with the name “[Case name] – Clone” will be created.


Cloning cases in Work Box (Lab)

- ① Click on the  icon.
- ② Duplicate all data to create a new case.
- ③ When cloning is complete, a case with the name “[Case name] – Clone” will be created.

Deleting cases (Lab and Clinic)

You can delete cases.

- ① Click on the  icon.
- ② Click “OK” to delete the case.



- ③ Deleted cases are moved to the Trash Box and can be recovered if necessary.

Inputting form information

You can enter treatment information based on the connected program of your selection.

The information you enter will be used as basic information for scanning, CAD, and ordering.

Inputting Form Information in Case Box

Entering dental information in the Scan tab

- ① Select whether you are going to input treatment information by teeth or by arch.
- ② Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the order list.
- ③ Click “Save” to save your registration.

Entering dental information in the CAD tab

- ① Select one of the following options that will define your scanning strategy: Teeth/
Flexible Multi-Die/Multi-Die.
- ② Select the tooth number/arch/die number.
- ③ Define the settings for detailed options and CAD operations.
- ④ Click “OK” to save the settings.
 - The selected information is automatically registered to the order list.
- ⑤ Click “Save” to save your registration.

Inputting Form Information in Work Box

Entering dental information on the Scan tab

- ① Select the scan program you want to use.
- ② Select one of the following options that will define your scanning strategy: Teeth/
Flexible Multi-Die/Multi-Die.
- ③ Methods may differ depending on the selected scan program.

	Medit Scan for Clinics	Medit Scan for Labs
Teeth	O	O
Arch	O	O
Flexible Multi-Die	X	O
Multi-Die	X	O


- ④ Select teeth, arch, or die numbers for the scanning process.
- ⑤ Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the list.
- ⑥ Click “Save” to save your registration.

Entering dental information in the CAD tab

- ① Select one of the following options that will define your scanning strategy: Teeth/lexible Multi-Die/Multi-Die.
- ② Select the tooth number/arch/die number.
- ③ Define the settings for detailed options and CAD operations.
- ④ Click “OK” to save the settings.
 - The selected information is automatically registered to the list.
- ⑤ Click “Save.”

Inputting Bridge Information

You can enter information for a bridge in case same information is applied to all selected teeth when selecting several teeth of the same type a bridge icon will appear by their side. You can turn the bridge on/off by clicking the icon.



Setting Scanning Options

Depending on the program you are working with, various scan options are provided.

The selected scan options are passed to each scan program.

Medit Scan for Labs

General options

	Adjacent teeth	Mark the prepared and adjacent teeth and pass this information to the scanning program.
	Antagonist	Automatically marks the antagonist and passes this information to the scanning program.
	Articulator	Displays options for selecting the articulator in the scanning program.
	Movement Marker	Allows to use the Jaw Motion scan strategy in the scanning program.

Options by teeth

	Scan body	Allows to set the strategy for scanning scan body in the scanning program.
	Pre-op Scan	Allows to set the strategy for scanning pre-op model in the scanning program.
	Wax-up Scan	Allows to set the strategy for scanning wax-ups in the scanning program.

Copying Form Information

You can copy the treatment information.

Copy dental information using context menu

- ① Select the tooth for which the tooth information is set.
- ② Right-click and select “Copy” from the context menu.
- ③ Select the target tooth.
- ④ Click  to complete the copy.

Copy dental information using shortcut keys

- ① Select the tooth for which the tooth information is set.
- ② Press and hold CTRL and select the target tooth.

Deleting treatment information

You can delete the dental information by clicking  in the order list.

No.	Type	Method	Material	Shade	
16	Inlay	-	Zirconia	A4	
25	Inlay	-	Gold	-	
26	Onlay	-	Gold	-	

Delete dental information using context menu

- ① Select the tooth for which the tooth information is set.
- ② Right-click to select “Delete” from the context menu.
- ③ Select the target tooth.
- ④ Click  to complete the delete.

Delete dental information using shortcut keys

- ① Select the tooth for which the tooth information is set.
- ② You can delete the selected information by pressing the keyboard “Del” button.

Adding Product Information to “Favorite”

You can quickly select a product by registering it as a favorite.

The list of registered favorites is managed as a template and can be shared through the export and import features.

Adding to Favorites

- ① Select the tooth number.
- ② Select Type – Method – Material – Shade.
 - Type – Method – Material – Shades can be registered as favorite even if not all of them are selected.
- ③ Click .
- ④ A combination of the selected parameters, excluding tooth number, is listed on your favorite search list.

Managing Favorite Templates

You can register, manage, and share different types of templates.


Add	Adds a new template.
Rename	Renames the current template.
Delete	Deletes the current template.
Export	Exports the current template to a file.
Import	Imports a template file.

Copying form information between Scan and CAD tabs

You can copy the treatment information entered between Scan tab and CAD tab.

From Scan tab to CAD tab

- ① If there is no information input in the CAD Tab, it automatically copies the information when you go from the Scan Tab to the CAD Tab.
- ② If there is information in the CAD tab, you will have to copy it manually.
 - Right-click on the CAD tab.
 - Select “Convert Form Information from Scan” from the context menu.




- The information in the Scan tab is copied to the CAD tab.

From CAD tab to Scan tab

To copy information, use the following methods:

- ① Right-click on the Scan tab.
- ② Select “Convert Form Information from CAD” from the context menu.



Note that only Type-Method information will be copied.

Scanning process

Based on the entered information, you proceed with scanning.

- ① Enter/save treatment information and select the necessary scan options.
- ② Click “Scan” button in the upper-right corner of the window.

CAD Process

Medit Link is integrated with exocad.

CAD	Registering treatment information	Confirmation of treatment information	Direct CAD Execution	CAD Results import back to Medit Link
exocad	○	○	○	○

- ① In Settings > Set Connected Program, you can link the CAD program of your choice from the list.
- ② Enter/save task information and design options in the CAD Tab.
- ③ Click “CAD” in the upper-right corner of the window.
- ④ The linked CAD program is run, and the data is sent automatically.
- ⑤ Save the results after performing CAD.



Medit Link for macOS does not support the function to link with CAD.

Get Result Files

Export scan result files

In addition to sharing results with partners through orders, Medit Link allows you to store your results in the desired file format using the export function.

Export files in case

Options **File List**

Export Files To

Folder Name

File Name

CAD Type

Mesh Format

Convert mesh files to the below formats

meditMesh OBJ PLY STL

Do not convert the attached files

Image Format

Convert images to the below format

JPG JPEG PNG BMP

Compress Files

Combine Individual Mesh

Export All Occlusal Relationship

Save Case Information as PDF [Preview](#)

[Close](#) [Export](#)


Export Files To	Specifies the path where you want to store the results.
Folder Name	<p>The system creates a new folder and stores all results in it.</p> <p>The default folder name is set by naming rule. But you can change the folder name directly.</p> <p>You can define the naming rule through Settings > Export Options > Rules for Folder Name</p>
File Name	<p>Sets the file name for the result. The default file name is set by naming rule. But you can change the file name directly.</p> <p>You can define the naming rule through Settings > Export Options > Rules for File Name</p>
CAD Type	<p>Select the CAD program to be used for your design program.</p> <p>Medit Link will adjust the files based on this selection, allowing to proceed with the design without any additional axis alignment.</p>
Mesh Format	Specifies the format of the mesh file.
Image Format	Specifies the format of the exported image files.
Compress Files	Compresses all results into a single file.
Combine Individual Mesh	Combines the result of maxilla or mandible with individual tooth. Useful when importing several files into a program is not possible.
Export All Occlusal Relationship	When the multiple occlusal results generated from the scan app, you can export all occlusal results at once.
Save Case Information as PDF	Save the case information as pdf with the result.

Export selected files


Medit Link provide the export function from the expand menu of file viewer. Use this, you can export the selected files easily.

How to export selected files

- ① Click  menu on the group or an entity in the file viewer.
- ② Click “Export.”



- ③ The export dialog shows and only the files selected that you want to export.



- ④ Click “Export.”


Export the Margin line

The Medit Scan for Clinics provide the function to create the margin line and Medit Link can transfer the margin line to the CAD program.

Export the margin Line as color

In general, Medit Link convert the margin line as color and export the result with maxilla or mandible together. You can create the margin line in their cad referencing the color information.

Example



 This feature is available for the format which include the color information such as meditMesh, OBJ, PLY. The margin line does not convert into STL.


 CAD type > 3shape does not support this feature.

Export the margin line as file

In case of exocad, Medit Link export the margin line to “~margin.xyz” which can be imported in exocad directly.

You can use this file in exocad without additional work.

Example



Export the case list as CSV

Medit Link provide the function to export the case list as CSV format.


You can export the case list with basic information and manage cases easily.

“Export to CSV” function is available in the following menu of Medit link application.

- In Box
- Work Box

How to export the case list

- Go to In Box or Work Box.
- Select the cases that you want to export.



- Click  then click “Export to CSV.”
- Define path and file name then save the result.

Download the files from the Medit Link Web

Medit Link provide the download function.

If you want to get the result files without install the Medit Link application, you can download the files from the Medit Link web service.




The data should be uploaded on the Medit Cloud to download the files from web.

The download function is available in the following menu of Medit Link web service

- Case Box
- Work Box
- Web Viewer

✓ Case Box/Work Box


You can find the  from the upper right.




✓ Web Viewer

When you share the case using Case Talk, you can set whether to allow downloading or not.

If you check “Grant Download Permission” and share the case, all users who connect to shared link can download the files.




In the Web Viewer, you can find the  from the upper right.



How to download

- ① Click .
- ② Select the mesh format from the download dialog.



- ③ Click .

Order Process and Order Management

You can easily put both acquired and attached data into an order to your partner.

Ordering

You can order cases to your partners from Case Box.

- ① Click “Order” button.
- ② Check the information required for your order.

Order

General Order List File List

Partner*

Wonhoonchoi_Lab_Dev

Wonhoonchoi_Lab_Dev Partner Details

seoul, KR
[02842] 35, Gaeunsa-gil, Seongbuk-gu, Seoul, Republic of Korea
1026088546

Case Name* Patient Name*

FaceScanTest's Case test

Share Patient Name

Requested Delivery Date Delivery Time

2022-04-06 11:45 AM AM 8:00 PM 8:00

Memo

Cancel OK

- ③ Select the partner
- ④ Check the product information from the Order List tab
 - You can decide whether to share patient name.

- ⑤ Check the files to order from the File List tab
- ⑥ Enter the requested delivery date.
 - The default value is set to 7 days after the order date.
- ⑦ You can memos to your partner to outline some specific details.
- ⑧ When the multiple occlusal results generated from the scan app, you can send all occlusal results at once by checking “Order All Occlusal Relationship.”
- ⑨ Proceed with your order.
 - **Skip payment:** Proceed with the order without using the payment system provided by Medit Link. Available when ordering from a partner with the “Contracted” status.
 - **Pay:** Utilize the payment system provided by Medit Link to proceed with the order. Available when ordering from a contracted partner.
 - **OK:** We proceed with the order without any conditions regarding the price of the product. Available when ordering from a partner with the “Partnership” status.




The following services are not available on Medit Link China.

- Skip payment
- Pay


Confirm Order – In Box

- ① Confirm Order
 - Accept the order upon receiving it.
 - Accepted orders will be forwarded to the Work Box.



Cancel/Reject Order – Order Box/In Box

- ① Cancelling Order – Order Box
 - You can cancel an order before your partner accepts it.
 - Canceled cases can be re-ordered.



- ② Reject order – In Box
 - You can reject a received order.
 - Rejected cases can be re-ordered.



Checking the ordered case

You can check the details of the case you ordered from Order Box in clinic account. You can also check the case you received from the In Box in lab account.

Status	Case Name	Patient Name	Order Date	Requested Delivery Date	Clinic Name	Order ID	Price
Pending	test's Case	test	2022-09-29 오전 8:13	2022-10-01 오전 10:30	Wonhoonchoi-Clinic	10413603	
Pending	test's Case - Clone	test	2022-09-22 오후 7:03	2022-09-24 오전 10:30	Wonhoonchoi-Clinic	10374886	
Pending	test's Case - Clone	test	2022-09-22 오후 7:01	2022-09-24 오전 10:30	Wonhoonchoi-Clinic	10374864	
Accepted	test's Case	test	2022-09-13 오후 2:43	2022-09-15 오전 10:30	Wonhoonchoi-Clinic	10314293	
Pending	test's Case	test	2022-09-05 오전 7:31	2022-09-07 오전 10:30	Wonhoonchoi-Clinic	10268220	
Case Completed	Worked with 2.3.6 - Clone - Clone - Clone - Clone	test	2021-07-01 오전 8:51	2021-07-08 오전 8:50	Wonhoonchoi-Clinic	8579990	500 KRW
Case Completed	test's Case - Clone - Clone - Clone - Clone - Clone	test	2021-08-31 오후 4:10	2021-09-07 오후 4:10	Wonhoonchoi-Clinic	8730155	2,000 KRW
Case Completed	test's Case memo	test	2022-04-08 오후 11:...	2022-04-15 오후 7:45	Wonhoonchoi-Clinic	9569817	
Pending	test's Case	test	2022-05-10 오후 9:21	2022-05-13 오후 8:00	Wonhoonchoi-Clinic	9710148	
Pending	test's Case	test	2022-05-10 오후 9:17	2022-05-13 오후 8:00	Wonhoonchoi-Clinic	9710120	
Pending	test's Case	test	2022-05-10 오후 9:03	2022-05-13 오후 8:00	Wonhoonchoi-Clinic	9710039	
Case Completed	test1's Case	***	2022-04-19 오후 5:15	2022-04-22 오후 8:00	Wonhoonchoi-Clinic	9608665	
Case Completed	test's Case	***	2022-04-21 오후 2:13	2022-04-24 오후 8:00	Wonhoonchoi-Clinic	9621238	

- ① Order information

The information that can be included in the order is as following: Order ID/ patient name/account name/scanning date/order date/delivery date/delivery due date/memo

You can check the details of the ordered product and verifying the delivery information.
- ② Order List

Order list displays the details of the ordered product.
- ③ Shipping

Enter/check the shipping information.

Managing Delivery Date

You can set the expected delivery date when placing an order, and this information will be forwarded to your partner.


- ① Set requested delivery date – Clinic Account.
 - When ordering from Case Box, you can enter the requested delivery date.

The screenshot shows the 'Order' form with the following details:


- Order** (Title)
- General** (Selected tab), Order List, File List
- Partner***: Wonhoonchoi_Lab_Dev
- Partner Details** (Link): Wonhoonchoi_Lab_Dev, seoul, KR, [02842] 35, Gaeunsa-gil, Seongbuk-gu, Seoul, Republic of Korea, 1026088546
- Case Name***: FaceScanTest's Case
- Patient Name***: test
- Share Patient Name
- Requested Delivery Date**: 2022-04-06
- Delivery Time**: 11:45 AM (Slider bar from AM 8:00 to PM 8:00)
- Memo**: (Empty text area)
- Buttons**: Cancel, OK

- The default value is 7 days after the date of the order.
- Enter the date/time you wish to receive the products and proceed with the order.
- You can set the time easily using slider bar.

- ② Set expected delivery date – Lab Account
 - You can enter the expected delivery date for the case ordered from In Box.



- ③ Check the requested delivery date with the partner and make adjustments using .
- ④ Enter the expected delivery date considering the actual working day. If required, you can leave a comment also.
- ⑤ Click “Apply.”



The expected delivery date and memo will be forwarded together with the order.


- You can check the date and memos in the Order Box.
- You can change the expected delivery date, and the changes will be forwarded to your partner.
 - You can check the case history using .

Job Completed – Work Box, In Box

Once all tasks are done, press “Done” to complete the case work.

✓ “Complete” a case from the case detail

You can find the  from the right top corner in the case. Click this menu, you can complete the work using “Done.”



✓ “Complete” the multiple cases from the case list


When you select the cases from the Work Box, the  is provided from the top of the list.

Click this, you can complete the multiple case using “Done” at once.

① Enter shipping information – In Box

- Enter the shipping information for the completed task and begin shipping.

- ② Enter shipping information.
- Click on the “Shipment” tab.



- ③ Enter shipping information:

Shipment Registration

Select a delivery method and enter the required information to register.
The information you enter is saved in the Shipment tab and can be checked at any time.

Shipment
 Self-Delivery

Company **Tracking Number**

Memo

Cancel
Register

- ④ Click the “Register” button to complete the shipment information entry.
 - The shipping information you entered will be forwarded to your partner.
- ⑤ Confirmation and Completion of Shipping – Order Box
 - After receiving the product proceed with the delivery confirmation.
 - The status of the case changes to Completed and the entire order process is completed.

Job Complete – Case Box, Order Box

After the case completes from the lab, the clinic needs to confirm the case to complete the whole order process.

If the lab completes their work, clinic can find the “Confirm Delivery” function in the Order Box.


Click this, whole order process completes, and case status will be changed to “Case Completed”



✓ Complete the case for none ordered case in the Case Box

For cases that are not ordered to the Lab, you can complete these cases from the Case Box.

- Go to Case Box.
- Select the case that want to complete.



- Click .
- Click “Case Complete.”