User Guide

i900 classic

Revision 2 (November 2025)

Document No.: UG-00003

CONTENTS

02	English
16	Български
30	中文(简体)
44	中文(繁體)
58	Hrvatski
72	čeština
86	Dansk
100	Dutch
114	Eesti keel
128	Suomi
142	Français
156	Deutsch
170	Ελληνικά
184	magyar
198	íslenska
212	Italiano
226	한국어
240	Latviešu
254	Lietuvis
268	Norsk
282	Português
296	Română
310	русский
324	Slovenčina
338	Slovenščina
352	Español
366	Svenska
380	Türk
394	українська

English

About This Guide	4
1 Introduction and Overview	4
1.1 Intended Use	4
1.2 Indication for Use	4
1.3 Contraindications	4
1.4 Qualifications of the Operator	4
1.5 Symbols	4
1.6 i900 classic Components Overview	5
1.6.1 Additional Components (Sold Separately)	6
1.7 Setting Up the i900 classic System	6
1.7.1 Basic Settings of i900 classic (Medit Plug & Scan)	6
2 Medit Scan for Clinics Overview	7
2.1 Introduction	7
2.2 Installation	7
2.2.1 System Requirements	7
2.2.2 Software Installation Guide	8
2.2.3 Medit Scan for Clinics User Guide	8
2.3 Error Messages	8
3 Maintenance	8
3.1 Calibration	8
3.1.1 How to Calibrate i900 classic	8
3.2 Cleaning, Disinfection, Sterilization Procedure	9
3.2.1 Reusable Tip	9
3.2.2 Mirror	9
3.2.3 Handpiece	9
3.2.4 Other Components	10
3.3 Disposal	10
3.4 Updates in Medit Scan for Clinics	10
4 Safety Guide	10
4.1 System Basics	10
4.2 Proper Training	10
4.3 In Case of Equipment Failure	11
4.4 Hygiene	11
4.5 Electrical Safety	11
4.6 Eye Safety	11
4.7 Explosion Hazards	11
4.8 Pacemaker and ICD Interference Risk	11
4.9 Cybersecurity Safety	11
5 Electro-Magnetic Compatibility Information	12
5.1 Electromagnetic Emissions	12
5.2 Electromagnetic Immunity	12
6 Specifications	15

About This Guide

Convention in This Guide

This guide uses various symbols to highlight important information to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

WARNING
The WARNING symbol indicates information that, if ignored, could result in a medium risk of personal injury.

CAUTION
The CAUTION symbol indicates safety information that, if ignored, could result in a slight risk of personal injury, property damage, or damage to the system.

TIPS
The TIPS symbol indicates hints, tips, and additional information for optimal operation of the system.

1. Introduction and Overview

1.1 Intended Use

The i900 classic system is an intraoral 3D scanner intended to record topographical characteristics of teeth and surrounding tissues digitally. The i900 classic system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations.

1.2 Indication for Use

The i900 classic system is for scanning the patient's intraoral features. Various factors (intraoral environment, operator's expertise, and laboratory workflow) may affect the final scan results when using the i900 classic system. The i900 classic system can also be used in full-arch scans, but various factors (intraoral environment, operator's expertise, and laboratory workflow) may affect the final results.

1.3 Contraindications

The i900 classic system is not intended to be used to create images of the internal structure of teeth or the supporting skeletal structure.






1.4 Qualifications of the Operator

- The i900 classic system is designed for use by individuals with professional knowledge in dentistry and dental laboratory technology.
- The user of the i900 classic system is solely responsible for determining whether or not this device is suitable for a particular patient's case and circumstances.
- The user is solely responsible for the accuracy, completeness, and adequacy of all data entered the i900 classic system and the provided software. The user should check the accuracy of the results and assess each individual case.
- The i900 classic system must be used in accordance with its accompanying User Guide.
- Improper use or handling of the i900 classic system will void its warranty. If you require additional information on the proper use of the i900 classic system, please contact your local distributor.
- The user is not allowed to modify the i900 classic system.

1.5 Symbols

No.	Symbol	Description
1		Serial number
2		Medical device
3		Date of manufacture
4		Manufacturer
5		Caution
6		Warning
7		Read the user guide
8		The official mark of the Europe Certificate
9		Authorized representative in the European Union
10		BF type of applied part
11		WEEE mark
12		Prescription use (U.S.A)
13		MET mark
14		AC
15		DC
16		Temperature limitation: -10 – 50°C (14 – 122°F)
17		Humidity limitation

18		Atmospheric pressure limitation
19		Fragile
20		Keep dry
21		This way up
22		Stacking more than fourteen layers is prohibited
23		Consult instructions for use
24		The official mark of the United Kingdom
25		Authorized representative in Switzerland
26		Authorized representative in the United Kingdom
27		Model Number
28		Quantity
29		Unique device identifier
30		Sterilized using ethylene oxide
31		Do not re-sterilize
32		Do not re-use
33		Use-by date
34		Do not use if package is damaged
35		Seven-Layer stacking prohibited
36		Keep away from sunlight

1.6 i900 classic Components Overview

No.	Item (Model Name)	Qty	Appearance
1	i900 classic Handpiece (MO4-i900c)	1ea	
2	i900 classic Handpiece Cover (MO1-HC1)	1ea	
3	Reusable Tip (Large) (MO1-RTL)	2ea	
4	Reusable Tip (Medium) (MO1-RTM)	2ea	
5	Calibration Tool (MO1-CT1)	1ea	
6	Practice Model	1ea	
7	Desktop Cradle (MO1-DC)	1ea	
8	Wall Mount Holder (MO1-WH1)	1ea	
9	Power Delivery Cable (2.5 m)	1ea	
10	Power Delivery Cable (2 m)	1ea	
11	USB Flash Drive (Installer of Medit Scan for Clinics included)	1ea	
12	User Guide	1ea	

1.6.1 Additional Components (Sold Separately)

No.	Item (Model Name)	Qty	Appearance
1	Reusable Tip (Large) (MO1-RTL)	4ea	
2	Reusable Tip (Medium) (MO1-RTM)	4ea	
3	Reusable Tip (Small) (MO1-RTS)	4ea	
4	i900 Series Single-Use Tip (MO1-ST1)	50ea	

- All components listed in 1.6 & 1.6.1 can be purchased separately.
- Availability of items for sale may vary depending on the medical device registration status in each country or region. Please contact Medit or your local distributor for the availability of specific items.
- For more information about the disposable tips, please refer to the i900 Series Single-Use Tip User Guide.


CAUTION

- Keep the practice model in a cool place away from direct sunlight. A discolored practice model may affect the results of the practice mode.
- Medit Scan for Clinics is included in the USB drive. This product is optimized for PC, and using other devices is not recommended. Do not use anything other than a USB port. It may cause malfunction or fire.

1.7 Setting Up the i900 classic System


1.7.1 Basic Settings of i900 classic (Medit Plug & Scan)

You can connect i900 classic directly to a PC with the power delivery cable.

Turning On the i900 classic

- ① Connect the i900 classic to your PC with the power delivery cable, and the scanner will automatically power on.
- ② When power is applied, the rear LED lights up.


Turning Off the i900 classic


When you disconnect the power delivery cable, the scanner switches off.

Desktop Cradle

Wall Mount Holder

2. Medit Scan for Clinics Overview

2.1 Introduction

Medit Scan for Clinics provides a user-friendly working interface to digitally record topographical characteristics of teeth and surrounding tissues using the i900 classic system.

2.2 Installation

2.2.1 System Requirements

Recommended System Requirements

	Windows OS		macOS
	Laptop	Desktop	Laptop/Desktop
CPU	Intel Core i7-13700H Intel Core i7-12700H AMD Ryzen 7 7735H AMD Ryzen 7 6800H	Intel Core i7-13700K Intel Core i7-12700K AMD Ryzen 7 7700X AMD Ryzen 7 5800X	M1 Pro (10-core CPU, 16-core GPU) M2 (8-core CPU, 10-core GPU) M2 Pro (10-core CPU, 16-core GPU)
RAM	32 GB		24 GB
Graphic	NVIDIA GeForce RTX 4060 (VRAM 8 GB or higher) NVIDIA GeForce RTX 3070 (VRAM 8 GB or higher) NVIDIA RTX A3000 (VRAM 8 GB or higher) * AMD Radeon is not supported.		
OS	Windows 10 64-bit Windows 11 (recommended for 12th Gen or later Intel Core processors)		Monterey 12 Ventura 13

Minimum System Requirements

	Windows OS		macOS
	Laptop	Desktop	Laptop/Desktop
CPU	Intel Core i5-13500H Intel Core i5-12500H AMD Ryzen 5 7535HS AMD Ryzen 5 6600H	Intel Core i5-13400 Intel Core i5-12400 AMD Ryzen 5 7500 AMD Ryzen 5 5600	M1 (8-core CPU, 7-core GPU) M2 (8-core CPU, 8-core GPU)
RAM	16 GB		16 GB
Graphic	NVIDIA GeForce RTX 4050 (VRAM 6 GB or higher) NVIDIA GeForce RTX 3060 (VRAM 6 GB or higher) NVIDIA RTX A2000 (VRAM 6 GB or higher) * AMD Radeon is not supported.		
OS	Windows 10 64-bit Windows 11 (recommended for 12th Gen or later Intel Core processors)		Monterey 12 Ventura 13

For accurate and up-to-date system requirements, please visit www.meditlink.com.

Use PC and monitor classified as class I and certified IEC 62368-1 (or IEC 60950-1), IEC 55032, IEC 55024.

- To prevent electric shock, connect the PC only to a properly grounded power source.
- Please ensure that the PC's power plug remains easily accessible at all times.

The device may not work when using cables other than the USB 3.0 cable provided by Medit. Medit is not responsible for any problems caused by cables other than the USB 3.0 cable provided by Medit. Be sure to use only the USB 3.0 cable included in the package.

2.2.2 Software Installation Guide

- 1 Connect the enclosed USB flash drive to a PC.
- 2 Run the installation file.
- 3 Select the setup language and click "Next."
- 4 Choose the installation path.
- 5 Read the "License Agreement" carefully, check "I agree to the License terms and conditions," and then click "Install."
- 6 The installation process may take several minutes. Please do not shut down the PC until the installation is complete.
- 7 After the installation is complete, restart the PC to ensure optimal program operation.



The installation will not be processed if the i900 classic system is connected to a PC. Please turn off the scanner before starting the installation.

2.2.3 Medit Scan for Clinics User Guide

Please refer to the User Guide of Medit Scan for Clinics: Medit Scan for Clinics > Menu > User Guide.

2.3 Error Messages

The Medit i900 classic scanner system presents error messages, enabling users to intuitively identify hardware and system conditions. Some messages include solutions to aid users in self-diagnosis and problem resolution. Error messages are presented in plain text and supplemented with images when needed to enhance comprehension.

CAUTION

If the instructions provided in the error message do not address the issue, please contact the local distributor from whom you acquired your scanner system or support@medit.com.

3. Maintenance

CAUTION

- Equipment maintenance should only be carried out by a Medit employee or Medit-Certified company or personnel.
- In general, users are not required to perform maintenance work on the i900 classic system besides calibration, cleaning, and sterilization. Preventive inspections and other regular maintenance are not required.

3.1 Calibration

Periodic calibration is required to produce precise 3D models. You should perform calibration when:

- The quality of the 3D model is not reliable or accurate when compared to previous results.
- Environmental conditions such as temperature have changed.
- The calibration period has expired. You can set the calibration period in the Menu > Settings > Calibration Period (Days).

The calibration panel is a delicate component.

Do not touch the panel directly. Check the calibration panel if the calibration process is not performed properly. If the calibration panel is contaminated, please contact your service provider.

If the calibration tool is long exposed to low temperatures below Medit's recommendations mentioned in the use and storage guidance, the rotation motion of the device may be compromised. In that case, forcibly turning the device may cause damage. To avoid it, keep the calibration device in the recommended temperature environment before use.

We recommend performing a calibration periodically.

You can set the calibration period via Menu > Settings > Calibration Period (Days). The default calibration period is 90 days.

3.1.1 How to Calibrate i900 classic

- 1 Turn on the i900 classic and launch the Medit Scan for Clinics.
- 2 Run the Calibration Wizard at the bottom of the main toolbar panel in Medit Scan for Clinics.
- 3 Prepare the calibration tool and the i900 classic handpiece.
- 4 Turn the dial of the calibration tool to the starting position.
- 5 Put the i900 classic handpiece into the calibration tool.
- 6 Click "Next" to start the calibration process.
- 7 If the calibration tool is mounted properly in the correct position, the system automatically acquires data.
- 8 When data acquisition is complete at the starting position, turn the dial to the next position.
- 9 Repeat the steps to the last position.
- 10 When data acquisition is complete at the last position, the system automatically calculates and shows the calibration results.

3.2 Cleaning, Disinfection, Sterilization Procedure

3.2.1 Reusable Tip

The reusable tip is the part that is positioned in the mouth of patients during scanning, and is reusable for a limited number of times. The tip needs to be cleaned and sterilized between patient uses to avoid cross-contamination.

Cleaning & Disinfection

- Prepare a cleaning solution.
 - » Dilute a general neutral detergent at a 1:100 ratio before use.
- Clean the reusable tips with a cleaning solution and a brush.
 - » Ensure the mirror of the tip is completely clean and stain-free. If the mirror appears stained or foggy, repeat the cleaning process.

CAUTION

- » The reusable tip has a complex structure, and automatic cleaning may not achieve complete cleaning; therefore, do not clean the reusable tip in an automatic washer.
- Rinse the reusable tips three times with purified water.
- Remove moisture with a paper towel and allow them to air-dry completely at room temperature for at least 80 minutes.
- Disinfect the reusable tips using a disinfectant containing 15% or less Isopropyl Alcohol (IPA) for 1 minute. Then, ensure they are thoroughly dried at room temperature for at least 5 minutes.
 - » Before using a disinfectant product, please refer to the product user manual for proper use.
 - » You can find the list of recommended disinfectants in the Medit Help Center at <http://support.medit.com/hc>.

Sterilization

- The tip should be cleaned manually using a disinfecting solution. After cleaning and disinfecting, inspect the mirror inside the tip to ensure there aren't any stains or smudges.
- Repeat the cleaning and disinfection process if necessary. Carefully dry the mirror with a paper towel.
- Insert the tip into a paper sterilization pouch and seal it, making sure that it is airtight. Use either a self-adhesive or heat-sealed pouch.
- Sterilize the wrapped tip in an autoclave with the following conditions:
 - » Sterilize in a gravity autoclave at 135°C (275°F) for 10 minutes and dry for 30 minutes.
 - » Sterilize in a pre-vacuum autoclave at 134°C (273.2°F) for 4 minutes and dry for 20 minutes.
- Use an autoclave program that dries the wrapped tip before opening the autoclave.
- Scanner tips can be re-sterilized up to 150 times. After reaching this limit, they must be disposed of according to the guidelines in the disposal section.
- Autoclave times and temperatures may vary depending on the autoclave type and manufacturer. For this reason, it may not be able to meet the maximum number of times. Please refer to the user's manual of the autoclave manufacturer you are using to determine whether the required conditions are met.

CAUTION

- The mirror of the scanner tip is a delicate optical component that should be handled with care to ensure optimal scan quality. Be careful not to scratch or smudge it as any damage or blemishes may affect the data acquired.
- Make sure to always wrap the tip before autoclaving. If you autoclave an exposed tip, it will cause stains on the mirror, which cannot be removed. Check the autoclave manual for more information.
- Tips that have been cleaned, disinfected, and sterilized must remain sterile until they are used on the patient.
- Medit is not responsible for any damage, including tip distortion, caused by cleaning, disinfection, or sterilization procedures that do not comply with the guidelines outlined above.

3.2.2 Mirror

The presence of impurities or smudges on the tip mirror may lead to poor scan quality and an overall poor scanning experience. In such situations, clean the mirror following the steps below:

- ① Disconnect the scanner tip from the i900 classic handpiece.
- ② Pour alcohol on a clean cloth or cotton-tipped swab and wipe the mirror. Make sure to use alcohol that is free of impurities or it may stain the mirror. You can use either ethanol or propanol (ethyl-/propyl alcohol).
- ③ Wipe the mirror dry using a dry, lint-free cloth.
- ④ Make sure the mirror is free of dust and fibers. Repeat the cleaning process as necessary.

3.2.3 Handpiece

After treatment, clean and disinfect all other surfaces of the i900 classic handpiece except for the scanner front (optical window) and end (air vent hole). Cleaning and disinfecting must be done with the device turned off. Use the device only after it is completely dry.

The recommended cleaning and disinfecting solution is denatured alcohol (ethyl alcohol or ethanol) – typically 60 – 70% Alc/Vol.

The general cleaning and disinfecting procedures are as follows:

- ① Turn off the power by unplugging the power delivery cable from the device.
- ② Clean the filter on the front end of the i900 classic handpiece.
 - » If alcohol is poured directly into the filter, it may seep inside the i900 classic handpiece and cause a malfunction.
 - » Do not clean the filter by pouring alcohol or cleaning solution directly into the filter. The filter must be gently wiped with a cotton or soft cloth moistened with alcohol.
 - » Do not wipe by hand or apply excessive force.
 - » Medit is not responsible for any damage or malfunction that occurs during cleaning that does not follow the guidelines above.
- ③ After cleaning the filter, put the cover on the front of the i900 classic handpiece.
- ④ Pour the disinfectant onto a soft, lint-free, and non-abrasive cloth.
- ⑤ Wipe the scanner surface with the cloth.
- ⑥ Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

CAUTION

- Do not clean the i900 classic handpiece when the device is turned on as the fluid may enter the scanner and cause malfunction.
- Use the device after it is completely dry.
- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

3.2.4 Other Components

- Pour the cleaning and disinfecting solution onto a soft, lint-free, and non-abrasive cloth.
- Wipe the component surface with the cloth.
- Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

CAUTION

- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

3.3 Disposal

CAUTION

- The scanner tip must be sterilized before disposal. Sterilize the tip as described in the section "3.2.1 Reusable Tip."
- Dispose of the scanner tip as you would any other clinical waste.
- Other components are designed to conform with the following directives:
RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)
WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

3.4 Updates in Medit Scan for Clinics

Medit Scan for Clinics automatically checks for updates when the software is in operation. If a new version of the software is released, the system will automatically download it.

4. Safety Guide

Please adhere to all the safety procedures as detailed in this User Guide to prevent human injury and equipment damage. This document uses the words WARNING and CAUTION when highlighting precautionary messages.

Carefully read and understand the guidelines, including all preventive messages as prefaced by the words WARNING and CAUTION. To avoid bodily injury or equipment damage, make sure to adhere strictly to the safety guidelines. All instructions and precautions as specified in the Safety Guide must be observed to ensure the system's proper functionality and personal safety.

The i900 classic system should only be operated by dental professionals and technicians who are trained to use the system. Using the i900 classic system for any purpose other than its intended usage as outlined in the section "1.1 Intended Use" may result in injury or damage to the equipment. Please handle the i900 classic system according to the guidelines in the safety guide.

Any serious incident involving the device should be reported to the manufacturer and the competent authority of the member state in which the user and patients are located.

4.1 System Basics

The i900 classic system is a high-precision optical medical device. Acquaint all the following safety and operating instructions before the installation, usage, and operation of i900 classic.

CAUTION

- If the product has been stored in a cold environment, give it time to adjust to the temperature of the environment before use.
- If used immediately, condensation may occur, which may damage the electronic parts inside the unit.
- Ensure that all components provided are free from physical damage. Safety cannot be guaranteed if there is any physical damage to the unit.
- Before using the system, check that there are no issues such as physical damage or loose parts. If there is any visible damage, do not use the product and contact the manufacturer or your local representative.
- Check the i900 classic handpiece and its accessories for any sharp edges.
- When not in use, the i900 classic system should be kept mounted on a desk stand or a wall mount stand.
- Do not install the desk stand on an inclined surface.
- Do not place any object on the i900 classic system.
- Do not place the i900 classic system on any heated or wet surface.
- Do not block the air vents located at the rear of the i900 classic system. If the equipment overheats, the i900 classic system may malfunction or stop working.
- Do not spill any liquid on the i900 classic system.
- The i900 classic handpiece and other included components are made of electronic components. Do not allow any kind of liquid or foreign objects to enter.
- Do not pull or bend the cable connected to the i900 classic system.
- Carefully arrange all the cables so that you or your patient do not trip or get caught in the cables. Any pulling tension on the cables may cause damage to the i900 classic system.
- Always keep an eye on the product and your patient while using the product to check for abnormalities.
- Proceed with calibration, cleaning, disinfection, and sterilization in accordance with the contents of the user guide.
- If you drop the tips on the floor, do not attempt to reuse them. Discard the tip immediately as there is a risk that the mirror attached to the tip may have been dislodged.
- Due to its fragile nature, the tips should be handled with care. To prevent damage to the tip and its internal mirror, be careful to avoid contact with the patient's teeth or restorations.
- If the i900 classic system is dropped on the floor or if the unit is impacted, it must be calibrated before use. If the instrument is unable to connect to the software, consult the manufacturer or authorized resellers.
- If the equipment fails to operate normally, such as having issues with accuracy, stop using the product, and contact the manufacturer or authorized resellers.
- Install and use only approved programs to ensure the proper functionality of the i900 classic system.
- In the event of a severe accident involving the i900 classic system, notify the manufacturer and report it to the competent national authority of the country or region where the user and patient reside.
- If the PC with the software installed does not have security software or if there is a risk of malicious code intrusion into the network, the PC may be breached with malware (malicious software such as viruses or worms that damage your computer).
- The software for this product must be used in compliance with medical and personal information protection laws.
- To provide electrical insulation and maintain electrical safety, a coating is applied to insulate the device except on the areas where the USB ports are located.
- Check the product's markings, such as certification symbols or manufacturer details, under lighting of about 500 lux for clear visibility.

4.2 Proper Training

WARNING

Before using your i900 classic system on patients:

- You should have been trained to use the system or read and fully understand this User Guide.
- You should be familiar with the safe use of the i900 classic system, as detailed in this User Guide.
- Before use or after changing any settings, the user should check that the live image is displayed properly in the camera preview window of the program.

4.3 In Case of Equipment Failure

WARNING

If your i900 classic system is not working properly or if you suspect that there is a problem with the equipment:

- Remove the device from the patient's mouth and discontinue use immediately.
- Disconnect the device from the PC and check for errors.
- Contact the manufacturer or authorized resellers.
- Modifications to the i900 classic system are prohibited by law as they may compromise the safety of the user, patient, or a third party.

4.4 Hygiene

WARNING

For clean working conditions and patient safety, ALWAYS wear clean surgical gloves when:

- Making any direct contact with the tip.
- Using the i900 classic system on patients.
- Touching the i900 classic system.

WARNING

The i900 classic system and its optical window should be always kept clean. Before using the i900 classic system on a patient, be sure to:

- Sterilize the i900 classic system as described in the section "3.2 Cleaning, Disinfection, Sterilization Procedure."
- Use a sterilized tip.

4.5 Electrical Safety

WARNING

- The i900 classic system is a Class I device.
- The i900 classic system only uses RF energy internally. The amount of RF radiation is low and does not interfere with surrounding electromagnetic radiation.
- There is a risk of electric shock if you attempt to access the inside of the i900 classic system. Only qualified service personnel should access the system.
- Do not place liquids such as beverages near the i900 classic system and avoid spilling any liquid on the system.
- Never spill liquid of any kind on the i900 classic system.
- Condensation due to changes in temperature or humidity can cause moisture buildup inside the i900 classic system, which may damage the system. Before connecting the i900 classic system to a power supply, be sure to keep the i900 classic system at room temperature for at least two hours to prevent condensation. If condensation is visible on the product surface, the i900 classic should be left at room temperature for more than 8 hours.
- The i900 classic system doesn't have a physical power button. Turn off the scanner by unplugging the power delivery cable.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A).
- If used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services.
- Avoid pulling on the power delivery cables used with the i900 classic system.
- Do not touch the connectors of the device and the patient simultaneously.

4.6 Eye Safety

WARNING

- The i900 classic system projects a bright light from its tip during scanning.
- The bright light projected from the tip of the i900 classic is not harmful to the eyes. However, you should not look directly at the bright light nor aim the light beam into the eyes of others. Generally, intense light sources can cause eyes to become brittle and the likelihood of secondary exposure is high. As with other intense light source exposure, you may experience a temporary reduction in visual acuity, pain, discomfort, or visual impairment, all of which increase the risk of secondary accidents.
- There is an LED that emits UV-C wavelengths inside the i900 classic handpiece. It is irradiated only inside the i900 classic handpiece and does not go outside. The blue light visible inside the i900 classic handpiece is for guidance, not UV-C light. It is harmless to the human body.
- The UV-C LED operates with a wavelength of 270 – 285 nm.
- Disclaimer for Risks Involving Patients with Epilepsy The Medit i900 classic should not be used on patients that have been diagnosed with epilepsy due to the risk of seizures and injury. For the same reason, dental staff who have been diagnosed with epilepsy should not operate the Medit i900 classic.

4.7 Explosion Hazards

WARNING

- The i900 classic system is not designed to be used near flammable liquids, gases, or in environments with high oxygen concentrations.
- There is a risk of explosion if you use the i900 classic system near flammable anesthetics.

4.8 Pacemaker and ICD Interference Risk

WARNING

- Implantable Cardioverter Defibrillators (ICDs) and pacemakers may have interference due to some devices.
- Maintain a moderate distance from the patient's ICD or pacemaker when using the i900 classic system.
- For more information on peripherals used with i900 classic, check the respective manufacturer's manuals.

4.9 Cybersecurity Safety

- If a cybersecurity incident occurs, stop using the scanner and software immediately. Power off the scanner and log out of the software.
- Promptly report the incident to our support team through email, phone, or other available means of contact. Please refer to the last page of the User Guide for contact information.
- When reporting an incident, please provide as much information as possible, including the time of occurrence and any unusual behavior you noticed. This information will assist us in resolving the issue quickly.

5. Electro-Magnetic Compatibility Information

5.1 Electromagnetic Emissions

The i900 classic system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i900 classic system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Emission			
Emission Test	Compliance	Electromagnetic Environment - Guidance	
RF Emissions CISPR 11	Group 1	The i900 classic uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class A	The i900 classic is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
Harmonic Emissions IEC 61000-3-2	Class A		
Voltage Fluctuations / Flicker Emissions	Complies		

WARNING

This i900 classic system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the i900 classic or shielding the location.

5.2 Electromagnetic Immunity

Guidance 1

The i900 classic system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i900 classic system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with a synthetic material, relative humidity of at least 30% is recommended.
Electrical Fast Transient / Burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	The mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV differential mode	± 0.5 kV, ± 1 kV differential mode	The mains power quality should be that of a typical commercial or hospital environment.
	± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV, ± 2 kV common mode	
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC 61000-4-11	0% Ur (100% dip in Ur) for 0.5/1 cycles 70% Ur (30% dip in Ur) for 25/30 cycles 0% Ur (100% dip in Ur) for 250/300 cycles	0% Ur (100% dip in Ur) for 0.5/1 cycles 70% Ur (30% dip in Ur) for 25/30 cycles 0% Ur (100% dip in Ur) for 250/300 cycles	The mains power quality should be that of a typical commercial or hospital environment. If the user of the i900 classic system requires continued operation during power mains interruptions, it is recommended that the i900 classic system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Fields (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a location in a typical commercial or hospital environment.
Proximity Magnetic Fields in the Frequency Range 9 kHz to 13.56 MHz Immunity IEC 61000-4-39	8 A/m 30 kHz CW modulation	8 A/m 30 kHz CW modulation	Resistance to magnetic fields was tested and applied only to surfaces of enclosures or accessories accessible during intended use.
	65 A/m 134.2 kHz PM 2.1 kHz	65 A/m 134.2 kHz PM 2.1 kHz	
	7.5 A/m 13.56 MHz PM 50 kHz	7.5 A/m 13.56 MHz PM 50 kHz	

NOTE: Ur is the main voltage (AC) prior to the application of the test level.

Guidance 2

Recommended Separation Distances Between Portable and Mobile Communication Equipment and the i900 classic			
Rated Maximum Output Power of the Transmitter [W]	Separation Distance According to the Frequency of Transmitter [M]		
	IEC 60601-1-2:2014		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$		80 MHz to 2.7 GHz $d = 2.0 \sqrt{P}$
0.01	0.12		0.20
0.1	0.38		0.63
1	1.2		2.0
10	3.8		6.3
100	12		20

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Guidance 3

The i900 classic system is intended for use in the electromagnetic environment specified below. The customer or the user of the i900 classic system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz Outside ISM Bands amateur	3 Vrms	Portable and mobile RF communications equipment should not be used closer to any part of the Ultrasound System, including cables, than the recommended separation distance. This is calculated using the equation applicable to the frequency of the transmitter. Recommended Separation Distance (d): $d = 1.2 \sqrt{P}$ IEC 60601-1-2:2007 $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80 MHz to 2.5 GHz IEC 60601-1-2:2014 $d = 2.0 \sqrt{P}$ 80 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
	6 Vrms 150 kHz to 80 MHz In ISM Bands amateur	6 Vrms	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

NOTE 3: The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

• **Guidance 4**

The i900 classic system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the i900 classic system. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity					
Immunity Test	Band ¹⁾	Service ¹⁾	Modulation	IEC 60601 Test Level	Compliance Level
Proximity Fields from RF Wireless Communications IEC 61000-4-3	380 – 390 MHz	TETRA 400	Pulse Modulation 18 Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460; FRS 460	FM ±5 kHz Deviation 1 kHz sine	28 V/m	28 V/m
	704 – 787 MHz	LTE Band 13, 17	Pulse Modulation 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM 800;900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	Pulse Modulation 18 Hz	28 V/m	28 V/m
	1700 – 1990 MHz	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217 Hz	28 V/m	28 V/m
	2400 – 2570 MHz	Bluetooth; WLAN 802.11b/g/n; RFID 2450; LTE Band 7	Pulse Modulation 217 Hz	28 V/m	28 V/m
	5100 – 5800 MHz	WLAN 802.11a/n	Pulse Modulation 217 Hz	9 V/m	9 V/m

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

 **WARNING**

- Use of the i900 classic adjacent to or on other equipment must be avoided as it may result in improper operation. If this use is necessary, it is advisable that this and the other equipment be observed to verify that they are operating normally.
- The use of accessories, transducers, and cables other than those specified or provided by the Medit of the i900 classic could result in high electromagnetic emissions or reduced electromagnetic immunity from this equipment and result in improper operation.

¹⁾ For some services, only the uplink frequencies are included.

6. Specifications

Model Name	MO4-i900c	
Trade Name	i900 classic	
Packing Unit	1 set	
Rating	5 V ⁺ 3 A	
Classifications for Protection Against Electric Shock	Class I, Type BF Applied Parts (Reusable Tip)	
* This product is a medical device.		
Handpiece (including medium-size tip)		
Dimension	223.4 x 36.7 x 35.3 mm (L x W x H)	
Weight	165 g	
Reusable Tip		
Dimension - Full Tip	Large	36.1 x 34.1 x 90.8 mm (W x H x L)
	Medium	36.1 x 34.1 x 90.4 mm (W x H x L)
	Small	36.1 x 34.1 x 90.3 mm (W x H x L)
Dimension - Tip Head	Large	26.9 x 19.7 mm (W x H)
	Medium	22.4 x 16.3 mm (W x H)
	Small	18.36 x 13.1 mm (W x H)
Calibration Tool		
Dimension	160 x 48.5 mm (H x Ø)	
Weight	205 g	
Operating, Storage, and Transport Conditions		
Operating Condition	Temperature	18 – 28°C (64.4 – 82.4°F)
	Humidity	20 – 75% relative humidity (non-condensing)
	Air Pressure	800 – 1,100 hPa
Storage Condition	Temperature	-10 – 50°C (14 – 122°F)
	Humidity	20 – 80% relative humidity (non-condensing)
	Air Pressure	800 – 1,100 hPa
Transport Condition	Temperature	-10 – 50°C (14 – 122°F)
	Humidity	20 – 80% relative humidity (non-condensing)
	Air Pressure	620 – 1,200 hPa
Emission limits per environment		
Environment	Hospital environment	
Conducted and radiated RF Emissions	CISPR 11	

EU REP

Meditrial Srl
Via Po 9 00198, Rome Italy

CH REP

Meditrial Europe Ltd
Bahnhofstrasse 23 6300 Zug, Switzerland

UK REP

Meditrial UK Ltd
19 The Circle, Queen Elizabeth Street, London, SE1 2JE, United Kingdom

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Connection failed.

Failed to connect to the cameras. Please follow the instructions below and try again.

- Press the "Retry" button. Device will automatically reboot and try to connect again.
- Turn off the device and disconnect all cables from the PC. Connect cables again and restart the device.
- If connection problem persists, please install the program again.
- If connection problem persists, please contact your service provider.

Cancel Retry

Turn Off the External Light

External light was detected while scanning.

This scan might have proceeded with an external light source, which would increase the chances of noisy and inaccurate scanned data.

We recommend turning the external light off before the next scan.

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