# User Guide

***i700 wireless***

**Revision 5 (November 2025)**

Document No.: ME-UG-7024

## CONTENTS

02	English
22	العربية
42	中文(简体)
62	Français
82	Deutsch
102	Italiano
122	Português
142	русский
162	Español
182	한국어

**English**

<b>About This Guide</b> .....	4
<b>1 Introduction and Overview</b> .....	4
1.1 Intended Use .....	4
1.2 Indication for Use .....	4
1.3 Contraindications .....	4
1.4 Qualifications of the Operating User .....	4
1.5 Symbols .....	4
1.6 i700 wireless Components Overview .....	5
1.6.1 Wireless Hub Package (*Sold Separately) .....	6
1.7 Setting Up the i700 wireless System .....	6
1.7.1 Basic Settings of i700 wireless .....	6
1.7.2 How to Use the Wireless Hub .....	8
1.7.3 How to Use the Battery .....	8
1.7.4 How to Use the Battery Charger .....	8
1.7.5 How to Mount the Handle .....	9
1.7.6 How to Place on the Desktop Cradle .....	9
1.7.7 How to Install the Wall Mount Holder .....	9
<b>2 Medit Scan for Clinics Overview</b> .....	10
2.1 Introduction .....	10
2.2 Installation .....	10
2.2.1 System Requirements .....	10
2.2.2 Software Installation Guide .....	10
2.2.3 Medit Scan for Clinics User Guide .....	10
<b>3 Maintenance</b> .....	11
3.1 Calibration .....	11
3.1.1 How to Calibrate i700 wireless .....	11
3.1.2 Auto Calibration Tool (Sold Separately) .....	11
3.2 Cleaning, Disinfection, Sterilization Procedure .....	11
3.2.1 Reusable Tip .....	11
3.2.2 Mirror .....	12
3.2.3 Handpiece .....	12
3.2.4 Other Components .....	12
3.3 Disposal .....	12
3.4 Battery Storage .....	13
3.5 The Battery Usage Precautions and Disposal Guide .....	13
3.6 Updates in Medit Scan for Clinics .....	13
<b>4 Safety Guide</b> .....	14
4.1 System Basics .....	14
4.2 Proper Training .....	14
4.3 In Case of Equipment Failure .....	15
4.4 Hygiene .....	15
4.5 Electrical Safety .....	15
4.6 Eye Safety .....	15
4.7 Explosion Hazards .....	16
4.8 Pacemaker and ICD Interference Risk .....	16
4.9 Cybersecurity Safety .....	16
<b>5 Electro-Magnetic Compatibility Information</b> .....	16
5.1 Electromagnetic Emissions .....	16
5.2 Electromagnetic Immunity .....	16
<b>6 Wireless Compliance Information</b> .....	19
6.1 FCC Compliance Statement .....	19
6.2 IC Compliance Statement .....	19
6.3 KC Compliance Statement .....	19
6.4 TELEC (Japan) Compliance Statement .....	19
<b>7 Specifications</b> .....	20

## About This Guide

### Convention in This Guide

This User Guide uses various symbols to highlight important information to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

#### WARNING

The WARNING symbol indicates information that, if ignored, could result in a medium risk of personal injury.

#### CAUTION

The CAUTION symbol indicates safety information that, if ignored, could result in a slight risk of personal injury, property damage, or damage to the system.

#### TIPS

The TIPS symbol indicates hints, tips, and additional information for optimal operation of the system.

## 1. Introduction and Overview

### 1.1 Intended Use

The i700 wireless system is an intraoral 3D scanner intended to record topographical characteristics of teeth and surrounding tissues digitally. The i700 wireless system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations.

### 1.2 Indication for Use

The i700 wireless system is for scanning the patient's intraoral features. Various factors (intraoral environment, operator's expertise, and laboratory workflow) may affect the final scan results when using the i700 wireless system.

### 1.3 Contraindications



The i700 wireless system is not intended to be used to create images of the internal structure of teeth or the supporting skeletal structure.





### 1.4 Qualifications of the Operating User

#### CAUTION

- The i700 wireless system is designed for use by individuals with professional knowledge in dentistry and dental laboratory technology.
- The user of the i700 wireless system is solely responsible for determining whether or not this device is suitable for a particular patient's case and circumstances.
- The user is solely responsible for the accuracy, completeness, and adequacy of all data entered the i700 wireless system and the provided software. The user should check the accuracy of the results and assess each individual case.
- The i700 wireless system must be used in accordance with its accompanying User Guide.
- Improper use or handling of the i700 wireless system will void its warranty. If you require additional information on the proper use of the i700 wireless system, please contact your local distributor.
- The user is not allowed to modify the i700 wireless system.

### 1.5 Symbols

No	Symbol	Description
1		Serial number
2		Medical device
3		Date of manufacture
4		Manufacturer
5		Caution
6		Warning
7		Read the user guide
8		The official mark of the Europe Certificate
9		Authorized representative in the European community
10		BF type of applied part
11		WEEE mark
12		Prescription use (U.S.A)
13		MET mark
14		AC
15		DC
16		Temperature limitation
17		Humidity limitation
18		Atmospheric pressure limitation
19		Fragile

20		Keep dry
21		This way up
22		Seven-Layer stacking prohibited
23		Consult instructions for use
24		The official mark of the United Kingdom
25		Authorized representative in Switzerland
26		Authorized representative in the United Kingdom
27		Model Number
28		Quantity
29		Unique device identifier
30		Non-sterile
31		Stand-by

### 1.6 i700 wireless Components Overview

No	Item	Qty	Appearance
1	i700 wireless Handpiece	1ea	
2	Wireless Hub	1ea	
3	Rechargeable Battery	3ea	
4	Extension Cable Battery	1ea	
5	Battery Charger	1ea	
6	Handpiece Cover	1ea	
7	Reusable Tip	4ea	
8	Small Tip (*Sold separately)	4ea	
9	Handle	1ea	
10	Calibration Tool	1ea	
11	Practice Model	1ea	
12	Wrist Strap	1ea	
13	Desktop Cradle	1ea	
14	Wall Mount Holder	1ea	
15	USB 3.0 Cable (C to A)	1ea	
16	Power Delivery Cable (C to C) – 1.0 Meter	1ea	
17	Power Delivery Cable (C to C) – 3.0 Meter (*Sold Separately)	1ea	
18	Medical Adapter for Wireless Hub	1ea	

19	Medical Adapter for Battery Charger (Same as above)	1ea	
20	Power Cord	2ea	
21	USB Flash Drive (Installer of Medit Scan for Clinics included)	1ea	
22	User Guide	1ea	
23	Monitor Mount	1ea	
24	Handle-Fit Cradle (*Sold Separately)	1ea	
25	Auto Calibration Tool (1 rechargeable battery included) (*Sold separately)	1ea	

- All components in the list can be purchased separately.
- The availability of items may vary depending on the medical device registration status in each country or region. Please contact Medit or your local distributor to check the availability of specific items.

**CAUTION**

- Keep the practice model in a cool place away from direct sunlight. A discolored practice model may affect the results of the practice mode.
- The strap is designed specifically for the weight of the i700 wireless and should not be used with other products.
- Medit Scan for Clinics is included in the USB drive. This product is optimized for PC, and using other devices is not recommended. Do not use anything other than a USB port. It may cause malfunction or fire.
- The auto-calibration tool is not included in the i700 wireless package. It can be bought separately from your place of purchase.


**1.6.1 Wireless Hub Package (\*Sold Separately)**

No	Item	Qty	Appearance
1	Wireless Hub	1ea	
2	USB 3.0 Cable (C to A)	1ea	
3	Power Delivery Cable (C to C) - 1.0 Meter	1ea	
4	Monitor Mount	1ea	
5	Medical Adapter for Battery Charger	1ea	
6	Power Cord	1ea	

**1.7 Setting Up the i700 wireless System**

**1.7.1 Basic Settings of i700 wireless**

Connect the i700 wireless Wireless Hub (1)

① Connect the USB 3.0 cable (C to A) to the wireless hub.

② Connect the adapter to the wireless hub.

③ Connecting the power cord to the adapter automatically powers the wireless hub.


④ Connect the USB 3.0 cable's A port connector to a PC.

※ The USB port is used only for signal transmission.


#### Connect the i700 wireless Wireless Hub (2)

① Connect the power delivery cable (C to C) to the wireless hub.

② Connect the power delivery cable to a PC.

※ 9 Vdc power is supplied to the USB port.

#### Turning On the i700 wireless

① Insert the battery into the i700 wireless handpiece and press the power button.

② When the power is supplied, the LED on the upper part of the i700 wireless handpiece turns blue.

③ Three LEDs on the bottom of the i700 wireless handpiece indicate the remaining battery level.

- 3 lights: 80 – 100%
- 2 lights: 50 – 79%
- 1 light: 20 – 49%
- 1 blinking light: 1 – 19%
- LED lights turned off: 0%

#### Turning Off the i700 wireless

Press and hold the power button at the bottom of the i700 wireless handpiece for 3 seconds.  
If the rechargeable battery is removed without turning the device off, it may shorten the lifespan of the i700 wireless and the battery.

### 1.7.2 How to Use the Wireless Hub

The i700 wireless is a wireless device that operates through a wireless module. For this reason, the i700 wireless handpiece has a transmitter and the wireless hub has a receiver. The i700 wireless system uses two types of frequencies to transmit data and control the i700 wireless handpiece.

Power is supplied by connecting the adapter cable to the wireless hub power port. The power turns off when the adapter port is removed.

When the i700 wireless is turned on, it will automatically try to pair with the wireless hub. Both the i700 wireless handpiece and wireless hub must be powered and placed in a range of each other to pair. When pairing is in progress, the LED on top of the i700 wireless handpiece flashes. When pairing is complete, the LED lights up. You can start scanning when the pairing is complete.

- i700 wireless uses two antenna modules: 60 GHz for data transmission and 2.4 GHz for the control. The actual frequency depends on the local regulations.
- The actual operational range is up to 5 m, and it may vary depending on the environment.
- 60 GHz Frequency: 57 – 64 GHz
- 2.4 GHz Frequency: 2.4 – 2.5 GHz

Power is supplied from the connected PC without using the adapters when connected with a power delivery cable. In this case, the PC must be turned on. If the power delivery cable is removed from the PC, the wireless hub will be automatically turned off and the connection status, such as pairing, will be initialized.

### 1.7.3 How to Use the Battery

- Rechargeable Battery
  - » Li-ion
  - » 3.6 V, 3100 mAh, 11.6 Wh
  - » Insert the battery at the bottom of the i700 wireless handpiece. Insert the battery into the i700 wireless handpiece with the terminal aligned correctly.

- » Battery life lasts up to 1 hour.
- » Battery life may vary depending on the user environment and the battery cycle count.

- Extension Cable Battery
  - » The extension cable battery is a cable-type wired battery that is used instead of a rechargeable battery.
  - » No charging is required as power is supplied through the cable.
  - » Connect the battery-shaped connector to the i700 wireless handpiece and connect the cable to the battery charger terminal.

- » Press the power button on the bottom of the i700 wireless handpiece to supply power.

### 1.7.4 How to Use the Battery Charger

- Connect the adapter to the battery charger power port to supply power. The power turns off when the adapter port is removed.

- Insert the rechargeable battery into the battery charger with the charging terminals oriented correctly.

- The battery charger is only for rechargeable batteries. It takes up to 2 hours and 30 minutes to be fully charged and the charging time may vary depending on the user environment and the battery cycle count.

The LED on the charger blinks in blue when the battery is charging. When fully charged, the LED lights up in blue.

If the battery is not correctly inserted into the battery charger, the LED on the charger will blink in red. In this case, remove the battery from the charger, check both terminals of the battery and the charger for foreign objects, wipe them gently with a soft cloth, and then insert the battery again.

### 1.7.5 How to Mount the Handle

The i700 wireless body is equipped with a wireless signal transmitter located where the i700 wireless logo is. Depending on your experiences and habits, you may hold the area where the transmitter is mounted. Covering the transmitter area may interfere with the wireless hub's communication. Therefore, a handle is provided to mount on the i700 wireless handpiece for a more comfortable hand-gripping method.

- ① Turn over the i700 wireless handpiece to find the silicone body.
- ② Remove the silicone with your hand.

- ③ Fully fasten the bolts of the handle to the handle mounting hole on the i700 wireless body.
- ④ Tighten it clockwise using the knob on the handle.

- ⑤ You can now use it by holding the handle. If you want to remove the handle, proceed in the reverse order of this instruction.

### 1.7.6 How to Place on the Desktop Cradle

Without the Handle

With the Handle

### 1.7.7 How to Install the Wall Mount Holder

## 2. Medit Scan for Clinics Overview

### 2.1 Introduction

Medit Scan for Clinics provides a user-friendly working interface to digitally record topographical characteristics of teeth and surrounding tissues using the i700 wireless system.

### 2.2 Installation

#### 2.2.1 System Requirements

##### Minimum System Requirements

	Windows		macOS
	Laptop	Desktop	Laptop/Desktop
CPU	Intel Core i5 – 13500H Intel Core i5 – 12500H	Intel Core i5 – 13400 Intel Core i5 – 12400	M1 (8-core CPU, 7-core GPU) M2 (8-core CPU, 8-core GPU) M3 (8-core CPU, 8-core GPU)
RAM	16 GB		16 GB
Graphic	NVIDIA GeForce RTX 4050 (VRAM 6 GB or higher) NVIDIA GeForce RTX 3060 (VRAM 6 GB or higher) NVIDIA RTX A2000 (VRAM 6 GB or higher) *AMD Radeon is not supported.		-
OS	Windows 10 64-bit Windows 11 (recommended for Intel Core 12th Gen or newer processors)		Monterey 12 Ventura 13

##### Recommend System Requirements

	Windows		macOS
	Laptop	Desktop	Laptop/Desktop
CPU	Intel Core Ultra 7 155H Intel Core i7 – 13700H Intel Core i7 – 12700H	Intel Core i7 – 13700K Intel Core i7 – 12700K	M1 Pro (10-core CPU, 16-core GPU) M2 (8-core CPU, 10-core GPU) M2 Pro (10-core CPU, 16-core GPU) M3 (8-core CPU, 10-core GPU) M3 Pro (11-core CPU, 14-core GPU)
RAM	32 GB		24 GB
Graphic	NVIDIA GeForce RTX 4060 (VRAM 8 GB or higher) NVIDIA GeForce RTX 3070 (VRAM 8 GB or higher) NVIDIA RTX A3000 (VRAM 8 GB or higher) *AMD Radeon is not supported.		-
OS	Windows 10 64-bit Windows 11 (recommended for Intel Core 12th Gen or newer processors)		Monterey 12 Ventura 13

For accurate and up-to-date system requirements, please visit [www.meditlink.com](http://www.meditlink.com).

Use PC and monitor certified IEC 60950, IEC 55032, IEC 55024.

The device may not work when using cables other than the USB 3.0 cable provided by Medit. Medit is not responsible for any problems caused by cables other than the USB 3.0 cable provided by Medit. Be sure to use only the USB 3.0 cable included in the package.

#### 2.2.2 Software Installation Guide

- ① Connect the enclosed USB flash drive to a PC.
- ② Run the installation file.
- ③ Select the setup language and click "Next."
- ④ Choose the installation path.
- ⑤ Read the "License Agreement" carefully, check "I agree to the License terms and conditions," and then click "Install."
- ⑥ The installation process may take several minutes. Please do not shut down the PC until the installation is complete.
- ⑦ After the installation is complete, restart the PC to ensure optimal program operation.

The installation will not be processed while the i700 wireless system is connected to a PC. Please turn off the scanner before starting the installation.

#### 2.2.3 Medit Scan for Clinics User Guide

Please refer to the User Guide of Medit Scan for Clinics: Medit Scan for Clinics > Menu > User Guide.

### 3. Maintenance

#### CAUTION

- Equipment maintenance should only be carried out by a Medit employee or Medit-Certified company or personnel.
- In general, users are not required to perform maintenance work on the i700 wireless system besides calibration, cleaning, and sterilization. Preventive inspections and other regular maintenance are not required.

#### 3.1 Calibration

Periodic calibration is required to produce precise 3D models. You should perform calibration when:

- The quality of the 3D model is not reliable or accurate when compared to previous results.
  - Environmental conditions such as temperature have changed.
  - The calibration period has expired.
- You can set the calibration period in the Menu > Settings > Calibration Period (Days).

The calibration panel is a delicate component. Do not touch the panel directly. Check the calibration panel if the calibration process is not performed properly. If the calibration panel is contaminated, please contact your service provider.

We recommend performing a calibration periodically. You can set the calibration period via Menu > Settings > Calibration Period (Days). The default calibration period is 14 days.

##### 3.1.1 How to Calibrate i700 wireless

- ① Turn on the i700 wireless and launch the Medit Scan for Clinics.
- ② Run the Calibration Wizard from Menu > Settings > Calibration.
- ③ Prepare the calibration tool and the i700 wireless handpiece.
- ④ Turn the dial of the calibration tool to position **1**.
- ⑤ Put the i700 wireless handpiece into the calibration tool.
- ⑥ Click "Next" to start the calibration process.
- ⑦ When the calibration tool is mounted properly in the correct position **1**, the system will automatically acquire the data.
- ⑧ When data acquisition is complete at position **1**, turn the dial to the next position.
- ⑨ Repeat the steps for positions **2** - **8** and the **LAST** position.
- ⑩ When data acquisition is complete at the **LAST** position, the system will automatically calculate and show the calibration results.

##### 3.1.2 Auto Calibration Tool (Sold Separately)

The i700 wireless auto calibration tool accessory can be purchased separately. This handy auto calibration tool will automatically perform calibration by calibrating the i700 wireless handpiece without having to turn the calibration dial. Please refer to Medit Scan for Clinics for details.

#### 3.2 Cleaning, Disinfection, Sterilization Procedure

##### 3.2.1 Reusable Tip

The reusable tip is the part that is positioned in the mouth of patients during scanning, and is reusable for a limited number of times. The tip needs to be cleaned and sterilized between patient uses to avoid cross-contamination.

##### Cleaning & Disinfection

- Prepare a cleaning solution.
  - » Dilute a general neutral detergent at a 1:100 ratio before use.
- Clean the reusable tips with a cleaning solution and a brush.
  - » Ensure the mirror of the tip is completely clean and stain-free. If the mirror appears stained or foggy, repeat the cleaning process.

#### CAUTION

- » The reusable tip has a complex structure, and automatic cleaning may not achieve complete cleaning; therefore, do not clean the reusable tip in an automatic washer.
- Rinse the reusable tips three times with purified water.
- Remove moisture with a paper towel and allow them to air-dry completely at room temperature for at least 80 minutes.
- Disinfect the reusable tips using a disinfectant containing 15% or less Isopropyl Alcohol (IPA) for 1 minute. Then, ensure they are thoroughly dried at room temperature for at least 5 minutes.
  - » Before using a disinfectant product, please refer to the product user manual for proper use.
  - » You can find the list of recommended disinfectants in the Medit Help Center at <http://support.medit.com/hc>.

##### Sterilization

- The tip should be cleaned manually using a disinfecting solution.
- After cleaning and disinfecting, inspect the mirror inside the tip to ensure there aren't any stains or smudges.
- Repeat the cleaning and disinfection process if necessary. Carefully dry the mirror with a paper towel.
- Insert the tip into a paper sterilization pouch and seal it, making sure that it is airtight. Use either a self-adhesive or heat-sealed pouch.
- Sterilize the wrapped tip in an autoclave with the following conditions:
  - » Sterilize in a gravity autoclave at 135°C (275°F) for 10 minutes and dry for 30 minutes.
  - » Sterilize in a pre-vacuum autoclave at 134°C (273.2°F) for 4 minutes and dry for 20 minutes.
- Use an autoclave program that dries the wrapped tip before opening the autoclave.
- Scanner tips can be re-sterilized up to 150 times. After reaching this limit, they must be disposed of according to the guidelines in the disposal section.
- Autoclave times and temperatures may vary depending on the autoclave type and manufacturer. For this reason, it may not be able to meet the maximum number of times. Please refer to the user manual of the autoclave manufacturer you are using to determine whether the required conditions are met.

#### CAUTION

- The mirror of the scanner tip is a delicate optical component that should be handled with care to ensure optimal scan quality. Be careful not to scratch or smudge it as any damage or blemishes may affect the data acquired.
- Make sure to always wrap the tip before autoclaving. If you autoclave an exposed tip, it will cause stains on the mirror, which cannot be removed. Check the autoclave manual for more information.
- Tips that have been cleaned, disinfected, and sterilized must remain sterile until they are used on the patient.
- Medit is not responsible for any damage, including tip distortion, caused by cleaning, disinfection, or sterilization procedures that do not comply with the guidelines outlined above.

### 3.2.2 Mirror

The presence of impurities or smudges on the tip mirror may lead to poor scan quality and an overall poor scanning experience. In such situations, clean the mirror following the steps below:

- ① Disconnect the scanner tip from the i700 wireless handpiece.
- ② Pour alcohol on a clean cloth or cotton-tipped swab and wipe the mirror. Make sure to use alcohol that is free of impurities or it may stain the mirror. You can use either ethanol or propanol (ethyl-/propyl alcohol).
- ③ Wipe the mirror dry using a dry, lint-free cloth.
- ④ Make sure the mirror is free of dust and fibers. Repeat the cleaning process as necessary.

### 3.2.3 Handpiece

After treatment, clean and disinfect all other surfaces of the i700 wireless handpiece except for the scanner front (optical window) and end (air vent hole). Cleaning and disinfecting must be done with the device turned off. Use the device only after it is completely dry.

The recommended cleaning and disinfecting solution is denatured alcohol (ethyl alcohol or ethanol) – typically 60 – 70% Alc/Vol.

The general cleaning and disinfecting procedures are as follows:

- ① Turn off the device using the power button.
- ② Unplug all the cables from the wireless hub.
- ③ Clean the filter on the front end of the i700 wireless handpiece.
  - » If alcohol is poured directly into the filter, it may seep inside the i700 wireless handpiece and cause a malfunction.
  - » Do not clean the filter by pouring alcohol or cleaning solution directly into the filter.  
The filter must be gently wiped with a cotton or soft cloth moistened with alcohol. Do not wipe by hand or apply excessive force.
  - » Medit is not responsible for any damage or malfunction that occurs during cleaning that does not follow the guidelines above.
- ④ Put the cover on the front of the i700 wireless handpiece after cleaning the filter.
- ⑤ Pour the disinfectant onto a soft, lint-free, and non-abrasive cloth.
- ⑥ Wipe the scanner surface with the cloth.
- ⑦ Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

#### CAUTION

- Do not clean the i700 wireless handpiece when the device is turned on as the fluid may enter the scanner and cause malfunction.
- Use the device after it is completely dry.
- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

### 3.2.4 Other Components

- Pour the cleaning and disinfecting solution onto a soft, lint-free, and non-abrasive cloth.
- Wipe the component surface with the cloth.
- Dry the surface with a clean, dry, lint-free, and non-abrasive cloth.

#### CAUTION

- Chemical cracks may appear if improper cleaning and disinfecting solutions are used during cleaning.

### 3.3 Disposal

#### CAUTION

- The scanner tip must be sterilized before disposal. Sterilize the tip as described in the section “3.2.1 Reusable Tip & Small Tip – Sterilization.”
- Dispose of the scanner tip as you would any other clinical waste.
- Other components are designed to conform with the following directives:  
RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)  
WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

### 3.4 Battery Storage

- Put it in a pack or a box and store indoors in a cool environment out of direct sunlight.
- Store the battery in a dry place with an ambient temperature of from  $-20^{\circ}\text{C}$  to  $+30^{\circ}\text{C}$  ( $-4^{\circ}\text{F}$  to  $+86^{\circ}\text{F}$ ).
- If not used for a long time, the battery self-discharge may accelerate and slip into sleep mode.
- To minimize the effect of deactivation, store the battery pack in between  $+10^{\circ}\text{C}$  -  $+30^{\circ}\text{C}$  ( $+50^{\circ}\text{F}$  to  $+86^{\circ}\text{F}$ ).
- When charging for the first time after long-term storage, the battery capacity may decrease due to pack deactivation. Restore the battery through several full charge/discharge cycles.
- If the battery is stored for more than 6 months, it should be charged at least once every 6 months to prevent the decrease in shelf life due to self-discharge.

#### CAUTIONS: Battery Safety Specifications

Safety Specifications		
Overcharge	Sense Voltage	$4.225\text{ V} \pm 0.020$
	Steady Voltage	$4.025\text{ V} \pm 0.03$
	Sense Delay	$1.0\text{ s} \pm 0.2$
Over-Discharge	Sense Voltage	$2.50\text{ V} \pm 0.035$
	Steady Voltage	$2.90\text{ V} \pm 0.50$
	Sense Delay	$64\text{ ms} \pm 12.8$
Overcurrent	Sense (Charge)	$10.0\text{ A} +5.0 / -4.0$
	Sense Delay	$8.0\text{ ms} \pm 1.6$
	Sense (Discharge)	$10.0\text{ A} +4.4 / -3.8$
	Sense Delay	$8.0\text{ ms} \pm 1.6$
Current consumption at operating mode		Max. $150.0\text{ }\mu\text{A}$

The possible safety specifications are determined by the propulsion control module (PCM) in the parts list.

### 3.5 The Battery Usage Precautions and Disposal Guide

#### CAUTION

- Be sure to fully understand how to replace the battery before use.
- Use a charger suitable for the specified voltage and current.
- Do not attempt to reverse the charge. Reverse charging may increase the gas pressure inside the battery and cause the battery to leak.
- Do not attempt to recharge a fully charged battery. Repeated overcharging could result in battery performance degradation and overheating.
- Charging efficiency decreases at temperatures above  $+40^{\circ}\text{C}$  ( $+104^{\circ}\text{F}$ ).
- Do not short-circuit the positive (+) and negative (-) terminals with metallic objects such as metal wires, necklaces, or chains.
- To avoid malfunction or damage, do not drop or throw the battery.
- Do not deform a battery by excessive pressure.
- Do not directly solder anything onto the battery pack.
- Do not allow children to change batteries without adult supervision.
- Do not dispose of or throw the battery into a fire. Heat can cause a battery explosion and a fire.
- Separate them from each other when disposing of the secondary batteries with different electrochemical systems.
- Dispose of the battery by discharging it to prevent heat from short-circuiting.
- Battery disposal methods may vary by country and region. Dispose of used batteries under local laws and regulations.

### 3.6 Updates in Medit Scan for Clinics

Medit Scan for Clinics automatically checks for updates when the software is in operation. If a new version of the software is released, the system will automatically download it.

## 4. Safety Guide

Please adhere to all the safety procedures as detailed in this User Guide to prevent human injury and equipment damage. This document uses the words WARNING and CAUTION when highlighting precautionary messages.

Carefully read and understand the guidelines, including all preventive messages as prefaced by the words WARNING and CAUTION. To avoid bodily injury or equipment damage, make sure to adhere strictly to the safety guidelines. All instructions and precautions as specified in the Safety Guide must be observed to ensure the system's proper functionality and personal safety.

The i700 wireless system should only be operated by dental professionals and technicians who are trained to use the system. Using the i700 wireless system for any purpose other than its intended usage as outlined in the section "1.1 Intended Use" may result in injury or damage to the equipment. Please handle the i700 wireless system according to the guidelines in the safety guide.

Any serious incident involving the device should be reported to the manufacturer and the competent authority of the member state in which the user and patients are located.

### 4.1 System Basics

The i700 wireless system is a high-precision optical medical device. Acquaint all the following safety and operating instructions before the installation, usage, and operation of i700 wireless.

#### CAUTION

- The USB 3.0 cable connected to the wireless hub is the same as a standard USB cable connector. However, the device may not operate normally if a standard 3.0 USB cable is used with i700 wireless.
- The wireless hub is designed specifically for the i700 wireless and should not be used with any other device.
- This device has been assessed for conformity for use in business environments. When used in residential environments, there is a risk of radio wave interference.
- If the product has been stored in a cold environment, give it time to adjust to the temperature of the environment before use. If used immediately, condensation may occur, which may damage the electronic parts inside the unit.
- Ensure that all components provided are free from physical damage. Safety cannot be guaranteed if there is any physical damage to the unit.
- Before using the system, check that there are no issues such as physical damage or loose parts. If there is any visible damage, do not use the product and contact the manufacturer or your local representative.
- Check the i700 wireless handpiece and its accessories for any sharp edges.
- When not in use, the i700 wireless system should be kept mounted on a desk stand or a wall mount stand.
- Do not install the desk stand on an inclined surface.
- Do not place any object on the i700 wireless system.
- Do not place the i700 wireless system on any heated or wet surface.
- Do not block the air vents located at the rear of the i700 wireless system.
- If the equipment overheats, the i700 wireless system may malfunction or stop working.
- The battery of the i700 wireless can only be compatible with the i700 wireless system.
- Do not touch the charging terminals of the rechargeable battery with your hands or other instruments.
- If the charging terminal of the rechargeable battery is damaged, do not use it and contact the manufacturer or regional manager.
- If the shape of the rechargeable battery has been deformed by dropping or physical damage, never use it, and contact the manufacturer or regional manager.
- Extension cable battery is not designed for charging. Do not attempt to charge it with the battery charger.
- Use only the extension cable battery provided by the manufacturer.
- Do not spill any liquid on the i700 wireless system.
- The i700 wireless handpiece and other included components are made of electronic components. Do not allow any kind of liquid or foreign objects to enter.
- Do not pull or bend the cable connected to the i700 wireless system.
- Carefully arrange all the cables so that you or your patient do not trip or get caught in the cables. Any pulling tension on the cables may cause damage to the i700 wireless system.
- Always place the plug of the power cord of the i700 wireless system in an easily accessible location.
- Always keep an eye on the product and your patient while using the product to check for abnormalities.
- Proceed with calibration, cleaning, disinfection, and sterilization in accordance with the contents of the user guide.
- If you drop the i700 wireless tip on the floor, do not attempt to reuse it.
- Discard the tip immediately as there is a risk that the mirror attached to the tip may have been dislodged.
- Due to its fragile nature, the i700 wireless tips should be handled with care. To prevent damage to the tip and its internal mirror, be careful to avoid contact with the patient's teeth or restorations.
- If the i700 wireless system is dropped on the floor or if the unit is impacted, it must be calibrated before use.
- If the instrument is unable to connect to the software, consult the manufacturer or authorized resellers.
- If the equipment fails to operate normally, such as having issues with accuracy, stop using the product, and contact the manufacturer or authorized resellers.
- Install and use only approved programs to ensure the proper functionality of the i700 wireless system.
- In the event of a severe accident involving the i700 wireless system, notify the manufacturer and report it to the competent national authority of the country where the user and patient reside.
- If the PC with the software installed does not have security software or if there is a risk of malicious code intrusion into the network, the PC may be breached with malware (malicious software such as viruses or worms that damage your computer).
- The software for this product must be used in compliance with medical and personal information protection laws.
- This high-precision 3D scanner is particularly sensitive to impacts. Handle with care to prevent performance degradation from shocks, including accidental drops.

### 4.2 Proper Training

#### WARNING

Before using your i700 wireless system on patients:

- You should have been trained to use the system or read and fully understand this User Guide.
- You should be familiar with the safe use of the i700 wireless system, as detailed in this User Guide.
- Before use or after changing any settings, the user should check that the live image is displayed properly in the camera preview window of the program.

### 4.3 In Case of Equipment Failure

#### WARNING

If your i700 wireless system is not working properly or if you suspect that there is a problem with the equipment:

- Remove the device from the patient's mouth and discontinue use immediately.
- Disconnect the device from the PC and check for errors.
- Remove the rechargeable battery from the i700 wireless system.
- Contact the manufacturer or authorized resellers.
- Modifications to the i700 wireless system are prohibited by law as they may compromise the safety of the user, patient, or a third party.

### 4.4 Hygiene

#### WARNING

For clean working conditions and patient safety, ALWAYS wear clean surgical gloves when:

- Handling and replacing the tip.
- Using the i700 wireless system on patients.
- Touching the i700 wireless system.

#### WARNING

The i700 wireless system and its optical window should be always kept clean. Before using the i700 wireless system on a patient, be sure to:

- Sterilize the i700 wireless system as described in the section "3.2 Cleaning, Disinfection, Sterilization Procedure."
- Use a sterilized tip.

### 4.5 Electrical Safety

#### WARNING

- The i700 wireless system is a Class I device.
- The i700 wireless system includes the i700 wireless handpiece, wireless hub, battery charger, and rechargeable battery altogether.
- To prevent electric shock, the i700 wireless system must only be connected to a power source with a protective earth connection. If you cannot insert the i700 wireless-supplied plug into the main outlet, contact a qualified electrician to replace the plug or outlet. Do not try to circumvent these safety guidelines.
- Do not use a grounding-type plug connected to the i700 wireless system for any other purpose than its intended use.
- The i700 wireless system only uses RF energy internally.
- The amount of RF radiation is low and does not interfere with surrounding electromagnetic radiation.
- There is a risk of electric shock if you attempt to access the inside of the i700 wireless system.
- Only qualified service personnel should access the system.
- Do not connect the i700 wireless system to a regular power strip or extension cord, as these connections are not as safe as grounded outlets. Failure to adhere to these safety guidelines may result in the following hazards:
  - » All connected equipment's total short circuit current may exceed the limit specified in EN/IEC 60601-1.
  - » The impedance of the ground connection may exceed the limit specified in EN/IEC 60601-1.
- Do not place liquids such as beverages near the i700 wireless system and avoid spilling any liquid on the system.
- Never spill liquid of any kind on the i700 wireless system.
- Condensation due to changes in temperature or humidity can cause moisture buildup inside the i700 wireless system, which may damage the system. Before connecting the i700 wireless system to a power supply, be sure to keep the i700 wireless system at room temperature for at least two hours to prevent condensation. If condensation is visible on the product surface, the i700 wireless should be left at room temperature for more than 8 hours.
- You should only disconnect the i700 wireless system from the power supply via its power cord or rechargeable battery.
- When disconnecting the power cord or the rechargeable battery, hold the surface to remove it.
- Before disconnecting, make sure to turn off the power on the device using the power switch on the handpiece.
- The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If used in a residential environment (for which CISPR 11 Class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services.
- Use only batteries provided for use with the i700 wireless. Other batteries could damage the i700 wireless system.
- Avoid pulling on the communication cables, power cables, etc. used with the i700 wireless system.
- Use only medical adapters provided for use with the i700 wireless. Other adapters could damage the i700 wireless system.
- Do not touch the connectors of the device and the patient simultaneously.

### 4.6 Eye Safety

#### WARNING

- The i700 wireless system projects a bright light from its tip during scanning.
- The bright light projected from the tip of the i700 wireless is not harmful to the eyes. However, you should not look directly at the bright light nor aim the light beam into the eyes of others. Generally, intense light sources can cause eyes to become brittle and the likelihood of secondary exposure is high. As with other intense light source exposure, you may experience a temporary reduction in visual acuity, pain, discomfort, or visual impairment, all of which increase the risk of secondary accidents.
- There is an LED that emits UV-C wavelengths inside the i700 wireless handpiece. It is irradiated only inside the i700 wireless handpiece and does not go outside. The blue light visible inside the i700 wireless handpiece is for guidance, not UV-C light. It is harmless to the human body.
- The UV-C LED operates with a wavelength of 270 – 285 nm.
- Disclaimer for Risks Involving Patients with Epilepsy  
The Medit i700 wireless should not be used on patients that have been diagnosed with epilepsy due to the risk of seizures and injury. For the same reason, dental staff who have been diagnosed with epilepsy should not operate the Medit i700 wireless.

#### 4.7 Explosion Hazards

##### WARNING

- The i700 wireless system is not designed to be used near flammable liquids, gases, or in environments with high oxygen concentrations.
- There is a risk of explosion if you use the i700 wireless system near flammable anesthetics.
- The rechargeable battery used with the i700 wireless is designed with safety devices included.
- The rechargeable battery shall not be exposed to excessive heat such as sunshine or the like.
- Failure to observe this could result in the battery explosion. Please be careful with the storage and maintenance of the battery.
- The battery charger is designed to adjust the current supply after charging is complete. But if it is not going to be used for a long time, disconnect the power to the battery charger or remove the fully charged battery from the battery charger.

#### 4.8 Pacemaker and ICD Interference Risk

##### WARNING

- Implantable Cardioverter Defibrillators (ICDs) and pacemakers may have interference due to some devices.
- Maintain a moderate distance from the patient's ICD or pacemaker when using the i700 wireless system.
- For more information on peripherals used with i700 wireless, check the respective manufacturer's manuals.

#### 4.9 Cybersecurity Safety

- If a cybersecurity incident occurs, stop using the scanner and software immediately. Power off the scanner and log out of the software.
- Promptly report the incident to our support team through email, phone, or other available means of contact. Please refer to the last page of the User Guide for contact information.
- When reporting an incident, please provide as much information as possible, including the time of occurrence and any unusual behavior you noticed. This information will assist us in resolving the issue quickly.

## 5. Electro-Magnetic Compatibility Information

### 5.1 Electromagnetic Emissions

The i700 wireless system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Emission		
Emission Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The i700 wireless uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class A	The i700 wireless is suitable for use in all establishments. This includes domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations / Flicker Emissions	Complies	

##### WARNING

This i700 wireless system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the i700 wireless or shielding the location.

### 5.2 Electromagnetic Immunity

#### • Guidance 1

The i700 wireless system is intended for use in the electromagnetic environment as specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be made of wood, concrete, or ceramic tiles. If floors are covered with a synthetic material, relative humidity of at least 30% is recommended.
Electrical Fast Transient / Burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	The mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0.5 kV, ± 1 kV differential mode	± 0.5 kV, ± 1 kV differential mode	The mains power quality should be that of a typical commercial or hospital environment.
	± 0.5 kV, ± 1 kV, ± 2 kV common mode	± 0.5 kV, ± 1 kV, ± 2 kV common mode	

Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Input Lines IEC61000-4-11	0% Ur (100% dip in Ur) for 0.5/1 cycles 70% Ur (30% dip in Ur) for 25/30 cycles 0% Ur (100% dip in Ur) for 250/300 cycles	0% Ur (100% dip in Ur) for 0.5/1 cycles 70% Ur (30% dip in Ur) for 25/30 cycles 0% Ur (100% dip in Ur) for 250/300 cycles	The mains power quality should be that of a typical commercial or hospital environment. If the user of the i700 Wireless system requires continued operation during power mains interruptions, it is recommended that the i700 wireless system be powered from an uninterruptible power supply or a battery.
--------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Power Frequency Magnetic Fields (50/60Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a location in a typical commercial or hospital environment.
------------------------------------------------------------	--------	--------	-----------------------------------------------------------------------------------------------------------------------------------

Proximity Magnetic Fields in the Frequency Range 9 kHz to 13.56 MHz Immunity IEC 61000-4-39	8 A/m 30 kHz CW modulation	8 A/m 30 kHz CW modulation	Resistance to magnetic fields was tested and applied only to surfaces of enclosures or accessories accessible during intended use.
	65 A/m 134.2 kHz PM 2.1 kHz	65 A/m 134.2 kHz PM 2.1 kHz	
	7.5 A/m 13.56 MHz PM 50 kHz	7.5 A/m 13.56 MHz PM 50 kHz	

NOTE: Ur is the main voltage (AC) prior to the application of the test level.

• **Guidance 2**

Recommended Separation Distances Between Portable and Mobile Communication Equipment and the i700 wireless			
Rated Maximum Output Power of the Transmitter [W]	Separation Distance According to the Frequency of Transmitter [M]		
	IEC 60601-1-2: 2014		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 2.7 GHz $d = 2.0 \sqrt{P}$	
0.01	0.12	0.20	
0.1	0.38	0.63	
1	1.2	2.0	
10	3.8	6.3	
100	12	20	

For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

• **Guidance 3**

The i700 wireless system is intended for use in the electromagnetic environment specified below. The customer or the user of the i700 wireless system should ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz Outside ISM Bands amateur	3 Vrms	Portable and mobile RF communications equipment should not be used closer to any part of the Ultrasound System, including cables, than the recommended separation distance. This is calculated using the equation applicable to the frequency of the transmitter.  <b>Recommended Separation Distance (d):</b> $d = 1.2 \sqrt{P}$ <b>IEC 60601-1-2:2007</b> $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 80 MHz to 2.5 GHz <b>IEC 60601-1-2:2014</b> $d = 2.0 \sqrt{P}$ 80 MHz to 2.7 GHz
	6 Vrms 150 kHz to 80 MHz In ISM Bands amateur	6 Vrms	
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:  

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

NOTE 3: The ISM (Industrial, Scientific, and Medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

• **Guidance 4**

The i700 wireless system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. Portable RF communications equipment should be used no closer than 30 cm (12 inches) to any part of the i700 wireless system. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity					
Immunity Test	Band <sup>1)</sup>	Service <sup>1)</sup>	Modulation	IEC 60601 Test Level	Compliance Level
Proximity Fields from RF Wireless Communications IEC 61000-4-3	380 – 390 MHz	TETRA 400	Pulse Modulation 18 Hz	27 V/m	27 V/m
	430 – 470 MHz	GMRS 460; FRS 460	FM ±5 kHz Deviation 1 kHz sine	28 V/m	28 V/m
	704 – 787 MHz	LTE Band 13, 17	Pulse Modulation 217 Hz	9 V/m	9 V/m
	800 – 960 MHz	GSM 800;900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	Pulse Modulation 18 Hz	28 V/m	28 V/m
	1700 – 1990 MHz	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation 217 Hz	28 V/m	28 V/m
	2400 – 2570 MHz	Bluetooth; WLAN 802.11b/g/n; RFID 2450; LTE Band 7	Pulse Modulation 217 Hz	28 V/m	28 V/m
	5100 – 5800 MHz	WLAN 802.11a/n	Pulse Modulation 217 Hz	9 V/m	9 V/m

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

**WARNING**

- Use of the i700 wireless adjacent to or on other equipment must be avoided as it may result in improper operation. If this use is necessary, it is advisable that this and the other equipment be observed to verify that they are operating normally.
- The use of accessories, transducers, and cables other than those specified or provided by the Medit of the i700 wireless could result in high electromagnetic emissions or reduced electromagnetic immunity from this equipment and result in improper operation.

<sup>1</sup> For some services, only the uplink frequencies are included.

## 6. Wireless Compliance Information

### 6.1 FCC Compliance Statement

Federal Communications Commission Interference Statement

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**FCC Caution:** Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

#### IMPORTANT NOTE:

**FCC Radiation Exposure Statement:**

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with a minimum distance of 20 cm (8 inches) between the radiator and your body.

### 6.2 IC Compliance Statement

This Class A digital apparatus complies with Canadian ICES-003.

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Any changes or modifications not expressly approved by the manufacturer could void the user's authority to operate this equipment.

This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

The device could automatically discontinue transmission in case of the absence of information to transmit or operational failure. Note that this is not intended to prohibit transmission of control or signaling information or the use of repetitive codes where required by the technology.

#### IMPORTANT NOTE:

**IC Radiation Exposure Statement**

This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment.

This equipment should be installed and operated with a minimum distance of 20 cm (8 inches) between the radiator and your body.

#### Transmit Antenna Notice

This radio transmitter has been approved by Innovation, Science and Economic Development Canada to operate with the antenna types listed below, with the maximum permissible gain indicated. Antenna types not included in this list that have a gain greater than the maximum gain indicated for any type listed are strictly prohibited from use with this device.

#### Antenna List

Model	Type	Maximum Gain (dBi)
Si16310	Patch Array Antenna	18 dBi
2450AT07A0100	Passive Antenna	1 dBi

### 6.3 KC Compliance Statement

This device has been assessed for conformity for use in business environments. When used in residential environments, there is a risk of radio wave interference.

### 6.4 TELEC (Japan) Compliance Statement

R 209-J00306  
R 209-J00282  
R 020-180117  
R 020-180116

## 7. Specifications

Model Name		MD-IS0300
Trade Name	i700 wireless	
Packing Unit	1 set	
Classifications for Protection Against Electric Shock	Class I, Type BF Applied Parts	
* This product is a medical device.		
Handpiece		
Dimension	312.7 x 43.9 x 47.4 mm (W x L x H)	
Weight	254.5 g	
Rating	4 V <sup>~</sup> , 4 A	
DC Adapter		
Model Name	ATM065T-P120	
Input Voltage	Universal 100 – 240 Vac, 50 – 60 Hz input, without any slide switch	
Output	12 V <sup>~</sup> , 5 A	
Case Dimension	119 x 60 x 36 mm (W x L x H)	
EMI	CE / FCC Class B, Conduction & Radiation met	
Protection	OVP (Over Voltage Protection)	
	SCP (Short Circuit Protection)	
	OCP (Over Current Protection)	
Protection Against Electric Shock	Class I	
Mode of operation	Continuous	
Rechargeable Battery		
Model Name	MD-IS0300REB	
Type	Lithium-ion	
Output	3.6 Vdc, 11.16 Wh	
Dimension	21.4 x 73.4 mm	
Weight	60 g	
Capacity	3,100 mAh	
Wireless Hub		
Input	12 V <sup>~</sup> , 5 A 9 V / 5 V <sup>~</sup> , 3 A	
Dimension	100 x 94.8 x 44.4 mm (W x L x H)	
Weight	181 g	
Battery Charger		
Input / Output	12 V <sup>~</sup> , 5 A	
Dimension	44.7 x 100 mm (H x Ø)	
Weight	247 g	
Calibration Tool		
Dimension	123.8 x 54 mm (H x Ø)	
Weight	220 g	
Auto Calibration Tool (*Sold Separately)		
Dimension	168.7 x 92.1 x 48.2 mm	
Weight	492 g	
Output	3.6 Vdc, 11.16 Wh (MD-IS0300ECB)	
Wireless Module		
60 GHz	Frequency Bands	HRP: 60.48 – 62.64 GHz MRP: 60.48 – 62.64 GHz LRP: 60.16 – 62.96 GHz
	Modulation Type	BPSK
	e.i.r.p.	HRP: 24.2 dBm MRP: 24.0 dBm LRP: 14.6 dBm
	Antenna Gain	18.0 dBi

2.4 GHz (Bluetooth LE)	Frequency Bands	2402 – 2480 MHz
	Channels	40
	Channel Bandwidth	1 MHz 2 MHz
	e.i.r.p.	9.8 dBm A-variant: 19.7 dBm N-variant: 19.8 dBm
	Modulation	GFSK
	Antenna Gain	A-variant: 1 dBi N-variant: 2.14 dBi
<b>Operating, Storage, and Transport Conditions</b>		
Operating Condition	Temperature	18 – 28°C (64.4 – 82.4°F)
	Humidity	20 – 75% relative humidity (non-condensing)
	Air Pressure	800 – 1,060 hPa
Storage Condition	Temperature	-5 – 45°C (23 – 113°F)
	Humidity	20 – 80% relative humidity (non-condensing)
	Air Pressure	800 – 1,100 hPa
Transport Condition	Temperature	-5 – 45°C (23 – 113°F)
	Humidity	20 – 80% relative humidity (non-condensing)
	Air Pressure	620 – 1,200 hPa

EC REP

**Meditrial Srl**  
Via Po 9 00198, Rome Italy

CH REP

**Meditrial Europe Ltd**  
Bahnhofstrasse 23 6300 Zug, Switzerland

UK REP

**Meditrial UK Ltd**  
19 The Circle, Queen Elizabeth Street, London, SE1 2JE, United Kingdom

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