User Manual

i500

Revision 6 (January 2019)

CONTENTS

02	English
26	한국어
46	Deutsch
70	español
94	français
118	italiano
142	Nederlands
166	polski
191	português
214	română
238	Ελληνικά
262	العربية
286	中文(简体)
306	中文(繁體)
326	日本語

English

1	About this guide	4
2	Introduction and overview	4
2.1	Intended use	4
2.2	Indication for use	4
2.3	Contraindications	5
2.4	Qualifications of the operating user	5
2.5	Symbol	5
2.6	i500 Components Overview	6
2.7	Setting up the i500	7
2.7.1	Basic settings of i500	7
2.7.2	Placing on Desktop Cradle	8
2.7.3	Installation of Wall Mount Holder	8
3	Image Acquisition Software Overview	8
3.1	Introduction	8
3.2	Installation	9
3.2.1	System Requirement	9
3.2.2	Installation Guide	9
4	Maintenance	10
4.1	Calibration	10
4.2	Cleaning and sterilization procedure	11
4.2.1	Reusable tip	11
4.2.2	Mirror	12
4.2.3	Handpiece	12
4.2.4	Other components	13
4.3	Disposal	13
4.4	Updates to Image Acquisition Software	13
5	Safety Guide	14
5.1	System Basics	14
5.2	Safety and Warning	15
5.2.1	Preventive Inspection Before Use of the System	15
5.2.2	Modification of the System	15
5.2.3	Approved Software Only	15
5.2.4	Proper Training	15
5.2.5	In Case of Equipment Failure	16
5.3	Mechanical Hazards	16
5.4	Explosion Hazards	16
5.5	Electrical Safety	16
5.6	Eye Safety	18
5.7	Pacemaker and ICD Interference Risk	18
5.8	Hygiene	18
5.9	Safety when scanning with i500	19
5.10	Protection from overheating	19
5.11	Technical Maintenance	19
6	Electro – Magnetic Compatibility Information	20
6.1	Electro-Magnetic Emissions	20
6.2	Electro-Magnetic Immunity	20
7	Specification	24

1 About this guide

Convention in this guide

This user guide uses various symbols to highlight important information so as to ensure correct usage, prevent injury to the user and others, and prevent property damage. The meanings of the symbols used are described below.

WARNING

The WARNING symbol indicates information that, if ignored, could result in medium risk of personal injury.

CAUTION

The CAUTION symbol indicates safety information that, if ignored, could result in slight risk of personal injury, property damage or damage to the system.

TIPS

The TIPS symbol indicates hints, tips and additional information for optimal operation of the system.

2 Introduction and overview

2.1 Intended use

The i500 system is a dental 3D scanner intended to be used to digitally record topographical characteristics of teeth and surrounding tissues. The i500 system produces 3D scans for use in computer-assisted design and manufacturing of dental restorations

2.2 Indication for use

The i500 system should be used on patients who require 3D scanning for dental treatments such as:

- Single custom abutment
- Inlays & Onlays
- Single Crown
- Veneer
- 3 Unit Implant Bridge
- Up to 5 Unit Bridge
- Orthodontics
- Implant Guide
- Diagnosis Model

2.3 Contraindications


The device is not intended to be used to create images of the internal structure of teeth or the supporting skeletal structure.

It is not intended to be used for cases with more than four (4) subsequent edentulous tooth positions.


2.4 Qualifications of the operating user

This device is designed for use by persons with professional knowledge in dentistry and dental laboratory technology. The user of this device is solely responsible for determining whether or not this device is suitable for a particular patient case and circumstances. The user is solely responsible for the accuracy, completeness and adequacy of all data entered into this device and the provided software. The user has to check the correctness and accuracy of the results and to assess each individual case. The i500 system must be used in accordance with its accompanying user guide. The user is not allowed to modify the i500 system. Improper use or handling of the i500 system will void its warranty, if any. If you require additional information on the proper use of the i500 system, please contact your local distributor.

2.5 Symbol

No	Symbol	Description
01		The serial number of the object
02		Date of manufacture
03		Manufacturer
04		Caution
05		Warning
06		Instruction for User Manual
07		The official mark of Europe Certificate
08		Authorized representative in the European Community
09		Type of applied part
10		WEEE Mark

2 Introduction and overview


11		Prescription use (U.S)
12		MET mark
13		AC
14		DC
15		Protective Earth (ground)

2.6 i500 Components Overview

No	Item	Qty	Appearance
01	i500 Handpiece + Power Hub	1ea	
02	i500 Handpiece Cover	1ea	
03	Reusable Tip	4ea	
04	Calibration Tool	1ea	
05	Desktop Cradle	1ea	
06	Wall Mount Holder	1ea	
07	USB 3.0 Cable	1ea	
08	Medical Adapter + Power Cord	1ea	
09	USB Memory (Pre-loaded with image acquisition software)	1ea	
10	User Guide	1ea	

2.7 Setting up the i500

2.7.1 Basic settings of i500

- ① Connect the USB 3.0 Cable to the Power Hub

- ② Connect Medical Adapter to the Power Hub

- ③ Connect Power Cord to the Medical Adapter

- ④ Connect Power Cord to a Power Source

- ⑤ Connect the USB 3.0 Cable to the PC


2 Introduction and overview

 Turn on the i500

Press the power button of i500 scanner

Wait until the USB connection indicator turns blue

 Turn off i500

Press and hold the power button of i500 for 3 seconds

2.7.2 Placing on Desktop Cradle

2.7.3 Installation of Wall Mount Holder

3 Image Acquisition Software Overview

3.1 Introduction

The image acquisition software provides a user-friendly working interface to digitally record topographical characteristics of teeth and surrounding tissues using the i500 scanner.

3.2 Installation

3.2.1 System Requirement

	Laptop	Desktop
CPU	Above Intel Core i7-8750H	Above Intel Core i7-8700K
RAM	Above 16 GB	Above 16 GB
Graphic	Above Nvidia Geforce GTX 1060	Above Nvidia Geforce GTX 1060
OS	Window 10 64 bit	

Use PC and monitor certified IEC 60950, IEC 55032, IEC 55024

3.2.2 Installation Guide

① Run Medit_iScan_X.X.X.X.exe

② Select the setup language then click "Next"

③ Select the installation path

④ Read the "License Agreement" carefully before checking "I agree to the License ~" then click Install

3 Image Acquisition Software Overview

- ⑤ It may take up to several minutes to finish the recommended installation process. Please do not shut down the PC until installation is completed.
- ⑥ After installation is completed, we recommend restarting the PC to ensure optimal program operation.

- ⑦ If the scanner is connected, please disconnect the scanner from the PC by removing the USB cable.

4 Maintenance

4.1 Calibration

The calibration process is essential in producing precise 3D models. You should perform the calibration process periodically.

Calibration is required when :

- The quality of the 3D model is not reliable or accurate as compared to previous results.
 - Environmental conditions such as temperature has changed.
 - Calibration period has expired.
- You can set the calibration period as detailed in Menu > Setting > Calibration Period(Days)

The calibration panel is delicate component. Do not touch the panel directly. When the calibration process does not perform properly, check the panel. If the calibration panel is contaminated, please contact your service provider.

We recommend performing the calibration process periodically. You can set the calibration period via Menu > Setting > Calibration Period(Days). The default calibration period is set to 14 days.

How to calibrate the i500

- Turn on the i500 and connect the device to the software.
- Run Calibration Wizard from Menu > Setting > Calibration.
- Prepare the Calibration Tool and i500
- Turn the dial of the calibration tool to the position. **1**
- Put the Handpiece into the calibration tool.
Click “Next” to start calibration process.
- When the handpiece is mounted in the correct position , the system will automatically acquire the data at the position **1**.
- When data acquisition is completed at the position **1** , turn the dial to the next position.
- Repeat the steps for positions **2** ~ **8** and **LAST** position.
- When data acquisition is completed at the **LAST** position, the system will automatically calculate and show the calibration result.

4.2 Cleaning and sterilization procedure

4.2.1 Reusable tip

The i500 reusable tip is the part of the scanner that is inserted into the patient’s mouth during scanning. The tip can be reused for a limited number of times, but it must be cleaned and sterilized between patients to avoid cross-contamination.

- Clean the tip manually using soap water and a brush. Inspect the mirror of the tip after cleaning. If the mirror appears stained, smudged or has a milky haze, repeat the cleaning process using a soft dish brush and soap water. Dry the mirror carefully with a paper towel.
- Put the tip in a paper sterilization pouch.
- The pouch should be sealed. Make sure it is airtight by using either a self-adhesive pouch or a heat-sealed pouch.

4 Maintenance

- Sterilize the wrapped tip in an autoclave with the following conditions:
 - » At 121°C(249.8°F) for 30 minutes and 15 minutes dry

Use an autoclave program that dries the wrapped tip before opening the autoclave.

CAUTION

- The mirror in the tip is a delicate optical component. It should be handled with care as its clean and undamaged surface is critical for scan quality. Be careful not to or smudge it before scanning a patient.
- Never autoclave a tip that is not wrapped as this will leave stains on the mirror which cannot be removed. See the autoclave manual for more information.
- New tips must be cleaned and sterilized / autoclaved before using them for the first time.
- The scanner tip can be re-sterilized up to 20 times and must be disposed afterwards as described in the next section on Disposal.
- Medit will not responsible for any damage including distortion, blacken etc at other condition.

4.2.2 Mirror

Using unclean tip mirror for scanning may cause bad quality of final result or bad scanning experience. In this situation, you can clean the mirror following steps

- Remove the scanner tip
- Apply alcohol to a clean cloth or cotton-tipped swab and wipe the tip mirror until it is clean. Use alcohol that is free of impurities, as those may stain the mirror. Ethanol or propanol (ethyl-/propyl alcohol) can be used.
- Wipe the mirror dry with a dry, lint-free cloth to avoid leaving dust and fiber on the mirror.
- Make sure there are no dust or fibers on the mirror and repeat the above process until it disappears.

4.2.3 Handpiece

After treatment, clean and disinfect all other surfaces of the handpiece except scanner front (optical window) and end (air vent hole).

Cleaning and disinfecting must be done with the device turned off. Use the device after it is completely dry.

Recommended disinfection solution:

Denatured alcohol (aka. ethyl alcohol or ethanol) – typically 60-70% Alc/Vol.

Following is the general cleaning and disinfecting procedure.

- Turn off the device using power button.
- Unplug all cables from the power hub.
- Attach the handpiece cover to the front of scanner.
- Apply the disinfection solution to a soft lint-free non abrasive cloth.
- Wipe the scanner surface with a cloth dampened with disinfectant.
- Remove liquid with a clean, lint-free dry, non-abrasive cloth.

**CAUTION**

- Do not clean the handpiece with the device is turned on. The fluid may enter the scanner and cause malfunction.
- Use the device after it is completely dry.

4.2.4 Other components

- Apply the disinfection solution to a soft lint-free non abrasive cloth.
- Wipe the component surface with a cloth dampened with disinfectant.
- Remove liquid with a clean, lint-free dry, non-abrasive cloth.

4.3 Disposal**CAUTION**

- The scanner tip must be sterilized before disposal. Sterilize the tip as described above.
- Dispose the scanner tip as you would other clinical waste.
- Other components are designed to conform with following Directives.
 - » RoHS, Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment. (2011/65/EU)
 - » WEEE, Waste Electrical and Electronic Equipment Directive. (2012/19/EU)

4.4 Updates to Image Acquisition Software

The Image acquisition software automatically checks for updates when the software is running.

When a new version of the software is available, the system will automatically download the new version.

5 Safety Guide

Please adhere to all safety precautions as detailed in this user guide to prevent human injury and equipment damage. This document uses the words WARNING and CAUTION when highlighting precautionary messages.

The i500 system should only be operated by dental professionals and technicians who are trained to use the system. Using the i500 system for any purpose other than its intended usage as outlined in section "2.1 Intended Use" may result in injury or damage to equipment. Please handle the i500 system according to the guidelines in the safety guide.

WARNING

Please read the safety guide carefully. Failure to follow the safety precautions as highlighted in the guide may result in injury or damage to equipment.

5.1 System Basics

CAUTION

- The USB 3.0 cable connector to the Power Hub is the same as a regular USB cable connector. However, the device may not operate normally unless the regular 3.0 USB cable is used with the i500.
- The connector provided by the Power Hub is designed specifically for the i500 and should not be used with any other device.
- If the product has been stored in a cold environment, give it time to adjust to the temperature of the environment before use. If used immediately, condensation may occur which may damage the electronic parts inside the unit.
- The calibration process should be carried out if the temperature of the environment has changed.
- The calibration panel is delicate. Do not touch the panel directly. When the calibration process does not perform properly, check the panel. If the calibration panel is contaminated, please contact your service provider.
- Before using the system, check that there are no issues such as physical damage, loose parts, and wear. If there is any visible damage, do not use the product and contact the manufacturer or your local representative.
- Check the i500 body and its accessories for any sharp edges.
- Always keep an eye on the product and your patient while using the product to for abnormalities.
- If the equipment fails to operate normally, such as having issues with accuracy, stop using the product and contact the manufacturer or authorized resellers.

5.2 Safety and Warning

Carefully read and understand the guidelines, including all precautionary messages as prefaced by the words WARNING and CAUTION. To avoid bodily injury or equipment damage, make sure to adhere to the safety guidelines. All instructions and precautions as specified in the Safety Guide must be observed to ensure proper functionality of the system and personal safety.

5.2.1 Preventive Inspection Before Use of the System

CAUTION

Ensure that all components provided are free from physical damage. Safety cannot be guaranteed if there is any physical damage on the unit.

5.2.2 Modification of the System

CAUTION

- Modifications to the i500 system are prohibited by law, as they may compromise the safety of the user, patient or a third party.

5.2.3 Approved Software Only

CAUTION

- Install and use only approved programs to ensure proper functionality of the i500 system.

5.2.4 Proper Training

WARNING

- Before using your i500 system on your patients:
 - » You should have been trained on how to use the system, or you should have read and fully understood this user guide.
 - » You should be familiar with the safe use of the i500 system as detailed in this user guide.
 - » Before using or after changing any settings, the user should check that the live image is displayed properly in the camera preview window of the program.

5 Safety Guide

5.2.5 In Case of Equipment Failure

WARNING

- If your i500 system is not working properly, or if you suspect that there is a problem with the equipment:
 - » Remove the device from the patient and discontinue use immediately.
 - » Disconnect the device from the PC and check for errors.
 - » Contact the manufacturer or authorized resellers.
 - » The i500 system should not be dismantled arbitrarily.

5.3 Mechanical Hazards

WARNING

- If you drop the i500 tip on the floor, do not attempt to reuse it. Discard the tip immediately as there is a risk that the mirror attached to the tip may have been dislodged.
- If the i500 unit is dropped on the floor or if the unit is impacted, it must be calibrated before use. If the instrument is unable to connect to the software, consult the manufacturer or authorized resellers.
- When not in use, the i500 should be kept mounted on a desk stand or wall mount stand.
- Do not install the desk stand on an inclined surface.
- Carefully arrange all cables so that you or your patient do not trip or get caught in the cable. Any pulling tension on the cables may cause damage to the i500 body.

5.4 Explosion Hazards

WARNING

- The i500 system is not designed to be used near flammable liquids or gases, or in environments with high oxygen concentrations.
- There is a risk of explosion if you use the i500 system near flammable anesthetics.

5.5 Electrical Safety

WARNING

- The i500 system is a Class 1 device.

-
- To prevent electric shock, the i500 system must only be connected to a power source with a protective earth connection. If you are unable to insert the i500-supplied plug into the main outlet, contact a qualified electrician to replace the plug or outlet. Do not try to circumvent these safety guidelines.
 - The i500 system only uses RF energy internally. The amount of RF radiation is low and does not interfere with surrounding electromagnetic radiation.
 - There is a risk of electric shock if you attempt to access the inside of the i500 system. Only qualified service personnel should access the system.
 - Do not connect the i500 system to a regular power strip or extension cord as these connections are not as safe as grounded outlets. Failure to adhere to these safety guidelines may result in the following hazards:
 - » The total short circuit current of all connected equipment may exceed the limit specified in EN / IEC 60601-1.
 - » The impedance of the ground connection may exceed the limit specified in EN / IEC 60601-1.
 - Do not place liquids such as beverages near the i500 system and avoid spilling liquid on the system.
 - Condensation due to changes in temperature or humidity can cause moisture buildup inside the i500 unit, which may damage the system. Before connecting the i500 system to a power supply, be sure to keep the i500 at room temperature for at least 2 hours to prevent condensation. If condensation is visible on the product surface, the i500 should be left at room temperature for more than 8 hours.
 - You should only disconnect the i500 system from the power supply via its power cord.
 - The radiation characteristics of the i500 system makes it suitable for use in industry and hospitals. (CISPR 11 class A). If the i500 system is used in a residential environment (CISPR 11 class B), it may not provide adequate protection from radio frequency communications.
 - Before disconnecting the power cord, make sure to turn off the power on the device using the power switch on the main unit.
 - Only use the power adaptor supplied together with the i500. Using other power adaptors may result in damage to the system.

**CAUTION**

- Avoid pulling on the communication cables, power cables, etc. used in the i500 system.

5 Safety Guide

5.6 Eye Safety

 **WARNING**

- The i500 system projects a bright light from its tip during scanning.

The bright light projected from the tip of the i500 is not harmful to the eyes. However, you should not look directly at the bright light nor aim the light beam into the eyes of others. Generally, intense light sources can cause eyes to become brittle, and the likelihood of secondary exposure is high. As with other intense light sources exposure, you may experience temporary reduction in visual acuity, pain, discomfort, or visual impairment, which increases the risk of secondary accidents.

 **WARNING**

<Disclaimer for risks involving patients with epilepsy>

The Medit i500 should not be used on patients that have been diagnosed with epilepsy due the risk of seizures and injury. For the same reason, dental staff who has been diagnosed with epilepsy should not work with the Medit i500.

5.7 Pacemaker and ICD Interference Risk

 **WARNING**

- Do not use the i500 system on patients with pacemakers and ICD devices.
- It is forbidden to use the i500 system for patients using pacemakers due to the risk of interference.

5.8 Hygiene

 **WARNING**

- For clean working conditions and patient safety, ALWAYS wear clean surgical gloves when:
 - » Handling and replacing the tip.
 - » Using the i500 on patients.
 - » Touching the i500 system.
- The i500 main unit and its optical window should be kept clean at all times.
- Before using the i500 on a patient, be sure to:
 - » disinfect the i500 system
 - » use a sterilized tip

5.9 Safety when scanning with i500

CAUTION

- The i500 system is a high-precision optical measuring device. Any shock to the system caused by dropping or bumping the device may cause serious damage. Make sure to follow these guidelines:
 - » Always mount the i500 main unit on its desktop cradle or wall mount holder when not in use.
 - » Do not place any object on the i500 body.
 - » Do not pull or bend the cable connected to the i500.
 - » Do not spill any liquid on the i500 body.
 - » Do not place the i500 on any heated or wet surface.
 - » Always be sure to hold the i500 body securely when dismounting from its holder or cradle, as well as during scanning.
 - » Always place the power cord of the i500 system in an easily accessible location.
- Due to its fragile nature, i500 tips should be handled with care. To prevent damage to the tip and its internal mirror, be careful to avoid contact with a patient's teeth or restorations.
- Do not scan for more than 10 minutes at a time.

5.10 Protection from overheating

WARNING

- Do not block the air vents located at the rear of the i500 system. If the equipment overheats, the i500 system may malfunction or stop working.

5.11 Technical Maintenance

CAUTION

- Equipment maintenance should be carried out only by a MEDIT employee or a MEDIT-certified company or personnel.
- In general, the user does not need to perform maintenance work on the i500 system except for calibration, cleaning, disinfection and sterilization. Preventive inspections and regular maintenance are not required.

6 Electro – Magnetic Compatibility Information

6.1 Electro-Magnetic Emissions

This EUT is intended for use in the electromagnetic environment as specified below. The customer or the user of the EUT should ensure that it is used in such an environment.

RF Emissions CISPR 11 – Group 1

The EUT uses RF energy only for its internal functions.

Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

RF Emissions CISPR 11 – Class A

The EUT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Immunity Test	Compliance	Electromagnetic Environment - Guidance
Harmonic emissions IEC 61000-3-2	A	The EUT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ Flicker emissions	Complies	The EUT is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

6.2 Elector-Magnetic Immunity

This EUT is intended for use in the electromagnetic environment as specified below. The customer or the user of the i500 system should ensure that it is used in such an environment.

Electrostatic discharge (ESD) IEC 61000-4-2

Floors should be of wood, concrete or ceramic tiles. If floors are covered with a synthetic material, the relative humidity should be at least 30%.

IEC 60601-1-2 Test Level	Compliance Level
± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 15 kV air

Electrical fast transient/burst IEC 61000-4-4

Mains power quality should be that of a typical commercial or hospital environment.

IEC 60601-1-2 Test Level	Compliance Level
± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency

Surge Line-to-line IEC 61000-4-5

Mains power quality should be that of a typical commercial or hospital environment.

IEC 60601-1-2 Test Level	Compliance Level
± 0.5 kV, ± 1 kV	± 1 kV

Surge Line-to-ground IEC 61000-4-5

Mains power quality should be that of a typical commercial or hospital environment.

IEC 60601-1-2 Test Level	Compliance Level
± 0.5 kV, ± 1 kV, ± 2 kV	± 2 kV

Voltage dips IEC 61000-4-11

Mains power quality should be that of a typical commercial or hospital environment. If the user of the EUT image intensifier requires continued operation during power mains interruptions, it is recommended that the EUT image intensifier be powered from an uninterruptible power supply or a battery.

IEC 60601-1-2 Test Level	Compliance Level
0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°

6 Electro – Magnetic Compatibility Information

Voltage interruptions IEC 61000-4-11

Mains power quality should be that of a typical commercial or hospital environment. If the user of the EUT image intensifier requires continued operation during power mains interruptions, it is recommended that the EUT image intensifier be powered from an uninterruptible power supply or a battery.

IEC 60601-1-2 Test Level	Compliance Level
0 % UT; 250/300 cycle	0 % UT; 250/300 cycle

RATED power frequency magnetic fields (50/60Hz) IEC 61000-4-8

RATED power frequency magnetic fields (50/60Hz) IEC 61000-4-8 Power frequency magnetic fields should be at levels characteristic of a location in a typical commercial or hospital environment.

IEC 60601-1-2 Test Level	Compliance Level
30 A/m	30 A/m

Conducted RF IEC 61000-4-6

Portable and mobile RF communications equipment, including cables, should be used no closer to any part of the EUT than the recommended separation distance as calculated using the equation below, according to the frequency of the transmitter.

Recommended separation distance

$$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$$

$$d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$$

$$d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$$

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.

Interference may occur in the vicinity of equipment marked with the following symbol:

IEC 60601-1-2 Test Level	Compliance Level
3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

Radiated RF IEC 61000-4-3

Portable and mobile RF communications equipment, including cables, should be used no closer to any part of the EUT than the recommended separation distance as calculated using the equation below, according to the frequency of the transmitter.

Recommended separation distance

$$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$$

$$d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$$

$$d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$$

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.

Interference may occur in the vicinity of equipment marked with the following symbol:

IEC 60601-1-2 Test Level	Compliance Level
3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz

7 Specification

Model Name	i500
Rating	+9V = 4A
DC Adapter	
Model name	ATM036T-P090
Input voltage	Universal 100~240 Vac / 50~60 Hz input, without any slide switch Vacv/ 50~60 Hz input, without any slide switch
Output	+9V / 0~4A
Case dimension	100 x 50 x 33mm (W x L x H)
EMI	CE / FCC Class A, Conduction & Radiation met
Protection	OVP (Over Voltage Protection)
	SCP (Short Circuit Protection)
	OCP (Over Current Protection)
Protection against Electric shock	Class I
Mode of operation	Continuous
Handpiece	
Dimension	264 x 44 x 54.5mm (W x L x H)
Weight	280g
Applied part	Type BF
Power Hub	
Dimension	109.5 X 37 X 19.8 mm (W x L x H)
Weight	80g
Calibration Tool	
Dimension	165 x 55mm (H x Ø)
Weight	280 g

Operating & Storage conditions		
Operating condition	Temperature	18°C to 28°C
	Humidity	20 to 75% relative humidity (non-condensing)
	Air pressure	800 hPa to 1100 hPa
Storage condition	Temperature	-5°C to 45°C
	Humidity	20 to 80% relative humidity (non-condensing)
	Air pressure	800 hPa to 1100 hPa
Transport condition	Temperature	-5°C to 45°C
	Humidity	20 to 80% relative humidity (non-condensing)v20 to 80% relative humidity (non-condensing)
	Air pressure	620 hPa to 1200 hPa
Emission limits per environment		
Environment	Hospital environment	
Conducted and radiated RF EMISSIONS	CISPR 11	
Harmonic distortion	See IEC 61000-3-2	
Voltage Fluctuations and flicker	See IEC 61000-3-3	

EC REP EU representative

MERIDIUS MEDICAL EUROPE LIMITED
 Unit 3D, North Point House, North point Business Park, New Mallow Road
 CORK, T23AT2P, Ireland, +353 212066448

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