

USER GUIDE

Medit Link




Table of Contents

Medit Link

Introduction	6
Medit Link	6
Service Overview	6
Cloud Service	7
Workflow Management Platform with Ordering System	7

Medit Link > Installation

System Requirements	8
System Requirements for Windows	8
System Requirements for macOS	8
System Requirements for Web Service	8
Recommended for Browsers	8
Minimum for Mobile Devices	8
Installation on Windows	9
Installation on macOS	11

Medit Link > Medit Link Web

Overview	15
Dashboard	16
Dashboard Overview	16
Medit Insights	16
Subscribe Medit Insights	17
Unsubscribe Medit Insights	18
Case Box/Work Box	19
Case Filtering	19
Filter by status	19
Filter by duration	19
Filter by Search	20
Case Details	20
Order Box/In Box	22
Order Details	22
Tabs in Order Box/In Box	22
Case Talk	23
Creating a Case Talk	23
Setting the sharing permission type and sharing options	24
Case Talk details	26
Web Viewer	26
Annotation	27
Comments	28
Real-time screen sharing	29
Data Tree	29
Data Tree Controls	29
3D Model Controls	30
Partners	32
Search for Partners	32
My Partner List	32
Partner Invitation	33
Membership	34
Membership	34
My Plan	34
Cloud Storage	35

Subscribe	36
Unsubscribe	37
Billing Info	39
How to add a payment method	39
How to change a payment method	41
Redeem	42
Register Redeem Code	42
Activate Redeem Code	43
Payment	44
Payment Details (Clinic & Lab Accounts)	44
Payment Account Information (Lab Account Only)	44
My Info	45
My Account	45
My Organization (Admin Account Only)	45
My Scanner	46
My Price List (Lab Admin Account Only)	46
Other Menu	47
Notice & Update	47
Medit Information	47
Download Medit Link App	47

Medit Link > Medit Link App

Overview	49
Box Overview	49
Title Bar	51
My Information	51
Notifications	51
Case Synchronization Manager	52
Upload	53
Download	53
Network Status	54
Menu	55
Dashborad	55
Case Box/Work Box	55
Tabs in Case Box/Work Box	56
Functions in Case Box/Work Box	56
Order Box/In Box	57
Tabs in Order Box/In Box	58
Case Talk	58
App Box	58
Trash Box	59
Settings	59
Patient Management	60
Add New Patient	60
Delete Patient	60
Merge Patients	61
Notice & Update	62
Go to Medit Link Web	62
Medit Information	62
Settings	63
General	63
Storage Management	64
Set Period	64
Local Storage	64
Cloud Storage	65
Set Connected Program	65

Export Options	65
Case Converting Tool	67
How to Use Case Converting Tool	67
Proxy Options	69
About	69
User Guide	70
File Viewer	71
Data Tree	71
Model Control Tools	71
3D Model Control	71
2D Image Control	72
Function Tools	72
Image Capture and Memo	73
Ortho Capture	73
Fill Holes	74
Trimming	74
Optimize Occlusion Alignment	77
Mandibular Movement	78
Manage Multi Occlusion Data in File Viewer	79
How to manage multiple results	79
View Options	81
Status Filters	81
Table Items	82
View Type	83
Gallery View	84
List View	85
Grouping	85
Sorting	86
Filtering and Search	88
Filter by Date	88
Advanced Search Capabilities	88

Medit Link > Workflow

Create/Edit/Clone/Delete Case	90
Creating Cases	90
Create Case in Case Box (Clinic)	90
Create Case in Work Box (Lab)	91
Editing Cases	91
Edit Case in Case Box (Clinic)	91
Edit Case in Work Box (Lab)	92
Cloning Cases	92
Clone Case in Case Box (Clinic)	92
Clone Case in Work Box (Lab)	93
Deleting Case	93
Delete Case (Clinic and Lab)	93
Form Information	94
Registering Form Information in Case Box	94
Enter on the Form Tab	94
Enter on the CAD tab	94
Registering Form Information in Work Box	94
Enter on the Form Tab	94
Enter on the CAD tab	95
Connecting Teeth into Bridge	95
Setting Scanning Options	95
Medit Scan for Labs	95
General options	95

Options by teeth	95
Adding to Favorite	96
Add to Favorites	96
Manage Favorite Templates	96
Work on Scanning/CAD	97
Scan Process	97
CAD Process	97
Get Result Files	98
Export Scan Result Files	98
Export Selected Files	98
Export Margin Line	99
Export Margin Line as a Color	100
Export Margin Line as a File	100
Export Case List as CSV	101
How to export the case list	101
Download Files from Medit Link Web	102
Case Box/Work Box	102
Web Viewer	102
How to download	103
Order/Mange Orders	105
Ordering Cases	105
Confirm Order (In Box)	106
Cancel Order (Order Box)	106
Reject Order (In Box)	106
Checking Ordered Cases	106
Entering Shipment Information (In Box)	107
Managing Delivery Date	107
Setting Requested Delivery Date (Clinic account)	107
Completing Case	109
Complete (Work Box/In Box)	109
Complete a case from the case detail	109
Complete the multiple cases from the case list	110
Complete (Case Box/Order Box)	110
Complete a non-ordered case in the Case Box	110
Communicate Using Tags	112
Add Tags When Creating a New Case	112
Add Tags When Ordering	113
Check Tags on Case Details	113

Introduction

Medit Link

Medit Link is a collaborative tool designed to facilitate communication between clinics and labs to provide an optimized workflow for each party.

It offers differentiated features based on the type and permissions of each account. The program provides a user-friendly user interface for managing the workflow, from scanning with oral and model scanners, ordering cases, and designing and producing prostheses.

Medit Link consists of a web service and a computer application.

Service Overview

		Web	App
		• Register and manage the account • Search • Search and manage partners • Share cases via Web Viewer with a sharable link • Access cases made from the App • Manage orders	• Enter form information • Scan • Place/receive orders and make payments • Link with CAD • Manage patients
Clinic	Admin	Search and manage partners and your organization, manage payment information, and share cases. - Dashboard - Case Box - Order Box - Case Talk - Partners - Membership - Payment - My Info	Enter form information, scan, link with CAD, place orders, make payments, and manage patients. - Dashboard - Case Box - Order Box - Case Talk - App Box - Trash Box - Patient Management
	Member	- Case Box - Case Talk - My Info	- Case Box - Case Talk - App Box - Trash Box
Lab	Admin	Search and manage partners and your organization, manage payment information, and share cases. - Dashboard - In Box - Work Box - Case Talk - Partners - Membership - Payment - My Info	Enter form information, scan, link with CAD, and receive/manage orders. - Dashboard - In Box - Work Box - Case Talk - App Box - Trash Box
	Member	- Work Box - Case Talk - My Info	- Work Box - Case Talk - App Box - Trash Box

iNote

The following service is not available on Medit Link China.

- Payment

Cloud Service

Medit Link supports the automatic upload/download of databases and files through AWS (Amazon Web Services).

The data stored in the cloud is automatically synchronized to the computers of multiple users in the organization so that they can access the same data from any computer at any time. Cloud-based backup and recovery keep all your data safe. It keeps your data safe even in case of data loss due to sudden computer failure.

The cloud service allows partnered dentists and labs to quickly share data without the hassle of emailing or manually transferring the scans. All data including patient information is securely encrypted, stored, and transmitted. Medit complies with medical laws and regulations like HIPAA and GDPR.

Workflow Management Platform with Ordering System

Medit Link offers its own ordering and payment system.

- Medit Link allows you to create contracts between dentists and dental laboratories and share data quickly and easily.
- You can check and manage the ordering status.

Credit card-backed system works with payment modules of online payment platforms, such as Stripe (www.stripe.com) and I'mPort (<http://www.iamport.kr>), making payment services easy and convenient.

iNote

The following service is not available on Medit Link China.

- Payment

System Requirements

iNote

Please refer to the system requirements for Medit Scan for Clinics and Medit Scan for Labs if you want to use the scanning software with your scanner.

System Requirements for Windows

	Laptop & Desktop
CPU	Intel Core i5 2.6 GHz or higher
RAM	16 GB or higher
Graphics	NVIDIA GeForce GTX 760 (2GB) or higher / or equivalent AMD video card
OS	Windows 8 64 Bit (unavailable in 32 Bit) or higher

System Requirements for macOS

	Mac
Processor	M1 with 8-core, 7-core GPU, 16-core Neural Engine
RAM	16 GB or higher
OS	Monterey 12

System Requirements for Web Service

Recommended for Browsers

	Browsers
PC	Chrome, Firefox, Safari, Edge (Chromium)
Android	Chrome, Firefox
iOS	Safari, Chrome, Firefox


Minimum for Mobile Devices

	Android	iOS
Device	Galaxy S7 equivalent or higher	iPhone 7 equivalent or higher
OS	Android 7.1 Nougat	iOS 11.0
RAM	3 GB or higher	2 GB or higher
AP	Above Snapdragon 710	A10 Fusion


Installation on Windows

Medit Link is installed as a package with scan acquisition programs - Medit Scan for Clinics and Medit Scan for Labs.


1. Run Medit_LinkSetup.exe.
2. Select a language for installation.
3. Select the folder to install Medit Link. Please read and agree to the license terms and conditions and click "INSTALL."



4. Make sure that the scanner cable is disconnected from the computer before proceeding to the next step. The installation may not proceed if the scanner is connected to your PC.
5. The program will be installed automatically and may take a few minutes. Do not turn off or unplug your PC until the installation is complete.



6. Please restart your PC after the installation is complete. You can proceed without restarting the PC if you do not have the scanner.




iNote

If you do not restart your PC after the installation is complete, the scanner may not function properly.


Installation on macOS

Medit Link is installed as a package with Medit Scan for Clinics.


1. Run Medit_Link_x.x.x.pkg and click "Continue."




2. Read the End User License Agreement and click "Continue."



3. Select the disk to install and click "Continue."




- Click the "Change Install Location..." to change the installation destination.




- Click "Install" and type your Mac password to allow installation if requested.




6. Wait until the installation is completed.



7. Click "Close" to exit the installation window.



Overview




The screenshot displays the 'Case Box' interface of the Medit Link application. The interface includes a sidebar with navigation options like Dashboard, Case Box, Order Box, Case Talk, Partners, Membership, Payment, and My Info. The main content area shows a table of cases with the following columns: Status, Case Name, Patient Name, Form Information, Last Modified Date, Scanning Date, Order Date, and Lab Name. The table contains 10 rows of data, each with a 'Share' button. The status of each case is indicated by a colored pill: Scan Completed (blue), Form Only (orange), Need Processing (green), Scan Completed (blue), CAD (blue), Ordered (yellow), and Case Completed (green).

Status	Case Name	Patient Name	Form Information	Last Modified Date	Scanning Date	Order Date	Lab Name
Scan Completed	mandibular movement - Clone -...	Test Occlusion	13, 14, 15, 16	2023-04-06 16:48	2023-04-03 10:11	-	-
Form Only	Occlusion Test's Case	Occlusion Test	14, 15, 16, 17	2023-04-06 14:40	-	-	-
Ordered	Occlusion Test's Case	Occlusion Test	21, 22, 23, 24, 25, 26	2023-04-06 14:13	2023-04-04 14:17	2023-04-06 14:13	WonhoonChoi_Lab...
Case Completed	mandibular movement - Clone	Test Occlusion	13, 14, 15, 16	2023-04-06 11:55	2023-04-03 10:11	-	-
Case Completed	test_initial's Case	test_initial		2023-04-06 11:55	-	-	-
Case Completed	mandibular movement	Test Occlusion	13, 14, 15, 16	2023-04-06 11:55	2023-04-03 10:11	-	-
Case Completed	COL-20532 test	test	12, 13, 14, 15, 16	2023-04-06 11:53	-	-	-
Case Completed	Order-test	test	23, 24, 25, 26, 27	2023-04-05 18:36	-	2023-03-28 16:46	WonhoonChoi_Lab...
Form Only	test_initial's Case	test_initial		2023-04-04 13:38	-	-	-
Form Only	test_initial's Case	test_initial		2023-03-31 08:18	-	-	-

Dashboard

Dashboard Overview

You can view the statistics for the entire work done on Medit Link in visualized graphs and charts – both in the Medit Link App and Web.



Graphs and charts in the dashboard show statistics for the entire work done in Medit Link.

Scanner	Provide information about the active scanners, including warranty, scanner usage time, calibration interval, and calibration elapsed time.
Cloud Storage	Display the total capacity of the cloud and the amount in use.
Statistics	Provide statistics on cases based on job status.
Order Status	Provide statistics on cases based on order status.
Graph	Graph job information by period.
Tooth Type	Provide statistics on the number of scans and orders by treatment information.


Medit Insights

Medit Insights is a reporting service that helps you analyze your usage statistics.

iNote

The following service is not available on Medit Link China.

- Medit Insights



Subscribe Medit Insights

You can receive a report email on the second day of every month. No additional sign-up process is required. Just enable "Subscribe Medit Insights."

If you don't use the email address you used for the Medit Link account, you can set another email address to receive Medit Insights.

Settings

Receive the report via


☒ Medit Link account
 ☐ Another email address

Enter another email address to receive the report email.

The report will be sent on the 2nd day of every month.

Cancel
OK

On the second of every month, you will receive an email with statistics as shown below.




Banner	Medit provides useful tips and information here. * If you don't want to see the banner, uncheck the Marketing and Promotion Agreement check box in My Info > My Account.
Statistical Information	The statistics information is displayed, such as scanning time, created cases, most ordered lab, etc.
Medit Link Download	The button provides the download link of the Medit Link.
Dashboard	All statistical information is based on your Dashboard information. You can see more details in Dashboard.

Unsubscribe Medit Insights

Disable "Subscribe Medit Insights" in Dashboard or click "Unsubscribe" in the received email.

Case Box/Work Box


The Case Box (Clinic) or Work Box (Lab) shows all cases the user has created on the Medit Link application. You can see files via a web browser without installing a separate program.



Case Filtering

Filter by status

You can filter cases based on their status.



- The user can click the status they want to view. The browser will display only cases with the chosen status.
- The user can select multiple statuses at once.

The following status type is provided based on account type:

Case Box (Clinic)	All / Form Only / Need Processing / Scan Completed / CAD / Ordered / Case Completed
Work Box (Lab)	All / Form Only / Need Processing / Scan Completed / CAD / Done


Filter by duration

You can filter and display cases that were last modified within your set period.

1. Click the duration filter.



2. Select the range.




- All: Shows all cases.
- Custom Range: Specify start and end dates to display cases modified within that period.
- Today: Show cases last modified today.
- Yesterday: Show cases last modified with yesterday's date.
- Last 7 days: Show cases last modified during the last 7 days, excluding today.
- Last 30 days: Show cases last modified during the last 30 days, excluding today.
- Last 90 days: Show cases last modified during the last 90 days, excluding today.
- This Month: Show this month's cases.
- Previous Month: Show the previous month's cases.

The case list will be updated in real time with the cases created or last modified within the selected duration.

Filter by Search


You can leverage search capabilities to search for the required information easily and quickly.



Search options based on clinic/lab account:

Case Box (Clinic)	Case Name, Patient Name, Patient ID, Lab Name
Work Box (Lab)	Case Name, Patient Name, Clinic Name, Order ID

Case Details



The following case management tools are provided in Case Box/Work Box.

	Case Status	Show the status of the case.
	Go to Detail	Provide a shortcut for the same case in another box (Case Box ↔ Order Box; In Box ↔ Work Box).
	Download	Download files in the case to the local PC.
	Case Information	Review case details and save them as a PDF or print them out.
	Memo & Tags	Add a memo or tags when registering or ordering a case. Use them to communicate any additional or supplementary information about the case to your partners.
	Share	Create a link for sharing. Users with the shared link can access the 3D data in a browser such as Chrome, Firefox, or Safari without installing additional programs.

Order Box/In Box

The Order Box (Clinic) or In Box (Lab) shows ordered cases on the Medit Link application.

Order Details

The screenshot shows the Medit Link application interface. The top bar includes a search bar with 'test (test)' and a user profile 'wonhoon choi'. The sidebar on the left contains navigation icons for Dashboard, Case Box, Order Box, Case Talk, Partners, Membership, Payment, and My Info. The main content area is divided into two sections: 'Order' and 'File Viewer'. The 'Order' section is active and shows detailed information for an order with ID 5995726. The 'Order List' section displays a table of ordered products.

No.	Type	Method	Material	Shade	Price
14	Inlay	-	Multi-Layered Zl...	None	0.00
15	Inlay	-	Multi-Layered Zl...	None	0.00
16	Inlay	-	Multi-Layered Zl...	None	0.00
23	Inlay	-	Multi-Layered Zl...	None	0.00
24	Inlay	-	Multi-Layered Zl...	None	0.00
25	Inlay	-	Multi-Layered Zl...	None	0.00
26	Inlay	-	Multi-Layered Zl...	None	0.00
27	Inlay	-	Multi-Layered Zl...	None	0.00
28	Onlay	-	Glass Ceramic	A3.5	0.00
Total					USD 0.00

1	Show case name and patient name.
2	Provide functions for order management (accepting/canceling; delivering/receiving).
3	Show available tabs: Order, File Viewer
4	Show functions for managing cases and files: case status, shortcuts for detailed information, case information, memo & tags, etc.
5	Show the detailed information for the order including the order number, patient name, partner name, scanning date, order date, delivery date, and memo & tags.
6	Show the ordered product information and delivery status.


Tabs in Order Box/In Box

Order	File Viewer
Check the detailed information for the order, such as order ID, patient name, account name, scanning date, order date, delivery date, memo, and tags.	View files generated by scan/CAD operations. The file viewer displays and manages the acquired or attached 3D models and 2D images.

Case Talk

Cases created in the Medit Link can be shared with specific people you specify, anyone with the shared link or registered members of Medit Link.

The user can see a list of all shared cases through "Case Talk" in Medit Link App and Medit Link Web.




Shared cases can be viewed via a web browser without installing a separate program, enabling easy communication over the 3D data on various devices.


Creating a Case Talk

You can create a shared link to a case that was created in Case Box/Work Box.

1. Find a case you want to share and click the "Share" button on the right.



2. Select the files you want to share.



3. Set the sharing permission type and sharing options.
4. Click "Get Link" to create a shared link.

Setting the sharing permission type and sharing options

Share through Case Talk



Share With

Share with specific people

Create a shared link and share it with specific people.

Share via

Set the expiration date for the link. After the link expires, the people you shared it with will no longer have access to files. You can extend the expiration period up to 90 days from today.

Expiration Date

2022-06-28

You can extend the period up to 90 days from today.

PIN Code

☐ PIN Code

Set the PIN code using six digits for privacy protection.

☒ Share the patient's name

☒ Grant Download Permission Premium Feature

If you check this option, all users who have access to the shared link can download the files.

The premium feature is provided to all users for a limited time. Feature availability is subject to change depending on the membership plan.

☒ Allow comments

If you check this option, you will be able to communicate through the shared link.

Comment

Please enter your comment.

0/1024

By pressing the "Get Link" button, you agree that you understand personal information can be collected and provided to a third party (or shared in any public domain) ONLY with the consent of the person whose images on which you have worked. You also agree the consent of that person was truthfully obtained, and Medit is not responsible for any problems which may arise there after.


Close

Get Link

Share With	- Share with anyone Creates a shared link. You can share the files with anyone. - Share with specific people Creates and sends a shared link to specific people.
Share via	If you select "Share with specific people," you can share the link with specific people in various ways. - Email Address - SMS - Partners
Expiration Date	Set the expiration date for the link. After the link expires, the people you shared it with will no longer have access to files. You can extend the expiration period up to 90 days from today.
PIN Code	Set a 6-digit access code to ensure data security.
Share the patient's name	If you check this option, all people with the link will see the patient's name.
Grant Download Permission	If you check this option, all users with the link can download the files.
Allow comments	- If you check this option, you will be able to communicate through the shared link. - If you uncheck the option, the recipient can only view your comments but will not be able to add their own.
Comment	You can leave a message in the comment section.

Case Talk details

You can see files in the shared case. The following functions are available on the Case Talk details page.



	Download	Download all the files in this case to the local PC.
	Case Information	Show detailed information about the case.
	Share	Modify the share setting.
	Copy Link	Copy the shared link to the clipboard. The copied link shows the shared files on your web browser.
	Live	Share your screen with the connected users. The screen sharing is synchronized with the connected users.
	Attendees	Show all users who are accessing the shared link.
	Delete	Delete the shared link. The deleted shared link is no longer accessible.

iNote


The following services are not available on Medit Link China.

- Live
- Attendees

Web Viewer

You can view shared data through Case Talk in a web browser.

It is possible to communicate with the people who have access to the link using annotations and comments.



	Download	Download all the files in this case to the PC.
	Case Information	Show detailed information about the case.
	Share	Modify the share settings.
	Live	Share your screen with the connected users. The screen sharing is synchronized with the connected users.
	Attendees	Show all users who are connected to the shared link.
	Tooth Labeling	Change the teeth labeling system.


iNote

The following services are not available on Medit Link China.

- Live
- Attendees

Annotation

You can add annotations to a specified location on the 3D model. Anyone with the link can view the annotations you add.



iNote

The following service is not available on Medit Link China.

- Annotation

1. Find the menu icon in the top left corner and select the "Annotation" tab. Use the "New Annotation" button to add an annotation anywhere on the data with a click.
2. Type the message, then click "OK" to save it.

Create Annotation

Please enter your annotation.

Please enter your message

Cancel
OK

3. You can use annotations as a simple communication channel. Click the created annotation to enter a reply and hit "Send."

wonhoonchoi Beta wonhoonchoi
🗑️ ✕

Annotation 1

1 09-09 20:31

Send

Comments

You can communicate with users who have access to the link through the "Comment" tab in the Web Viewer menu.

iNote

The following service is not available on Medit Link China.

- [Comment](#)

Real-time screen sharing


You can share your screen in real-time.

iNote

The following service is not available on Medit Link China.

- [Live](#)

1. Click "Live" to share your Web Viewer screen with the connected users.





2. Click "OK" to start sharing your screen.
3. Other people who have accessed the link must accept the screen sharing.
4. The screen of the user who initiated screen sharing is displayed in sync with the screen of the person who received the request.

Data Tree





All data is grouped and displayed according to attributes in the data tree. You can show or hide data, adjust opacity, and so on.

Data Tree Controls

Show/Hide	Click the 'eye' icon to show or hide data from the view. Control data visibility as a group or individual file.
Opacity	Hover the mouse over a visible element to see the opacity slide bar. Adjust the slider to change the transparency of your data. 

Extended Menu	The extended menu is provided on the right side of each Data Tree element. It includes the following options:  • Show This Only: Hide all other data except for the current one. • Hide: Hide the current data. • Duplicate: Clone data within the current case. • Export: Save data to a local PC. • Delete: Delete the current data or data group.
---------------	--

3D Model Controls

	Pan	Move the model.
	Rotate	Rotate the model.
	Zoom In/Out	Zoom in and out on the model.
	Zoom Fit	Position the model in the center of the screen.
	Model Display Mode	Show data in real texture colors. 
		Show data in a single color. 
	Grid Settings (mm)	Show the grid in the background. 
		Overlay the grid over the model. 



Grid Off

Hide the grid.

Partners

You can search for potential partners to request a partnership and view the list of current partners and partners with contracts.


iNote

Services using Google Maps are not available on Medit Link China.

Search for Partners

You can search for Medit Link users based on their business name and/or address.

In the search results, you can send a partnership request to an organization you want to work with.




My Partner List

You can review your current partnerships and contracts with all your existing partner businesses in My Partner List.

How is a contract different from a partnership?

Once you partner with a business, you can request to make a contract with that organization.

- If your partner accepts your contract request, you can proceed with payment for the placed orders.



iNote

The following service is not available on Medit Link China.

- **Contract**

Partner Invitation

If you have an organization you want to work with, you can invite them to join Medit Link. The partnership will be automatically formed once they accept the invitation.

1. Click "Invite Partner" in the top left of your screen.
2. Enter the email address of the vendor you want to work with and send them an invitation with a partnership request.



Invite Partner

Invite your desired partner to join Medit Link. When the lab accepts your invitation, a partnership is automatically formed.

Email Address

Cancel
Send Invitation

3. Once the invited organization accepts the invitation and signs up for Medit Link, the partnership will be automatically established.

Membership

Medit Link's membership policy aims to establish a more extensive digital environment, gradually expanding to provide various services based on the cloud platform.

We aim to provide our customers with a more consistent and personalized experience with our new membership policy.

There are two membership plans available:

Standard

- Free standard plan for all signed-in users
- 1 GB cloud storage
- For new users

Premium

- Paid plan
- 10 TB cloud storage
- For users with many cases and various data
- One-month free trial

	Standard	Premium
Cloud Storage	1 GB	10 TB
Monthly Fee	Free of charge	\$0.99/month
* Cloud Storage does not calculate the space taken up by Raw Data.		

iNote

The membership plans will be available on Medit Link China. You can get the Premium plan for free until the official release of the feature.


Membership

The screenshot displays the Medit Link web interface. On the left is a sidebar with navigation icons. The main content area is titled 'Membership' and includes tabs for 'Membership', 'Billing Info', and 'Redeem'. The 'My Plan' section shows the 'Standard Free' plan with a list of benefits: 'Free plan for every member', '1 GB of cloud storage', and 'Efficient for new users'. A 'Subscribe' button is at the bottom. The 'Cloud Storage' section shows a donut chart indicating 353.21 MB / 1 GB usage. Below the chart is a table of storage usage:

Category	Usage
Raw Data	11.97 GB
Result Files & Attachments	353.21 MB


My Plan

In the left-side section, you can see what membership plan you currently have. Here you can also change your membership subscription status at any time.



Cloud Storage

In the right-side section, you can manage your cloud storage.



Set Period	Set the specific period to delete data created within that time frame. • All: See all files. • 1 Week Ago: See the files that have passed a week after creating. • 1 Month Ago: See the files that have passed a month after creating. • 6 Months Ago: See the files that have passed six months after creating. • 1 Year Ago: See the files that have passed a year after creating.
------------	---


Raw Data	The raw data does not take up space in the cloud storage, so it will not be in the pi-chart. But you can manage the raw data from here if needed, e.g., delete it to protect the patient's private information, etc.
Result & Attachment	The result and attachment files take up space in the cloud. If the cloud is full, you can secure more storage space by removing unnecessary files.
Delete	Remove selected files from the cloud storage to free up space.

Subscribe


If you are currently on the Standard plan or the Premium plan with a redeem code, you can subscribe to the Premium membership.

The payment is completed through Stripe service (www.stripe.com).

1. Click the "Subscribe" button.



2. Check the Premium plan description and click "Subscribe" again.



3. Check the price and conditions and input information necessary for the payment. Then, click "Subscribe."

← MEDIT Link TEST MODE

Subscribe to Premium membership for Clinic

\$0.00

Then \$0.99 per day

Medit Link Membership

Premium membership for Clinic	\$0.99
Billed daily	
Subtotal	\$0.99
1 MONTH FREE – for first subscription	
100% off for a day	-\$0.99
Total due today	\$0.00

Powered by stripe | Terms Privacy

Pay with card

Email

Card information

1234 1234 1234 1234

VISA

MM / YYCVC

Name on card

Country or region


South Korea

☐ Save information to pay faster next time

Subscribe

By confirming your subscription, you allow MEDIT SG PTE. LTD. to charge your card for this payment and future payments in accordance with their terms.


4. The payment may take several seconds. When done, you will see the confirmation message.




Unsubscribe

You can unsubscribe from the membership at any moment. After unsubscribing, you can still use the Premium membership for the remaining days of that paid month. Your plan will be changed to Standard from the next payment day, and the cloud storage will be limited to 1 GB. If the cloud storage capacity is exceeded, the files will not be backed up, and the order and share functions will not be available.

1. Click "Unsubscribe."



2. Click "Cancel Subscription."



3. To confirm that you want to unsubscribe, click "Cancel Subscription" again in the new pop-up window.

Confirm Unsubscribe

Your subscription will be canceled at the end of your billing period on 2021-09-10

Cancel Subscription Stay Subscribed

4. Your subscription has been canceled.
5. Your plan will be changed to Standard from the next payment day, and the cloud storage will be limited to 1 GB. If the cloud storage capacity is exceeded, the files will not be backed up, and the order and share functions will not be available.

Confirm Unsubscribe

Your subscription has been canceled.
 Your current plan can use up to 2021-09-10. After that, the cloud capacity is limited to 1GB.
 If the cloud capacity is exceeded, the files will not be backed up, and the order and share functions will not be available.

Close

Billing Info

You can check the payment history for membership subscriptions, including an invoice for each payment. You can also add or change the payment method, which will be applied starting from the next scheduled payment.

MEDIT Link

Dashboard

Case Box

Case Task

Partners

Membership

Payment

My Info

Home

Medit Notice

Terms & Conditions

Privacy Policy

Copyright © Medit Corp. All rights reserved.

Membership

Billing Info

Redeem

My Plan and Billing Information

Plan

Premium

Amount

0.99\$ USD / Monthly

Next Payment Date

2021-09-09

Payment Method

Visa *****4242 (11/21)

Stripe

Edit Payment Method

If you haven't subscribed yet?

Experience a differentiated premium.

Billing History

Payment Date	Membership	From	To	Payment Amount	Status
2021-09-09	Premium	2021-09-09	2021-09-10	0.99 \$ (USD)	Payment Completed
2021-09-09	Premium	2021-09-09	2021-09-09	0 \$ (USD)	Payment Completed

How to add a payment method

1. Click "Edit Payment Method."

MEDIT Link

Dashboard

Case Box

Case Task

Partners

Membership

Payment

My Info

Home

Medit Notice

Terms & Conditions

Privacy Policy

Copyright © Medit Corp. All rights reserved.

Membership

Billing Info

Redeem

My Plan and Billing Information

Plan

Premium

Amount

0.99\$ USD / Monthly

Next Payment Date

2021-09-09

Payment Method

Visa *****4242 (11/21)

Stripe

Edit Payment Method

If you haven't subscribed yet?

Experience a differentiated premium.

Billing History

Payment Date	Membership	From	To	Payment Amount	Status
2021-09-09	Premium	2021-09-09	2021-09-10	0.99 \$ (USD)	Payment Completed
2021-09-09	Premium	2021-09-09	2021-09-09	0 \$ (USD)	Payment Completed

2. Click "Add payment method."

39

Billing

CURRENT PLANS

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

PAYMENT METHOD

 **** 4242 Default Expires 11/2021

[+ Add payment method](#)

3. Input the required information and click "Add."

Add payment method

Card information

	Card number	MM / YY
---	-------------	---------

☒ Use as default payment method

Add

Go back

By adding a payment method, you agree to MEDIT SG PTE. LTD.'s
Terms of Service and Privacy Policy.

4. The newly added payment method will be set as default and used for the next payment.

Billing

CURRENT PLANS

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

PAYMENT METHODS

 **** 5556 Default Expires 11/2021 ×

 **** 4242 Expires 11/2021 ...

+ Add payment method

How to change a payment method

1. Click the expanded menu icon next to the payment method you want to change in the billing dialog window.

Billing

CURRENT PLANS

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

Premium membership for Clinic

\$0.99 per day

Your plan renews on September 10, 2021.

PAYMENT METHODS

 **** 5556 Default Expires 11/2021 ×

 **** 4242 Expires 11/2021 ...

+ Add payment method

2. Click "Make default" to set the selected card as the default payment method.

The screenshot shows a 'Billing' page with two main sections: 'CURRENT PLANS' and 'PAYMENT METHODS'. Under 'CURRENT PLANS', there are two identical entries for 'Premium membership for Clinic' at '\$0.99 per day', both renewing on September 10, 2021. The 'PAYMENT METHODS' section lists two Visa cards. The first card, ending in 5556, is marked as 'Default'. The second card, ending in 4242, has a red box highlighting the 'Make default' and 'Delete' buttons. Below the cards is a '+ Add payment method' link. The 'BILLING INFORMATION' section is partially visible at the bottom.

Redeem

Register the redeem code provided as a part of the promotion. You can use the Premium membership together with redeem code.

Register Redeem Code

1. Input redeem code and click "Register."

The screenshot shows the 'MEDIT Link' interface with the 'Redeem' tab selected. The 'Register Redeem Code' section features a text input field for the 'Redeem Code' and a blue 'Register' button. Below this is a table titled 'My Redeem Code' with columns: Promotion, Code Number, Dealer Name, Registration Date, From, To, and Status. The table currently displays 'No search results were found.' The left sidebar contains navigation icons for Dashboard, Case Box, Case Talk, Partners, Redeem, Payment, and My Info. The footer includes links for Home, Medit Notice, Terms & Conditions, Privacy Policy, and a copyright notice for Medit Corp.

2. Read the notice and choose "Register" or "Apply" in the pop-up dialog.
 - **Register:** Register the redeem code. You can activate and use the redemption code later after registration.
 - **Apply:** Register and activate the redemption code. The promotion discount applies immediately.

Register Redeem Code

You can register for the promotion or apply immediately.
 Register allows you to apply and activate when you need it, If applied immediately, it will start after the expiration of the active promotion.
 Periodic payments are stopped automatically.

Close
Register
Apply

Activate Redeem Code

You can register multiple redeem codes, but you need to activate the one you want to use.

1. Click the "Activate" button next to the code you want to use.
 If there is no code with the "Active" status, the one you choose will be applied immediately. If any code is already activated, its status will change to "Waiting," and the one you chose will be used instead. When the current redeem code expires, the waiting one will be automatically re-activated.

2. **Note**
 Only one code can be in "Waiting" status.

The screenshot shows the 'Medit Link' web application interface. The top navigation bar includes 'Membership', 'Billing Info', and 'Redeem'. The left sidebar contains various icons for navigation. The main content area is titled 'Redeem' and contains two sections: 'Register Redeem Code' and 'My Redeem Code'.

The 'Register Redeem Code' section has a text input field for the 'Redeem Code' and a 'Register' button.

The 'My Redeem Code' section displays a table with the following columns: Promotion, Code Number, Dealer Name, Registration Date, From, To, and Status. The table contains four rows of test data.

Promotion	Code Number	Dealer Name	Registration Date	From	To	Status
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Yoshida	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate
Test Promotion - 1 month	[Redacted]	Medit Link	2021-09-09	-	-	Registered Activate

Payment

iNote


The following service is not available on Medit Link China.

- Payment

Payment Details (Clinic & Lab Accounts)

You can check the payment and transaction history between you and your partners through the Payment Detail page.

Payment Account Information (Lab Account Only)

The screenshot shows the 'Payment Detail' page header. On the left, it says 'Payment Detail' and 'Check the details of payment information here.' Below that, it says 'Total Number of Transactions : 0 Case'. On the right, there is a date range selector showing '2020-08-12 ~ 2020-09-10', a dropdown menu for 'All Partners', and a print icon. There are also settings and help icons in the top right corner.

This feature is only available to the lab administrator account to allow them to enter information to get paid by your partners.

- You must enter your payment account information for your partners to place orders through the online payment system.

iNote

If you select a country where Stripe is supported, you will be directed to Stripe by clicking a button to sign up for Stripe.

You can select South Korea, then you can enter a local bank account number to receive payment, and KCP will process your payment.

My Info

My Account

Here you can review and edit the information you provided during registration and check your membership status.

[My Account](#) [My Organization](#) [My Scanner](#) [My Price List](#)

Account Information



1. JPG, GIF, or PNG only. Max size of 10 MB.

Email Address
coldgreentea@naver.com [Change Email Address](#)

Password
[Change Password](#)

Name
wonhoonchoi_Lab3

Language
한국어 (대한민국) ▼

Country
Korea, Republic of ▼

2. Your chosen language is applied to the email sent by Medit Link and the notification.

Subscription Information

[Subscribe](#)

Standard

Organization Information

[Edit Information](#) >

Organization Name
wonhoonchoi_Lab3

My Organization (Admin Account Only)

Here you can review and edit the information about your organization. This section is only available to administrator-type accounts.

[My Account](#) [My Organization](#) [My Scanner](#) [My Price List](#)

[Address Information](#) [Organization Information](#) [Additional Information](#)

Address Information (Mandatory) *

Address 1 *
23 Incheon-ro 22-gil, Seongbuk-gu, Seoul, Korea [View Location on Map](#)

1. Street Address, P.O. Box, Company Name, C/O
2. Please enter the correct address so labs can check your company's location on the map.

Address 2
Apartment, Suite, Unit, Building, Floor, etc.

1. Apartment, Suite, Unit, Building, Floor, etc.

Zip/Postal Code *
02835

City *
Seoul

Country *
Korea, Republic of ▼

State/Province *
Seoul

1. If you have a contract history, you cannot change the country information. Please contact Medit Link's operation team.

Office Phone *
+ 82 01026088546

Mobile Phone (Optional)
+ 82 01026088546

My Scanner

You can check the connected device information for the account.

My Account My Organization My Scanner My Price List									
No.	Type	Serial Number	Usage Time ⓘ	Last Connected	Last Calibrated	Warranty Start Date	Warranty Expiry Da...	Days Left ⓘ	
1	i500	AA1909103895	1d 3h 4m	-	2020-07-06 09:48	2019-09-30	2021-03-30	Purchase -738 Days	
2	Hybrid	1H16091L9003	3d 22h 31m	-	2020-06-24 12:08	2016-09-27	2018-09-26	Purchase -1654 Days	

My Price List (Lab Admin Account Only)

You can enter the products offered by your lab and the price information for each product.

My Account My Organization My Scanner **My Price List**

Teeth Arch

AED (₩)

All (99/557)

Inlay (8/12)

Onlay (8/12)

Coping (10/17)

Denture Coping (5/5)

Crown (15/26)

Pontic (23/36)

Post & Core (3/4)

Veneer (8/11)

Implant Crown (12/53)

Custom Abutment (2/3)

Surgical Guide (1/1)

Offset Substructure (4/5)

Partial Denture (8/7)

Diagnostic Wax-Up (8/1)

Custom Abutment + Implant Crown (0/159)

Custom Abutment + Coping (0/45)

Post & Core + Crown (0/100)

Post & Core + Coping (0/60)

☐	Name	Type	Method	Material	Price
☒	인레이 지르코니아	Inlay	-	Zirconia	\$6.83
☒	인레이 골드	Inlay	-	Gold	1.83
☒	인레이 프레스 세라믹	Inlay	-	Press Ceramic	2.83
☒	인레이 레진	Inlay	-	Resin	3.83
☒	인레이 하이브리드 세라믹	Inlay	-	Hybrid Ceramic	4.83
☒	인레이 글라스 세라믹	Inlay	-	Glass Ceramic	1.00
☒	인레이 컴퍼지트 레진	Inlay	-	Composite Resin	1.00
☒	인레이 e.max	Inlay	-	e.max	1.00

Save

iNote

- Once you contract with a partner, your partner can order your products based on this product list and price.
- If you only have a partnership with your partner, Medit Link does not support orders that require price information.

Other Menu

Notice & Update

You can check the important notices from Medit and program updates of Medit software.

Notice & Update

Notice
NFT Service Termina...
2023-04-19

Notice
End of Support Notic...
2023-04-14

Update
Release Note for Me...
2023-03-29

Notice
[Medit] We've updat...
2023-03-29

Update
Release Note for Me...
2023-03-17

Update
Medit Occlusion Ana...
2023-03-08

Notice
NFT Service Termination Notice
2023-04-19

To our users,

We want to thank all Medit Link users who have supported and used our NFT service.

Unfortunately, we must inform you that this service will be terminated following the upcoming program update.

Please see the details below for more information:

- **Feature in Question**
 - Share via NFT (in Medit Link Case Talk)
- **Effective Date**
 - Starting from the update in the first half of 2023 (scheduled for April 27)

※ The update schedule is subject to change.

We'll continuously work to provide better service.

Medit Link Team

Medit Information

You can get various information about using Medit Link. Access Medit's educational and informational resources, such as Medit Academy, Medit Academy YouTube, and Medit Help Center.

Need more help?

Medit Academy >

YouTube >

Help Center >


Download Medit Link App

You can download the latest version of the program installer from the web. This menu is helpful when you lose the installer provided in the product package, or you have an outdated installer.


47

iNote

Medit Link updates automatically when you have an internet connection.



Overview



1	Title Bar	Refer to Medit Link App > Title Bar for detailed information.
2	Menu	Refer to Medit Link App > Menu for detailed information.
3	Tools and Functions	Refer to Medit Link App > File Viewer and Medit Link App > View Options for detailed information.
4	Case List	Refer to Medit Link App > Filtering and Search for detailed information.

Box Overview

Workflow		Clinic Account	Lab Account
Step 1	Create Cases	Case Box: Create cases and easily search and view the cases with various filtering and search options. With the created case, you can directly link with scanning programs and CAD programs to work on them.	
Step 2	Order/Manage Cases	Order Box: Manage ordered cases. Functions for managing the process from ordering to completing cases are provided.	In Box: Manage cases of received orders. Functions for managing the process of accepting, rejecting, shipping, and completing cases are provided.
Step 3	Work on Scanning/CAD		Work Box: Manage cases that have been forwarded from In Box or that you have created. You can directly link with scanning programs and CAD programs to work on them.

Step 4	Share Cases	Case Talk: Share the created cases with restrictions, without restrictions, with specific people, or other registered members of Medit Link.
Step 5	Supplemental or Additional work	App Box: Utilize a variety of apps that complement your digital dentistry workflow. You can download and install Medit or third-party apps.
Step 6	Delete/Restore Cases	Trash Box: View, restore, or permanently delete cases from the Trash Box.

The table below summarises what Boxes are available in each type of clinic and lab account.

Account Type		Case Box	Order Box	In Box	Work Box	Case Talk	Trash Box	App Box
Clinic Accounts	Admin	O	O	X	X	O	O	O
	Member	O	X	X	X	O	O * Restore feature only	O
Lab Accounts	Admin	X	X	O	O	O	O	O
	Member	X	X	X	O	O	O * Restore feature only	O

Title Bar


The title bar provides quick access to features like network status, cloud synchronization, and notifications.



My Information

My information contains basic user information such as profile image, organization name, email address, and user name.

It also provides a logout option and shortcut to the My Account page in Medit Link Web.



1	Show the profile image.
2	Show the organization name, email address, and user name.
3	Log out of the current account.
4	Go to Medit Link Web while logged in with the current account. You can view, edit, and manage the account information you entered when you signed up.

Notifications

You can see notifications related to contracts, orders, and work processes.

Notification

1

☐ Display Unread Messages Only

Mark All As Read

2

Delivery
Wonhoonchoi_Lab_Dev has dispatched the ordered package.

2020-08-20 오후 5:20

Order
Wonhoonchoi_Lab_Dev has accepted the order.

2020-08-20 오후 5:18

Order
Wonhoonchoi_Lab_Dev has accepted the order.

2020-08-17 오후 4:05

Order
Wonhoonchoi_Lab_Dev has accepted the order.

2020-08-17 오후 4:05

Order
Wonhoonchoi_Lab_Dev has accepted the order.

2020-08-17 오후 4:05

Order
Wonhoonchoi_Lab_Dev has accepted the order.

2020-08-17 오후 4:05

Work
Wonhoonchoi_Clinice_Dev_Member registered a new case test's Case.

2020-08-14 오후 6:33

Work
Wonhoonchoi_Clinice_Dev_Member registered a new case test's Case.

2020-08-14 오전 10:15

Close

- | | |
|---|--|
| 1 | Check the box to show only the unread messages. |
| 2 | Click any notification to instantly move to the related case or webpage. |

Case Synchronization Manager

Case Synchronization Manager helps you to upload and download the cases in Medit Link. You can upload multiple cases to the cloud or download them to your local PC.

[illegible]

Upload

Total Cases / Selected Cases	Check the number of all/selected cases that need to be uploaded.
Upload All	Upload all cases to the cloud.
Cancel All	Cancel all upload requests.
	Upload the case to the cloud.
	Cancel uploading of the case.
	Retry uploading the case.

Download

Total Cases / Selected Cases	Check the number of all/selected cases that need to be downloaded.
Download All	Download all cases to the PC.
Cancel All	Cancel all download requests.
	Download the case to the PC.
	Cancel downloading of the case.

	Retry downloading the case.
---	-----------------------------

Network Status

Medit Link enables you to use all the features seamlessly when you are connected to a network.

To help you check the status of your network connection in real time, the network status icon is displayed in the title bar.

If the network connection is not good, you may experience the following restrictions:

- Some functions, such as ordering a case, creating a case, sharing via Case Talk, and cloud synchronization, will be restricted.
- The overall performance of the service can be reduced due to the impact of the network environment.

In this case, you can enter the Offline Mode to perform basic tasks and then use the features that require internet once you are connected.

iNote

The network symbol indicates the response speed between the current area and AWS (Amazon Web Service). It does not mean your network speed, so sometimes it may show the status as poor or bad, even though your network speed is fast enough.


If the other process uses the same network, the speed may show as poor or bad, and the overall performance of Medit Link could be bad. In this case, please do not share the network with other processes.

For example, the "Upload Raw Data" feature is a huge part of the network process. If the symbol shows poor or bad status, go to Settings > Upload Raw Data and turn "Upload Automatically" off.

Menu

Dashborad


You can view the statistics for the entire work done on Medit Link in visualized graphs and charts – both in the Medit Link App and Web.



INote

Please refer to the [Medit Link Web > Dashboard](#) for more information on the dashboard and Medit Insights.

Case Box/Work Box



1	Show case name and patient name.
2	Provide key features available: Order/Scan(Open)/Order buttons and Apps icons.
3	Show available tabs: Form, CAD, and File Viewer.
4	Show functions for managing cases and files.
5	Show the form registration UI for the selected CAD program.
6	Provide the type/method/material/shade information for users to choose for the selected tooth number.
7	Show the type/method/material/shade information applied to the selected teeth.

Tabs in Case Box/Work Box

Form	CAD	File Viewer
The form tab provides a UI for entering basic information for scanning work, such as type, method, material for teeth numbers, and scan options. The information you enter is passed to the connected scanning program and used as order information when placing an order with your partner.	The CAD tab provides a UI to set CAD options. The information you set up is passed to the connected CAD program.	View result files generated by scanning programs and CAD. The file viewer displays and manages the acquired or attached 3D models and 2D images.


Functions in Case Box/Work Box

The following functions are provided in Case Box/Work Box.

	Case Status	Shows the status of the case.
	Go to Detail	Provides a shortcut for the same case in another box (Case Box ↔ Order Box; In Box ↔ Work Box).
	Case Information	Review case details and save them as a PDF or print them out.
	Memo & Tags	Add a memo or tags when registering or ordering a case. Use them to communicate any additional or supplementary information about the case to your partners.
	Case History	Provides detailed information about the case history. History Last Modified Date 3/29/2022 1:57 PM Scanning Date 4/12/2021 2:06 PM Creation Date 3/22/2022 1:17 PM Case Size 435MB Close
	Delete	Delete the case. You cannot delete ordered cases or cases that were created offline.

	Clone	Create a duplicate of the current case in Case Box. All data will be copied to the cloned case. This feature helps replicate treatment information for the same patient or make an additional order with the same information or files.
	Attach	Attach files in various formats.
	Export	Export scan or CAD result files to the local PC. You can choose your preferred file format when exporting files. The exported files can be used for further work in CAD and CAM software.
	Share	Create a link for sharing. Users with the shared link can access the 3D data in a browser such as Chrome, Firefox, or Safari without installing additional programs. You can share the link freely with anybody or specific people only. A PIN code or a link expiration date can be set to secure private information.
	Submit Support Request	Land on a Medit Help Center page to submit the support request.

Order Box/In Box



1	Show case name and patient name.
2	Provide functions for order management (accepting/canceling; delivering/receiving).
3	Show available tabs: Order, File Viewer
4	Show functions for managing cases and files: case status, shortcuts for detailed information, case information, memo & tags, etc.
5	Show the detailed information for the order including the order number, patient name, partner name, scanning date, order date, delivery date, and memo & tags.
6	Show the ordered product information and delivery status.

Tabs in Order Box/In Box

Order	File Viewer
Check the detailed information related to the order, such as order number, patient name, account name, scanning date, order date, delivery date, memo, and tags.	View files generated by scan/CAD operations. The file viewer displays and manages the acquired or attached 3D models and 2D images.

Case Talk

Cases created in the Medit Link can be shared with specific people you specify, anyone with the shared link or registered members of Medit Link.

The user can see a list of all shared cases through Case Talk in Medit Link App and Medit Link Web.

iNote

Case Talk service requires an internet connection.

MEDIT Link 3.0.0
Wonhoon Choi_clinic

Case Talk

Case Talk Name	Patient Name	Form Information	Scanning Date	Shared Date ↓	Expiration Date
test's Case-262 - Clone_2022-08-05 18:00:05	test	14, 15, 16, 17	2022-07-28 08:27	2022-08-05 18:00	2022-11-03 23:59
Wonhoonchoi_real - Clone_2022-07-25 18:09:10	test1	Maxilla, Mandible	-	2022-07-25 18:09	2022-10-23 23:59
Wonhoonchoi_real - Clone_2022-07-25 18:02:36	test1	Maxilla, Mandible	-	2022-07-25 18:02	2022-10-23 23:59
Multi Occlusion - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Cl_2022-07-21 11:58:18	test	24, 25, 26	2021-05-14 09:27	2022-07-21 11:58	2022-10-19 23:59
Multi Occlusion - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Cl_2022-07-21 11:56:47	test	24, 25, 26	2021-05-14 09:27	2022-07-21 11:56	2022-10-19 23:59
Multi Occlusion - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Clone - Cl_2022-07-21 11:56:47	test	24, 25, 26	2021-05-14 09:27	2022-07-21 11:55	2022-10-19 23:59
test1's Case_2022-07-14 14:45:16	test1	Maxilla	2022-07-06 10:36	2022-07-14 14:45	2022-10-12 23:59


Shared cases can be viewed via a web browser without installing a separate program, enabling easy communication over the 3D data on various devices.

App Box

App Box provides various add-in programs that complement your basic scanning workflow and enhance productivity.

Medit Apps includes Medit-developed original software and integration apps with our partners from different sectors of digital dentistry such as 3D printing, other CAD/CAM software, etc.

Download available apps from Medit Apps and Partners App to create dental appliances, analyze scan data, assist you during a consultation, or streamline your work process.



Trash Box

The Trash Box temporarily stores the deleted cases. Any case in the Trash Box can be either restored or permanently deleted.

⚠ Caution

The data deleted from the Trash Box will be permanently deleted from your PC and the server and cannot be recovered.

Status	Case Name	Patient Name	Lab Name	Deletion Date	Deleted By	Size
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	816 KB
☐ Form Only	wonhoonchoi's Case	wonhoonchoi		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	10,204 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	3,628 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:13	Wonhoon Choi_clinic	0 KB
☐ Form Only	test's Case	test		2022-09-21 12:11	Wonhoon Choi_clinic	4,469 KB

Settings

Please refer to the [Medit Link App > Settings](#) for detailed information on the Settings window.

Patient Management

You can add, edit, delete, or merge patients on the Patient Management window.

Patient Management

Allows to select patients to delete or merge their information from different cases.
Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
Press 'Delete' to remove patients from the system.

All

Patients with No Cases ☐

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date	
☒	5-9_Scan_CAMOrder_Teeth2				3	2020-09-18	2020-09-18	
☒	5-9_Scan_CAMOrder_Teeth2				1	2020-09-18	2020-09-18	
☒	test				2	2020-09-17	2020-09-17	
☒	test				22	2019-07-26	2020-09-17	
☒	test				1	2020-09-17	2020-09-17	
☒	ASD Demo				1	2020-09-16	2020-09-16	
☒	ASD Demo				1	2020-09-16	2020-09-16	
☐	gfm				1	2020-09-15	2020-09-15	
☐	Medit Link 2.3.0			2020-9-24	2	2020-09-14	2020-09-14	
☐	test				1	2020-09-13	2020-09-13	
☐	测试记录				2	2020-09-13	2020-09-13	
☐	test3				1	2020-09-13	2020-09-13	
☐	test1				1	2020-09-13	2020-09-13	
☐	testq1231				26	2019-07-25	2020-09-11	

Add New Patient

1. Click "New Patient."

New Patient

Patient Name*

General

☐ Female ☐ Male

Birth Date

Patient ID

Phone Number

Email Address

Memo

Creating patient information involves consent from the patient.
Patient information can be shared with 3-parties on order.

[Terms & Conditions](#)

2. Fill in the required information in the New Patient dialog, then click "Save."

Delete Patient

iNote

You can delete a patient only when there are no cases for the patient.

1. Select the patient to delete from the list.

Patient Management

Allows to select patients to delete or merge their information from different cases.
 Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
 Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date		
☐	test				0	2019-08-28	2019-09-02		
☐	test				0	2019-08-28	2019-08-28		
☐	test				0	2019-08-28	2019-08-28		
☐	A.J.Abument2				0	2019-08-28	2019-08-28		
☐	Scarbody				0	2019-08-28	2019-08-28		
☐	Margin Line Creation				0	2019-08-28	2019-08-22		
☐	Abument Library				0	2019-08-28	2019-08-22		
☐	GuideMessage				0	2019-08-28	2019-08-16		
☒	margin				0	2019-08-28	2019-08-07		
☒	test2				0	2019-08-28	2019-08-07		
☒	suljadholf				0	2019-08-28	2019-08-07		
☒	Test				0	2019-08-28	2019-08-07		
☒	test				0	2019-08-28	2019-07-26		
☒	test				0	2019-08-28	2019-07-25		

Close

2. Click "Delete."

Merge Patients

You can merge all cases of multiple patients into one patient.

1. Select the patients you want to merge from the list.

Patient Management

Allows to select patients to delete or merge their information from different cases.
 Allows to merge together information of a single patient spread across multiple cases. After merging, all cases of the patient are combined together under the patient name you selected.
 Press 'Delete' to remove patients from the system.

All Search / Contains + New Patient

Delete Merge Patients Patients with No Cases

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date		
☐	5-9 /Scan_CAMOrder_test2				3	2020-09-18	2020-09-18		
☐	5-9 /Scan_CAMOrder_test2				1	2020-09-18	2020-09-18		
☒	test				2	2020-09-17	2020-09-17		
☒	test				22	2019-07-25	2020-09-17		
☒	test				1	2020-09-17	2020-09-17		
☒	ASD Demo				1	2020-09-16	2020-09-16		
☒	ASD Demo				1	2020-09-16	2020-09-16		
☒	gmt				1	2020-09-15	2020-09-15		
☐	Medit Link 2.3.3			2020-9-24	2	2020-09-14	2020-09-14		
☐	test				1	2020-09-13	2020-09-13		
☐	SCX-111111				2	2020-09-13	2020-09-13		
☐	test2				1	2020-09-13	2020-09-13		
☐	test1				1	2020-09-13	2020-09-13		
☐	test21231				26	2019-07-25	2020-09-11		

Close

2. Click "Merge Patients."

3. Select the target patient you would like to merge all cases into.

4. All cases are merged into the target patient and the patient information of the target patient maintains.

Merge all patient information to

Please select the patient you'd like to merge all data to.
After merging all data together, all information in other patient files will be deleted.

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Date Registered	Recent Visit Date
☒	5-9_jScan_CAMOrd...				1	2020-09-18	2020-09-18
☐	test				2	2020-09-17	2020-09-17
☐	test				22	2019-07-26	2020-09-17

Cancel
Apply

5. Click "Apply" to proceed with the merge.

Notice & Update

You can check the important notices from Medit and program updates of Medit software.

Go to Medit Link Web

You will land on the Medit Link Web under the account you are currently logged in with.

Medit Information

You can get various information about using Medit Link. Access Medit's educational and informational resources, such as Medit Academy, Medit Academy YouTube, and Medit Help Center.

Need more help?

×



Medit Academy

>



YouTube

>



Help Center

>

Settings

General

Provides options for overall system setup.

The screenshot shows the 'Settings' window with the 'General' tab selected. The left sidebar lists various settings categories: General, Storage Management, Set Connected Program, Export Options, Case Converting Tool, and Proxy Option. The main area displays the 'General' settings, which include:

- Save Folder:** A text field showing 'C:\MeditWork\Temp\LinkSave\'. There is a folder icon button to the right.
- Language Settings:** A dropdown menu currently set to 'English'.
- Teeth Labeling:** Two radio buttons: 'FDI' (selected) and 'US'.
- Notification:** A checkbox labeled 'Allow Notifications' which is checked.
- Upload Raw Data:** A checkbox labeled 'Upload Automatically' which is checked.
- Quick Scan:** A section with a description: 'Choose if you would like to have the scan stages below automatically appear in Medit Scan for Clinics without having to fill in the form information.' Below this is a checkbox labeled 'Pre-Op for Maxilla, Pre-Op for Mandible, Face Scan' which is unchecked.
- Viewing Options:** A checkbox labeled 'Advanced Rendering' which is checked. Below it is a slider for 'Adjust Color Texture'.
- Medit User Experience Improvement Program:** A section with a description: 'We value your feedback! Join the Medit User Experience Improvement Program and share your thoughts. The collected information will be used to improve the product and services. Medit Link does not collect any information that identifies individuals.' There is a small icon to the left of the text.

At the bottom right of the settings area is an 'Apply' button. At the bottom right of the entire window is a 'Close' button.

Save Folder	Specify a folder to store the data results generated in Medit Link and connected programs.	
Language Settings	Select the language for the program.	
Teeth Labeling	Select a teeth labeling system between FDI and US.	
Notifications	When checked, you will receive notifications for program updates.	
Upload Raw Data	When checked, the raw data is automatically uploaded to the cloud. It helps to back up the original data and perform seamless scanning when using Medit Link with the same account on another PC. The original data cloud backup only applies to data acquired in Medit Scan for Clinics.	
Quick Scan	Choose if you want to have the scan stages below automatically appear in Medit Scan for Clinics without having to fill in the form information. • Pre-op for Maxilla • Pre-op for Mandible • Face Scan	
Viewing Options	Advanced Rendering	When checked, the 3D model is displayed in a more vivid and realistic rendering.
	Adjust Color Texture	Adjust the brightness of the 3D model. * The resulting colors may slightly differ from Medit Link when you open the same data with other software.


Medit User Experience Improvement Program

Join the Medit User Experience Improvement Program.

- We use the collected information to improve our products and services.
- We do not collect any information that identifies individuals.

Storage Management

Check and manage the local and cloud storage space.



Set Period

Set the specific period to delete data created within the below time frame.

- All: See all files.
- 1 Week Ago: See the files that have passed a week after creating.
- 1 Month Ago: See the files that have passed a month after creating.
- 6 Months Ago: See the files that have passed six months after creating.
- 1 Year Ago: See the files that have passed a year after creating.

Local Storage

Manage the PC storage space.

Raw Data	Remove the selected raw data from the PC. If the raw data is already uploaded to the cloud, you can download and use it whenever you want.
Result & Attachment	Remove the result and attachment files from the PC, such as mesh, images, and videos. If these files are already uploaded to the cloud, you can download and use them whenever you want.
Others	See the files used to operate the system and other programs. It shows the storage of your PC. You cannot remove these files from Medit Link.
Delete	Remove selected files from the PC.


Cloud Storage

Manage the cloud storage space.

Raw Data	The raw data does not take up space in the cloud storage, so it will not be in the pi-chart. But you can manage the raw data from here if needed, e.g., delete it to protect the patient's private information, etc.
Result & Attachment	The result and attachment files take up space in the cloud storage. If the cloud is full, you can secure more storage space by removing unnecessary files.
Buy Storage	If you're on a Standard plan, subscribe to the Premium one to expand your cloud storage.
Delete	Remove selected files from the cloud storage to free up storage space.

Set Connected Program

Manage the programs that support integration in Medit Link.



You can connect the following programs:

	Clinic Account	Lab Account
Scan Program	Medit Scan for Clinics	Medit Scan for Clinics Medit Scan for Labs
CAD Program	exocad	exocad

iNote

Integration of the following services is not available for Medit Link for Mac.

- Medit Scan for Labs
- exocad

Export Options

Set how to export files from Medit Link.

Settings

General

Storage Management

Set Connected Program

Export Options

Case Converting Tool

Proxy Option

Export Options

CAD Type

Maestro

Export Files To

C:/Users/Medit-Wonhoonchoi/Desktop/test/

Mesh Format

Convert mesh files to the below formats

☐ medtMesh
☐ OBJ
☒ PLY
☐ STL

☒ Do not convert the attached files

Image Format

Convert images to the below format

☐

☐ JPG
☐ JPEG
☒ PNG
☐ BMP

Compress & Combine

☐ Compress Files
☐ Combine Individual Mesh

Options when same folder or file already exists.

☒ Always ask when exporting
☐ Export a new file and overwrite the existing file. The file currently stored will be lost.
☐ Create a folder and export the files

Rule for Folder Name

Creation Date

Scanning Date

Order Date

Patient Name

Case Name

Partner Name

Patient ID

Order ID

Folder Name

Case Name

Rule for File Name

Creation Date

Scanning Date

Order Date

Patient Name

Case Name

Partner Name

Patient ID

Order ID

File Name

Case Name

Apply

About

User Guide

Close

CAD Type	Select the CAD program you want to use for designing. The program will adjust the results based on your selection, allowing you to proceed with the design without any additional axis alignment.
Export Files To	Specify the path where you want to store the results.
Mesh Format	Specify the format of the final mesh file.
Image Format	Specify the format of the exported image files.
Compress & Combine	Compress Files: Compress all results into a single file. Combine Individual Mesh: Combine all results for the maxilla or mandible into one file. When using some programs that can't import and work with multiple files, you can merge them into one to work in those programs.
Options when the same folder or file already exists.	Choose the default action for cases when the same folder or file name exists in the target path. - Always ask when exporting. - Export a new file and overwrite the existing files. The file currently stored will be lost. - Create a folder and export the files
Rule for Folder Name	Define rules for creating a folder name when exporting files. The set rules will be applied to folder naming, and the exported file will be saved in the created folder. Below information can be used to set the rule. - Clinic account Case Creation Date, Scan Date, Order Date, Patient Name, Case Name, Partner Name, Patient ID, Order ID - Lab account Case Creation Date, Scan Date, Order Date, Patient Name, Case Name, Partner Name, Order ID

Rule for File Name	Define rules for the exported file name. Below information can be used to set the rule. - Clinic account Case Creation Date, Scan Date, Order Date, Patient Name, Case Name, Partner Name, Patient ID, Order ID - Lab account Case Creation Date, Scan Date, Order Date, Patient Name, Case Name, Partner Name, Order ID
--------------------	--

Case Converting Tool

Move cases between different accounts.

The screenshot displays the 'Case Converting Tool' interface. On the left is a sidebar with 'Settings' as the main header and several sub-items: General, Storage Management, Set Connected Program, Export Options, Case Converting Tool (highlighted), and Proxy Option. At the bottom of the sidebar are 'About' and 'User Guide' links.

The main area is divided into two panels: 'Export' and 'Import'.

Export Panel:

- At the top, there's a date range selector set to '2020-09-13 ~ 2020-09-20' and a dropdown menu set to 'All'. A search bar with the placeholder 'Search / Contains' is also present.
- Below this is a table titled 'Files to Export' with columns: Case Name, Patient Name, and Date Modified. It lists several files, some with checkboxes for selection.
- At the bottom of the export panel, there is a 'File Name' field containing 'CaseBox_2020-9-20_123057' and an 'Export Files To' field containing 'C:/Users/Medit-Wonhoonchoi/Documents/'.
- An 'Export' button is located at the bottom right of the export panel.


Import Panel:

- A note at the top states: 'A new case is created for every imported case. It contains case data, patient information and scan file, but not the order information.'
- Below the note is an 'Import from' field containing 'C:/Users/Medit-Wonhoonchoi/Documents/' and an 'Import' button.
- At the bottom is a table titled 'File List' with columns: File Name and Import Time. It lists several imported files and their corresponding import times.


A 'Close' button is located at the bottom right of the entire interface.

How to Use Case Converting Tool

1. Open Medit Link and click "Settings."




2. Choose the "Case Converting Tool" section in the left-side menu.



3. In the import section, click the folder icon to browse your PC for the .meditex file you want to import.


4. Click the "Import" button and check the imported file in the file list below.



- The imported case will also appear in your Case Box (Clinic account) or Work Box (Lab account).


Proxy Options

In some limited network environments, you can change network settings.



About

Show detailed information about Medit Link.



User Guide

Get the user guide for Medit Link.

Opening the user guide requires an Internet connection.

File Viewer

File Viewer groups and shows various data, including the data acquired from programs linked to Medit Link and data attached by users.


It provides basic management functions such as showing/hiding data, adjusting opacity, and deleting data.

It supports previews of files in the following formats:

- 3D data: STL, OBJ, PLY
- 2D images: PNG, JPG, JPEG, BMP, GIF, TIF, TIFF
- Video files: MKV, MP4, AVI, WMV

Data Tree

You can group and display all data for the case in the data tree. You can also show/hide data, adjust transparency, and delete data from the data tree.

	Show/Hide	You can show or hide the data.
	Opacity	The slide bar appears when you click one of the data on the data tree. You can adjust the opacity of the shown data by moving the slider bar. 
	More	Clicking the More icon shows additional features available.  • Show This Only: Hide all other data but display only the corresponding data. • Hide: Hide the data. • Duplicate: Copy the data and add it to the attachment group. • Export: Export the selected files. • Rename: Change the name of the data. • Delete: Delete the data.





Note

You can rename data when online.

Model Control Tools

3D Model Control

	Pan	Move the model.
	Rotate	Rotate the model.

	Zoom In/Out		Zoom in and out on the model.
	Zoom Fit		Position the model in the center of the screen.
	Model Display Mode	Texture On	Show data in real texture colors. 
		Texture Off	Show data in a single color. 
	Grid Settings (mm)	Grid On	Show the grid on the background. 
		Overlay On	Overlay the grid over the model. 
		Grid Off	Hide the grid.

2D Image Control


	Next/Previous Image	Move among images.
	Edit	Edit the image.
	Pinch/Paste Window	You can detach the image window and adjust it to the desired position. You can resize the window. When the window is restored to its original location, it returns to its default size.

Function Tools

Image Capture and Memo


Capture the current screen as a 2D image, draw on it, or add notes.

You can save it as a file with your notes and send it as a reference.



Ortho Capture


Capture the images in the current viewing direction and top/bottom/left/right from the current viewing direction.



Sync View	When on, you can rotate or zoom in/out the center image. The other images will be adjusted simultaneously as fixed viewpoints as you move the center image. When off, you can pan, rotate, or zoom in/out all images individually.
Save Each Image Separately	Save each image in high resolution.
Show Texture	Turn on/off the texture for the model.

Reassign Data


Change the data for maxilla/mandible with the Reassign Data feature. You can use the changed data to create and utilize the result images of Ortho Capture.



- **Scan Data:** Scan data for the case are all shown.
- **Maxilla:** The data you selected for the maxilla are displayed.
- **Mandible:** The data you selected for the mandible are displayed.


Fill Holes

Fill holes in 3D mesh data.



Trimming


Trimming provides various tools to edit data.



The selected 3D mesh data can be edited using different tools located at the bottom of the window.

Selection		Smart Teeth Selection	Automatically select all teeth of the arch leaving out gingiva parts. iNote This function is only available for the scan data that has been acquired by Medit Scan for Clinics with the "Use GPU" option on.
		Smart Single Tooth Selection	Automatically select the area of a single tooth leaving out gingiva parts. Click, press, and drag the mouse on the tooth.
		Polyline Selection	Select all entities within a polyline shape drawn on the screen.
		Brush Selection	Select all entities on a freehand-drawn path on the screen. Only the front face will be selected. The brush comes in three different sizes.
Adjustment of Selected Area		Shrink Selected Area	Reduce the selected area each time you press the button.
		Expand Selected Area	Expand the selected area each time you press the button.
		Invert Selected Area	Invert the selection.
Deselection		Selection/Deselection	When on, deselect the area using various tools.
		Clear All Selection	Clear all selected areas.
		Delete Selected Area	Delete the data from the selected area.

1. Select the area you want to trim with the area selection tools.




iNote


The Polyline trimming tool selects all entities within a polyline shape drawn on the screen, while the Brush trimming tool selects only the front face data.

2. You can invert, shrink or expand the selection with the following three tools.


- Invert Selected Area



- Shrink Selected Area




- Expand Selected Area



3. You can deselect the selected areas with the Delete Selected Area tool.





Optimize Occlusion Alignment

Optimize the occlusion alignment. You can choose the target data and adjust the occlusion alignment.




Compensation Level



Mid	Compensation Level  A slider control with a blue circle positioned approximately one-third of the way from the left end. The slider is flanked by two small icons of a tooth with a bracket.	 A 3D model of a lower dental arch (mandible) showing a dark grey base with several circular indentations representing teeth. The model is shown from a slightly elevated, front-facing perspective.
High	Compensation Level  A slider control with a blue circle positioned approximately three-quarters of the way from the left end. The slider is flanked by two small icons of a tooth with a bracket.	 A 3D model of a lower dental arch (mandible) showing a dark grey base with several circular indentations representing teeth. The model is shown from a slightly elevated, front-facing perspective.
Off	Compensation Level  A slider control with a blue circle positioned at the far left end. The slider is flanked by two small icons of a tooth with a bracket.	 A 3D model of a lower dental arch (mandible) showing a dark grey base with several circular indentations representing teeth. The model is shown from a slightly elevated, front-facing perspective.

Mandibular Movement

You can play the mandibular movements recorded in Medit Scan for Clinics to see the dynamic occlusal contacts between the upper and lower dentitions.



	Free Movement	Play the recording of free mandibular movement.
	Left Lateral Movement	Play the recording of the mandible moving to the left.
	Right Lateral Movement	Play the recording of the mandible moving to the right.
	Protrusive Movement	Play the recording of the mandible moving forward.

Manage Multi Occlusion Data in File Viewer

You can get multiple results of occlusion data from Medit Scan Application.

Medit Link provides the function to manage the results through the "File Viewer."

How to manage multiple results


When the user gets multiple occlusion results from Medit Scan for Clinics, Medit Link provides the option to select a result.

The list of results is listed by name, and you can select one item.

The selected model will be shown in the "Data Tree" and "Model Display Area."



You can collapse the list and change the item using the arrow icons.



View Options

You can change the View Options to customize the layouts of the case list in the Case Box/Work Box/In Box/Order Box.

View Options

 Status Filters

6 in use >

 Table Items

6 in use >

 View Type

List View >

 Grouping

None >

 Sorting

Last Modified Date >

Status Filters

Show or hide each filter representing the status of the case. The hidden status filters are not displayed on the screen in the Case Box/Work Box/In Box/Order Box.

You can easily find the case you need with the shown filters.

View Options

< Status Filters



Shown

Hide All

Form Only



Need Processing



Scan Completed



Hidden

Show All

CAD



Ordered




Case Completed



Cancel

Apply



- You can click the status you want to view to display only cases in that state.
- You can select multiple statuses at a time.

Available filtering options for each box:

Case Box	All/Form Only/Need Processing/Scan Completed/CAD/Ordered/Case Completed
Order Box	All/Pending/Canceled/Rejected/Accepted/Shipped/Case Completed
In Box	All/Pending/Canceled/Rejected/Accepted/Ready/Shipped/Case Completed
Work Box	All/Form Only/Need Processing/Scan Completed/CAD/Done
Trash Box	Not supported

Table Items

Select the items to display in the case list in the Case Box/Work Box/In Box/Order Box. You can change the order of table items or adjust the table column width in the list view to best optimize for the user environment.

View Options

< Table Items



Shown

Hide All

- ⋮ Status
- ⋮ Case Name
- ⋮ Patient Name
- ⋮ Form Information
- ⋮ Last Modified Date
- ⋮ Scanning Date

Hidden

Show All

- ⋮ Birth Date
- ⋮ Order Date
- ⋮ Lab Name

Cancel

Apply

MEDIT Link 3.0.0 Wonhoon Choi_clinic

Case Box

☒ All
 ☐ Form Only
 ☐ Need Processing
 ☐ Scan Completed
 ☐ CAD
 ☐ Ordered
 ☐ Case Completed

Total Count: 514 2018-01-01-2022-09-29 All / Contains

Status	Case Name	Patient Name	Form Information	Last Modified Date	Scanning Date
Form Only	test's Case	test	-	2022-09-29 오전 8:45	
Ordered	test's Case	test	14-Implant Crown / 15-Implant Crown / 16-Implant Crown	2022-09-29 오전 8:13	2022-09-29 오전 7:56
Form Only	test's Case	test	-	2022-09-28 오후 5:56	
Form Only	test's Case	test	-	2022-09-28 오전 10:41	
Scan Completed	test's Case	test	24-Inlay / 25-Inlay / 26-Inlay	2022-09-28 오전 9:27	
Scan Completed	Multi Occlusion - Clone - Clone - Clone ...	test	24-Inlay / 25-Inlay / 26-Inlay	2022-09-27 오후 7:03	2021-05-14 오전 9:27
Scan Completed	Preop 사용된 케이스 - Clone	straumann	15-Inlay / 26-Pontic	2022-09-27 오후 5:09	2021-11-18 오후 6:53

View Type

Change the view type between the Gallery View and List View on the Case Box.

MEDIT Link
3.0.2
Wonhoon Choi_clinic

All

Form Only

Need Processing

Scan Completed

CAD

Ordered

Case Completed

Total Count: 644

2018-01-01~2022-10-28

All / Contains

Dashboard	Case Box	Order Box	Case Talk
test	temporaries - Clone - Clone test Crown 25, 27 Pontic 26 Scan Completed	Preop 사용한 케이스 straumann Inlay 15 Pontic 26 Ordered	test's Case test Inlay 28 Crown 16 Custom Abutment + 21, 22, 23, 24, 25, 26 Scan Completed
No Registered Form	주문테스트 1 - Clone test Inlay 16, 15, 14, 23, 24, 2... Ordered	temporaries - Clone - Clone test Crown 25, 27 Pontic 26 Scan Completed	R2.4.0_TAK_Medit Scan for Clinics_Basic_with ... 2.4.0 Data_TAK Inlay 46 Crown 25, 27 Implant Crown 16, 36 2More Scan Completed
	temporaries	주문테스트 1 - Clone - Clone	test's Case

1

Medit Link2.3.0 Release

2

Medit Link 2.3.0

2020-09-24

3



Implant Crown Zirconia

11, 12, 13, 14,...



Inlay Zirconia

3, 4, 5, 6, 7

+

1

4



Crown Metal

18, 19, 20, 21,...

5

Scan Completed

2020-09-14 오전 7:43

84

List View

Dashboard

Case Box

Order Box

Case Talk

App Box

Trash Box

Case Box

Order Box

Case Talk

App Box

Trash Box

MEDIT Link3.0.2

Wonhoon Choi_clinic

+ New Case

AllForm OnlyNeed ProcessingScan CompletedCADOrderedCase Completed

Total Count: 644

2018-01-01~2022-10-28AllAll / Contains

Status	Case Name	Patient Name	Form Information	Last Modified Date	Scanning Date	Lab Name
Scan Completed	test's Case	test	-	2022-10-28 오전 7:32		
Scan Completed	temporaries - Clone - Clone - Clone	test	25-Crown / 26-Pontic / 27-Crown	2022-10-27 오후 9:29		
Form Only	test's Case	test	-	2022-10-27 오후 8:45		
Form Only	test's Case	test	-	2022-10-27 오후 8:39		
Scan Completed	test's Case	test	12-Inlay / 13-Inlay / 14-Inlay / 15-Inlay	2022-10-27 오후 8:31		
Form Only	test's Case	test	-	2022-10-27 오후 8:30		
Form Only	test's Case	test	-	2022-10-27 오후 7:33		
Form Only	test's Case	test	-	2022-10-27 오후 7:07		
Ordered	test's Case	test	13-Onlay / 14-Onlay / 15-Onlay	2022-10-27 오후 3:50		Wonhoonchoi-Lab
Scan Completed	케이스해모있고주문해모수경한케이스 - Clone	test	14-Inlay / 15-Inlay / 16-Inlay	2022-10-27 오전 11:28	2021-12-03 오후 2:37	
Scan Completed	파일리스트레스12 - Clone	tst	13-Inlay / 14-Inlay / 15-Inlay / 16-Inlay / 17-In...	2022-10-27 오전 10:53	2022-10-11 오후 8:12	
Scan Completed	test's Case	test	-	2022-10-27 오전 9:10		
Scan Completed	test's Case	test	-	2022-10-27 오전 8:59	2022-10-27 오전 8:54	

Grouping

Select a grouping option to classify and rearrange cases in the list by the selected option.

View Options

Status Filters3 in use >

Table Items6 in use >

View TypeList View >

GroupingNone v

None

Patient Name

Order Date

Scanning Date

SortingLast Modified Date >

For each box, the following grouping options are available:

Box	Grouping
Case Box	None
	Patient Name
	Scan Date
	Order Date
Order Box	None
	Scan Date
	Order Date

85

	Requested Delivery Date
In Box	None
	Order Date
	Requested Delivery Date
Work Box	None
	Requested Delivery Date

iNote

The grouping and sorting options are unavailable in the Trash Box and Case Talk.

Sorting

Select a sorting option to change the order of cases in the list according to the selected option.

View Options

Status Filters 3 in use >

Table Items 6 in use >

View Type List View >

Grouping None >

Sorting Last Modified Date ▾

☐ Patient Name
☒ Last Modified Date
☐ Scanning Date
☐ Order Date
☐ Lab Name

For each box, the following sorting options are available:

Box	Grouping	Sorting
Case Box	None	Patient Name / Last Modified Date / Scan Date/ Order Date / Lab Name
	Patient Name	Last Modified Date / Scanning Date/ Order Date / Lab Name
	Scan Date	Patient Name / Last Modified Date / Order Date / Lab Name
	Order Date	Patient Name / Last Modified Date / Scan Date/ Lab Name
Order Box	None	Patient Name / Scan Date / Order Date/ Requested Delivery Date / Lab Name
	Scan Date	Patient Name / Order Date / Requested Delivery Date / Lab Name
	Order Date	Patient Name / Scan Date / Requested Delivery Date / Lab Name

	Requested Delivery Date	Patient Name / Scan Date / Order Date / Lab Name
In Box	None	Patient Name / Order Date / Requested Delivery Date / Clinic Name
	Order Date	Patient Name / Scan Date / Requested Delivery Date / Clinic Name
	Requested Delivery Date	Patient Name / Order Date / Clinic Name
Work Box	None	Patient Name / Last Modified Date / Requested Delivery Date / Clinic Name
	Requested Delivery Date	Patient Name / Last Modified Date / Clinic Name

iNote

The grouping and sorting options are unavailable in the Trash Box and Case Talk.

Filtering and Search


Filter by Date

Filter and display cases last modified within the set duration.

1. Click the duration filter.

A date range filter button with a calendar icon and the text "2020-08-10 ~ 2020-09-09".

2. Select the duration from the list.


A duration filter interface. At the top is a dropdown menu with "Last 30 Days" selected. Below it are two date range inputs. The "From" input has a calendar icon and the date "2020-08-21". The "To" input has a calendar icon and the date "2020-09-20".

- All: Shows all cases.
- Custom Range: Specify the start and end dates to display cases last modified within the period.
- Today: Show the cases last modified today.
- Yesterday: Show the cases last modified yesterday.
- Past 7 days: Show the cases last modified during the previous 7 days excluding today.
- Past 30 days: Show the cases last modified during the previous 30 days excluding today.
- This Month: Show the cases last modified during this month.
- Previous Month: Show the cases modified during the last month.

The case list will be updated in real time with the cases created or last modified within the selected duration.

Advanced Search Capabilities

Leverage search capabilities to quickly find the information you want.

A search bar with a dropdown menu showing "All", a search icon, and the text "Search / Contains".

The following Searching Criteria are available for each box:

Case Box	All, Case Name, Patient Name, Patient ID, Form Information, Lab Name
Order Box	All, Case Name, Patient Name, Patient ID, Lab Name, Order ID
In Box	All, Case Name, Patient Name, Clinic Name, Order ID
Work Box	All, Case Name, Patient Name, Form Information, Clinic Name, Order ID
Case Talk	Case Talk Name (Contains), Patient Name (Equals)
Trash Box (Clinic)	All, Case Name, Patient Name, Lab Name, Deleted By

Trash Box
(Lab)

All, Case Name, Patient Name, Clinic Name, Deleted By

You can select the Searching Criteria and select between the Contains and Equals logic.

All

▼ All / Contains

🔍

Searching Criteria

☒ All

☐ Case Name

☐ Patient Name

☐ Patient ID

☐ Form Information

☐ Lab Name

☒ Contains


☐ Equals

Create/Edit/Clone/Delete Case

Creating Cases

Create Case in Case Box (Clinic)

1. Click the "New Case" button in the upper right corner.




The form is titled "New Case Information". It contains the following fields and elements:

- Patient Name***: A text input field with a "Patient List" icon (three horizontal lines with a magnifying glass) to its right.
- Case Name***: A text input field.
- Information**: A small icon with an 'i' followed by the text: "Creating patient information involves consent from the patient. Patient information can be shared with 3rd parties while placing an order."
- Tag**: A text area with the instruction: "You can add the entered tag with Enter or Spacebar."
- Footer**: A note stating "You can add up to 10 tags of 100 characters or less." followed by a link "Terms & Conditions". At the bottom are three buttons: "Cancel", "Register", and "Register & Scan".

2. Enter basic information, such as the patient's name and case name.
 - If you enter a name not on the patient list, a new patient will be added.
 - If you enter a name that is in the patient list, you can select it from the list.
 - If registering a case for an existing patient, click the Patient List icon to select an existing patient.



- In the Patient List dialog, you can register a new patient or edit the information of an existing patient.



The "Patient List" dialog shows a table of existing patients. At the top, there is a search bar with "All" selected and a "Search / Contents" label. A "+ New Patient" button is in the top right. The table has columns: Patient Name, Patient ID, Gender, Birth Date, Number of Cases, Registered Date, and Recent Visit Date. Each row has a radio button in the first column and an edit icon in the last column.

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Registered Date	Recent Visit Date	
☐	test				0	2020-08-28	2020-08-28	
☐	test				0	2020-08-28	2020-08-28	
☐	test				13	2020-08-23	2020-08-28	
☒	test				85	2020-08-23	2020-05-08	
☐	Worhosen Chis				1	2020-04-01	2020-04-01	
☐	0000				0	2020-08-18		
☐	13				0	2020-08-20		
☐	14				0	2020-08-20		
☐	13				0	2020-08-20		

At the bottom right are "Close" and "Select" buttons.

3. The case name will be filled in automatically, but you can still change it.
4. Click "Register & Scan" to complete the case registration and proceed to scan.
 - **Register**: You will be redirected to the main page. This is useful when registering multiple cases.

- **Register & Scan:** Move to the Form tab to register form information for the case.

Create Case in Work Box (Lab)

1. Click the "New Case" button in the upper right corner.

New Case Information

Clinic Name*

Patient Name*

Case Name*

Requested Delivery Date

2023-04-13

Tag

You can add the entered tag with Enter or Spacebar.

You can add up to 10 tags of 100 characters or less.

Cancel Register

2. Enter or select the clinic name.
3. Enter the patient's name and case name. The case name will be filled in automatically, but you can still change it.
4. Enter the expected delivery date considering the completion date of the work.
5. After registering the case, you will be redirected to the main page.

Editing Cases

Edit Case in Case Box (Clinic)

1. Click the "Edit" icon next to the case name.

Edit Case Information

Patient Name*

Test Occlusion

Test Occlusion 00010101 is selected.

Case Name*

mandibular movement - Clone - Clone

Creating patient information involves consent from the patient. Patient information can be shared with 3rd parties while placing an order.

Tag

You can add the entered tag with Enter or Spacebar.

You can add up to 10 tags of 100 characters or less.

[Terms & Conditions](#)

Cancel Apply

2. Click the "Search" icon to change the patient's name. Note that it cannot be edited directly.; instead, you can find a patient name you would like to replace the current one with and can edit the patient information from the list.

Patient List

Select a patient in the list to set the patient of the case. You can add or manage the patient information in Medit Link.

All Search / Contains  [+ New Patient](#)

	Patient Name	Patient ID	Gender	Birth Date	Number of Cases	Registered Date	Recent Visit Date	
☐	test				0	2020-08-28	2020-08-28	
☐	test				0	2020-08-28	2020-08-28	
☒	test				13	2020-08-23	2020-08-28	
☐	test				85	2020-08-23	2020-05-08	
☐	Wanhooon Choi				1	2020-04-01	2020-04-01	
☐	ttttt				0	2020-08-18		
☒	13				0	2020-08-20		
☐	14				0	2020-08-20		
☐	15				0	2020-08-20		

[Close](#) [Select](#)

3. Check the patient's name you want to replace the current one with.
4. To create a new patient, press the "New Patient" button.
5. Use the "Edit" button to change the registered information about the patient.
6. Click "Apply" to finish the process.

Edit Case in Work Box (Lab)

1. Click the "Edit" icon next to the case name.

Edit Case Information

Clinic Name*

test

Patient Name*

test

Case Name*


test's Case

Requested Delivery Date

 2022-04-18

Tag

You can add the entered tag with Enter or Spacebar.



● You can add up to 10 tags of 100 characters or less.

[Cancel](#) [Apply](#)

2. Click "Apply" after the editing is done.

iNote

You cannot edit the information of the cases received through the order in In Box.

Cloning Cases

Clone Case in Case Box (Clinic)

1. Click the "Clone" icon at the top.



2. You will be asked what kind of information you would like to clone to the new case.

Information
 Select a clone type.

[Cancel](#) [Form Data Only](#) [All Data](#)

3. Choose whether to replicate only form information only or all data.
4. When cloning is complete, a case named "[Case name] – Clone" will be created.

Clone Case in Work Box (Lab)

1. Click the "Clone" icon at the top.



2. Duplicate all data to create a new case.
3. When cloning is complete, a case named "[Case name] – Clone" will be created.

Deleting Case

Delete Case (Clinic and Lab)

1. Click the "Delete" icon.



2. Click "OK" to delete the case.

Warning
 Do you want to delete the case?

[Cancel](#) [OK](#)

3. Deleted cases are moved to the Trash Box and can be recovered if necessary.

Form Information

You can enter the form information based on the connected program and the treatment type.

The information you enter will be used as basic information for working on scanning or CAD programs and ordering.

Registering Form Information in Case Box

Enter on the Form Tab

1. Select whether you are going to input treatment information by teeth or by arch.
2. Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the order list.
3. Click "Save" to save your registration.

Enter on the CAD tab

1. Select one of the following options that will define your scanning strategy: Teeth/Flexible Multi-Die/Multi-Die.
2. Select the tooth number/arch/die number.
3. Define the settings for detailed options and CAD operations.
4. Click "OK" to save the settings.
 - The selected information is automatically registered to the order list.
5. Click "Save" to save your registration

Registering Form Information in Work Box

Enter on the Form Tab

1. Select the scan program you want to use.
2. Select one of the following options that will define your scanning strategy: Teeth/Flexible Multi-Die/Multi-Die.
3. Methods may differ depending on the selected scan program.

	Medit Scan for Clinics	Medit Scan for Labs
Teeth	☐	☐
Arch	☐	☐
Flexible Multi-Die	☒	☐
Multi-Die	☒	☐


4. Select teeth, arch, or die numbers for the scanning process.
5. Select Type – Method – Material – Shade.
 - If you entered information incorrectly, you go back to select it again.
 - The selected information is automatically registered to the list.
6. Click "Save" to save your registration.

Enter on the CAD tab

1. Select one of the following options that will define your scanning strategy: Teeth/ Flexible Multi-Die/Multi-Die.
2. Select the tooth number/arch/die number.
3. Define the settings for detailed options and CAD operations.
4. Click "OK" to save the settings.
 - The selected information is automatically registered to the list.
5. Click "Save."

Connecting Teeth into Bridge

You can connect multiple teeth into a bridge when the same material is applied to multiple adjacent teeth. Click the bridge icon beside the selected teeth to connect or disconnect them.



Setting Scanning Options

Depending on the program you are working with, various scan options are provided.

The selected scan options are passed to each scanning program.

Medit Scan for Labs

General options

	Adjacent Teeth	Mark the prepared and adjacent teeth and pass this information to the scanning program.
	Antagonist	Automatically mark the antagonist and passes this information to the scanning program.
	Articulator	Display options for selecting the articulator in the scanning program.
	Movement Marker	Allow using the "Jaw Motion" scan strategy in the scanning program.

Options by teeth

	Scan Body	Allow setting a strategy for scanning scan bodies in the scanning program.
	Pre-Op Scan	Allow setting a strategy for scanning a pre-op model in the scanning program.
	Wax-Up Scan	Allow setting a strategy for scanning wax-ups in the scanning program.

Adding to Favorite

You can quickly select a product by registering it as a favorite.

The list of registered favorites is managed as a template and can be shared through the export and import features.

Add to Favorites

1. Select the tooth number.
2. Select Type – Method – Material – Shade.
 - Type – Method – Material – Shades can be registered as favorites even if not all are selected.
3. Click the "Star" icon.

4. A combination of the selected parameters, excluding tooth number, is on your favorite search list.

Manage Favorite Templates

You can register, manage, and share different types of templates.

Add	Add a new template.
Rename	Rename the current template.
Delete	Delete the current template.
Export	Export the current template to a file.
Import	Import a template file.

Work on Scanning/CAD

Scan Process

Based on the entered information, you proceed with scanning.

1. Enter/save treatment information and select the necessary scan options.
2. Click the "Scan" button in the upper-right corner of the window.

CAD Process

Medit Link is integrated with exocad.

CAD	Registering treatment information	Confirmation of treatment information	Direct CAD execution	CAD results import back to Medit Link
exocad	O	O	O	O

1. In Settings > Set Connected Program, you can link the CAD program of your choice from the list.
2. Enter/save task information and design options in the CAD Tab.
3. Click "CAD" in the upper-right corner of the window.
4. The linked CAD program is run, and the data is sent automatically.
5. Save the results after performing CAD.

iNote

Medit Link for Mac does not support the function to link with CAD.

Get Result Files

Export Scan Result Files


In addition to sharing results with partners through orders, Medit Link allows you to store your results in the desired file format using the export function.

Export Files To	Specify the path where you want to store the results.
Folder Name	The system creates a new folder and stores all results in it. The default folder name is set by the naming rule. But you can change the folder name directly. You can define the naming rule on Settings > Export Options > Rules for Folder Name
File Name	Set the file name for the result. The naming rule sets the default file name. But you can change the file name directly. You can define the naming rule through Settings > Export Options > Rules for File Name
CAD Type	Select the CAD program to be used for your design program. Medit Link will adjust the files based on this selection, allowing you to proceed with the design without any additional axis alignment.
Mesh Format	Specify the format of the mesh file.
Image Format	Specify the format of the exported image files.
Compress Files	Compress all results into a single file.
Combine Individual Mesh	Combine the result of the maxilla or mandible with the individual tooth. Useful when importing several files into a program is not possible.
Export All Occlusal Relationship	When multiple occlusal results are generated from the scan app, you can export all occlusal results at once.
Save Case Information as PDF	Save the case information as a pdf with the result.


Export Selected Files

Medit Link provides an export function from the expanded menu in the file viewer. Use this, and you can export the selected files easily.

1. Click the "Menu" icon next to the data group or an entity in the data tree of the file viewer.
2. Then click "Export."



3. The export dialog shows only the files selected that you want to export.



4. Finally, click "Export."


Export Margin Line

Medit Scan for Clinics provides a function to create the margin line, and Medit Link can transfer the margin line to the CAD program.

Export Margin Line as a Color

In general, Medit Link converts the margin line as color and exports the result with maxilla or mandible together. You can create the margin line in their cad referencing the color information.

Example:



iNote
This feature is available for formats that include color information such as meditMesh, OBJ, and PLY. The margin line does not convert into STL.

iNote
CAD type > 3shape does not support this feature.


Export Margin Line as a File

In the case of Exocad, Medit Link export the margin line to "~margin.xyz" which can be imported into exocad directly.

You can use this file in exocad without additional work.

Example:

Maxilla	Margin Line
---------	-------------



Export Case List as CSV

Medit Link provides the function to export the case list in CSV format.

You can export the case list with basic information and manage cases easily.

The "Export to CSV" function is available in the following menu of the Medit Link application.

- In Box
- Work Box

How to export the case list

1. Go to In Box or Work Box.
2. Select the cases that you want to export.

MEDIT Link 3.0.0 Wonhoon Choi

In Box

✓ All Pending Accepted Case Completed View Options

Cancel Selection Delete Items Selected: 5 2018-01-01~2022-09-29 All / Contains

	Status	Case Name	Patient Name	Order Date	Requested Delivery Date	Clinic Name	Order ID	Price
☐	Pending	test's Case	test	2022-09-29 오전 8:13	2022-10-01 오전 10:30	Wonhoonchai-Clinic	10413603	
☐	Pending	test's Case - Clone	test	2022-09-22 오후 7:03	2022-09-24 오전 10:30	Wonhoonchai-Clinic	10374886	
☐	Pending	test's Case - Clone	test	2022-09-22 오후 7:01	2022-09-24 오전 10:30	Wonhoonchai-Clinic	10374864	
☐	Accepted	test's Case	test	2022-09-13 오후 2:43	2022-09-15 오전 10:30	Wonhoonchai-Clinic	10314293	
☒	Pending	test's Case	test	2022-09-05 오전 7:31	2022-09-07 오전 10:30	Wonhoonchai-Clinic	10268220	
☐	Case Completed	Worked with 2.3.6 - Clone - Clone - Clone - Clone	test	2021-07-01 오전 8:51	2021-07-08 오전 8:50	Wonhoonchai-Clinic	8579990	500 KRW
☒	Case Completed	test's Case - Clone - Clone - Clone - Clone - Clone	test	2021-08-31 오후 4:10	2021-09-07 오후 4:10	Wonhoonchai-Clinic	8730156	2,000 KRW
☒	Case Completed	test's Case memo	test	2022-04-08 오후 11:...	2022-04-15 오후 7:45	Wonhoonchai-Clinic	9569817	

3. Click the "Expanded menu" button at the top to select the "Export to CSV" option.



4. Define the path and file name then save the results.

Download Files from Medit Link Web

Medit Link provides a download function.

If you want to get the result files without installing the Medit Link application, you can download the files from the Medit Link web service.

iNote


The data should be uploaded on the Medit Cloud to download the files from the web.

The download function is available in the following menu of the Medit Link web service.

- Case Box
- Work Box
- Web Viewer

Case Box/Work Box

You can find the "Download" button at the top right.



Web Viewer

When you share the case using Case Talk, you can set whether to allow downloading or not.

If you check "Grant Download Permission" and share the case, all users who connect to the shared link can download the files.

Share through Case Talk

Share with Anyone

1 Create a shared link and share it with anyone.

Expiration Date: 2022-04-27

1 You can extend the period up to 90 days from today.

PIN Code: ☐ PIN Code

1 Set the PIN code using six digits for privacy protection.

☒ Share the patient's name

☒ **Grant Download Permission** *Premium Feature*

1 If you check this option, all users who have access to the shared link can download the files.

1 The premium feature is provided to all users for a limited time. Feature availability is subject to change depending on the membership plan.

☒ Allow comments

1 If you check this option, you will be able to communicate through the shared link.


Comment: Please enter your comment.

0/1024

By pressing the "Get Link" button, you agree that you understand personal information can be collected and provided to a third party (or shared in any public domain) ONLY with the consent of the person whose images on which you have worked. You also agree the consent of that person was truthfully obtained, and Medit is not responsible for any problems which may arise there after.

Close **Get Link**

In the Web Viewer, you can find the "Download" icon at the top right.



How to download

1. Click "Download."



2. Select the mesh format from the download dialog.

Download Premium Feature

The premium feature is provided to all users for a limited time.
Feature availability is subject to change depending on the membership plan.

Mesh Format

Convert mesh files to the below formats

☐ meditMesh ☐ OBJ ☒ PLY ☐ STL

i The selected format will be applied to the result files created in the Medit Link application.

i The attachment files will not be converted but will keep their original format.

Close **OK**

3. Then click "OK."

Order/Mange Orders

You can easily order a case with data acquired from connected programs and attached by the user.

Ordering Cases

You can place an order for your partners in Case Box.

1. Click the "Order" button.
2. Check the information required for your order.

The screenshot shows the 'Order' form with the following fields and options:

- Order** (Title)
- General** (Selected tab), **Order List**, **File List**
- Partner**: WonhoonChoi_Lab_Stage
- Partner Information**: WonhoonChoi_Lab_Stage, Seoul, KR, [02835] 23 Incheon-ro 22-gil, Seongbuk-gu, Seoul, Korea, 01026088546
- Case Name**: mandibular movement - Clone - Clone
- Patient Name**: Test Occlusion
- ☒ **Share patient information**
- Requested Delivery Date**: 2023-04-13
- Delivery Time**: 오전 8:15 (Range: AM 8:00 to PM 8:00)
- Memo**: (Empty text area)
- Tags**: (Empty text area with placeholder: "You can add the entered tag with Enter or Spacebar.")
- You can add up to 10 tags of 100 characters or less.
- ☒ **Include all occlusion relationships in the order**
- ☒ **Include Tags In The Order**
- Buttons**: Cancel, OK

3. Select a partner.
4. Check the product information from the Order List tab.
 - You can decide whether to share the patient's name.
5. Check the files to order from the File List tab.
6. Enter your desired delivery date.
 - The default value is set to 7 days after the order date.
7. You can add a memo and tags to your partner to outline specific details.
8. When multiple occlusal results are generated from the scan app, you can send all occlusal results at once by checking "Order All Occlusal Relationship."
9. Proceed with your order.
 - **Skip Payment**: Proceed with the order without using the payment system provided by Medit Link. Available when placing an order with a contracted partner.
 - **Pay**: Utilize the payment system provided by Medit Link to proceed with the order. Available when placing an order with a contracted partner.
 - **OK**: We proceed with the order without any conditions regarding the price of the product. Available when placing an order with a company you have partnered with

but do not have a contract with.


iNote

The following services are not available on Medit Link China.

- Skip payment
- Pay


Confirm Order (In Box)

- Accept the order upon receiving it.
- Accepted orders will be forwarded to the Work Box.




Cancel Order (Order Box)

- You can cancel an order before your partner accepts it.
- Canceled cases can be re-ordered.




Reject Order (In Box)

- You can reject a received order.
- Rejected cases can be re-ordered.



Checking Ordered Cases

You can check the details of the case you ordered from the Order Box in the clinic account. You can also check the case you received from the In Box in the lab account.



No.	Type	Method	Material	Shade	Price
23	Inlay	-	Zirconia	None Shade	1.00
24	Inlay	-	Zirconia	None Shade	1.00
25	Inlay	-	Zirconia	None Shade	1.00
26	Inlay	-	Zirconia	None Shade	1.00

Total USD 4.00

Order Info	The information that can be included in the order is as follows: order ID/patient name/account name/scan date/order date/delivery date/delivery due date/memo. You can check the details of the ordered product and delivery information.
Order List	The order list displays the details of the ordered product.
Shipment	Enter/check the shipping information.

Entering Shipment Information (In Box)

Enter the shipping information of the completed case and start shipping.

1. Click the Shipment tab.
2. Enter the shipment information and click "Register."

Shipment Registration

Select a delivery method and enter the required information to register.
The information you enter is saved in the Shipment tab and can be checked at any time.

☒ Shipment
 ☐ Self-Delivery

Company

Tracking Number

Memo

Cancel

Register

3. The entered information is sent to your partner.

Managing Delivery Date

You can set the expected delivery date when placing an order, and this information will be forwarded to your partner.

Setting Requested Delivery Date (Clinic account)

- When ordering from Case Box, you can enter the requested delivery date.

Order

General Order List File List

Partner
WonhoonChoi_Lab_Stage

WonhoonChoi_Lab_Stage
Seoul, KR
[02835] 23 Incheon-ro 22-gil, Seongbuk-gu, Seoul, Korea
01026088546

Case Name
mandibular movement - Clone - Clone

Patient Name
Test Occlusion

☒ Share patient information

Requested Delivery Date
2023-04-13

Delivery Time
오전 8:15 AM 8:00 PM 8:00

Memo

Tags
You can add the entered tag with Enter or Spacebar.

☒ Include all occlusion relationships in the order
☒ Include Tags In The Order

Cancel OK

- The default value is 7 days after the date of the order.
- Enter the date/time you wish to receive the products and proceed with the order.
- You can set the time easily using the slider.

Setting Expected Delivery Date (Clinic account)

- You can enter the expected delivery date for the case ordered from In Box.

MEDIT Link 3.0.0 Wonhoon Choi

test's Case / Wonhoonchoi-Clinic

Reject Confirm Order

Pending

Order ID
10413603

Patient Name
test

Partner Name
Wonhoonchoi-Clinic

Scanning Date
2022-09-29 오후 7:56

Order Date
2022-09-29 오후 8:13

Requested Delivery Date
2022-10-01 오후 10:30

Expected Delivery Date
2022-10-01 오후 10:30

Memo
-Scanbody Library-
16 : 0DENTIS_Custom_(NMR) ver 1.41 - Custom_Mini_Extra Long_16mm - N-HEX_T10
15 : 0DENTIS_Custom_(NMR) ver 1.41 - Custom_Mini_Extra Long_16mm - N-HEX_T10
14 : 0DENTIS_Custom_(NMR) ver 1.41 - Custom_Mini_Extra Long_16mm - N-HEX_T10

No.	Type	Method	Material	Shade	Price
16	Implant Crown	Cementation Ty...	Multi-Layered Zirconia	None Sh...	0
15	Implant Crown	Cementation Ty...	Multi-Layered Zirconia	None Sh...	0
14	Implant Crown	Cementation Ty...	Multi-Layered Zirconia	None Sh...	0

Total KRW 0

1. Check the requested delivery date with the partner and reschedule using the "Edit" button.
2. Enter the expected delivery date considering the actual working day. If required, you can leave a comment also.
3. Click "Apply."

March 2022 >

S	M	T	W	T	F	S
27	28	1	2	3	4	5
6	7	8	9	10	11	12
13	14	15	16	17	18	19
20	21	22	23	24	25	26
27	28	29	30	31	1	2
3	4	5	6	7	8	9

Delivery Time

11:45 AM AM 8:00 PM 8:00

Comment

Cancel Apply

4. The expected delivery date and the memo will be forwarded together with the order.
 - You can check the date and memos in the Order Box.
 - You can change the expected delivery date, and the changes will be forwarded to your partner.
5. You can check the case history using the "History" button.



Completing Case

Complete (Work Box/In Box)


Once all tasks are done, click "Done" to complete the case.

Complete a case from the case detail

You can find the "Menu" button on the right from the Open button.



Complete the work by clicking the "Menu" icon and selecting "Done."




Complete the multiple cases from the case list

When you select the cases from the Work Box, the expanded menu button is provided at the top of the list.



Complete multiple cases as follows:

1. Go to Work Box.
2. Select all cases to complete.



3. Click the above-mentioned "Expanded Menu" button.




4. Select "Done" from the list

Complete (Case Box/Order Box)

After the case completes from the lab, the clinic needs to confirm the case to complete the whole order process.

If the lab completes its work, the clinic can find the "Confirm Delivery" function in the Order Box.


Click this, the whole order process completes, and the case status will be changed to "Case Completed."



Complete a non-ordered case in the Case Box

For cases that are not ordered to the Lab, you can complete these cases from the Case Box.

1. Go to the Case Box.
2. Select the case that you want to complete.



3. Click the "Expanded Menu" button.



4. Choose "Case Completed."

Communicate Using Tags

You can add additional information for your case using tags and send them to your partners.

Add Tags When Creating a New Case

You can enter tags when you create a new case or edit case information.

New Case Information

Patient Name*

Case Name*

Creating patient information involves consent from the patient.
Patient information can be shared with 3rd parties while placing an order.

Tag

Doctorname × Specialmaterial × SpecialComment ×

You can add up to 10 tags of 100 characters or less.

[Terms & Conditions](#)

Cancel

Register

Register & Scan

You can easily search and select tags you've used before.

New Case Information

Patient Name*

Case Name*

Creating patient information involves consent from the patient.
Patient information can be shared with 3rd parties while placing an order.

Tag

t

t (new)

Suggestions

tet

tagtest

tag

tat2

Add Tags When Ordering

You can check and edit tags in the Order dialog while ordering.

Order

General Order List File List

Partner
WonhoonChoi_Lab_Stage

WonhoonChoi_Lab_Stage
Seoul, KR
[02835] 23 Incheon-ro 22-gil, Seongbuk-gu, Seoul, Korea
01026088546

Case Name
mandibular movement - Clone - Clone

Patient Name
Test Occlusion

☒ Share patient information

Requested Delivery Date
2023-04-13

Delivery Time
오전 8:15

Memo

Tags
Doctorname Specialmaterial SpecialComment

You can add up to 10 tags of 100 characters or less.

☒ Include all occlusion relationships in the order

☒ Include Tags In The Order

Cancel OK

Check Tags on Case Details

